

HALYARD*

Clear Procedure Drapes

CLEAR CHOICES. CLEAR COVERAGE. CLEAR SAVINGS.

CLEAR COVERAGE

HALYARD* Surgical Products are recognized for quality and innovation in protection for O.R. patients and staff. Now, in addition to our extensive line of surgical drapes and equipment covers, HALYARD* brings you a full selection of Clear Procedure Drapes.

CLEAR CHOICES

The right drapes are essential to maintaining a sterile field to protect surgical patients and reduce the risk of cross-contamination. HALYARD* Clear Procedure Drapes are:

- Transparent, impervious and drapable
- Strong, residue-free adhesive
- Non-linting release liners
- Low-glare matte finish

Usable alone or attachable to our non woven patient drapes for instrument or fluid management

CLEAR SAVINGS


Check out the following pages for details on our new range of Clear Procedure Drapes, and talk to your O&M HALYARD* Clinical Specialist about how your team can realize clear savings for your O.R.



Our
Commitment
to you

Better
Infection
Prevention


Clear Procedure Drapes



INCISE PROCEDURE DRAPE

- Clear Impervious Drapable Polyethylene Film
- Incise area
- 2.4" diameter circular fenestration



Sterile Code	Non-Sterile Code	Size/Description	Sterile Case Count	Non-Sterile Case Count
49050	49060NS	Clear Incise Drape w/ Fenestration, 40x40 cm	80	200



CLEAR U-DRAPES

- U-shaped Fenestration with Adhesive
- Clear, impervious, drapable polyethylene film
- Matte finish to reduce glare


Sterile Code	Non-Sterile Code	Size/Description	Sterile Case Count	Non-Sterile Case Count
49048	49058NS	Clear U-Drape, 130x120 cm	40	200
49049	49059NS	Clear Long U-Drape, 180x120 cm	40	200




TOWEL DRAPES

- Clear polyethylene film
- Adhesive strip
- Versatile, conformable drapes

Sterile Code	Non-Sterile Code	Size/Description	Sterile Case Count	Non-Sterile Case Count
49046	49056NS	Clear Towel Drape, 30 x45 cm	40	200
49047	49057NS	Clear Towel Drape, 60 x45 cm	40	300






49041 / 49051NS


FLUID COLLECTION POUCHES

- Collects and holds fluid during procedures
- Impervious Clear Drapable Polyethylene Film
- Adhesive Strips for Attachment
- Drainage Ports
- Filter Screen no 49042/49052NS Pouch to catch bone chips and tissue for specimen collection



49042 / 49052NS


Sterile Code	Non-Sterile Code	Size/Description	Sterile Case Count	Non-Sterile Case Count
49041	49051NS	Clear U-Pouch, 89x76 cm	20	50
49042	49052NS	Clear Irrigation Pouch, 60x50 cm	40	200



49043 / 49053NS


INSTRUMENT POUCHES

- Multiple Compartment ot hold Instruments
- Adhesive strips



49044 / 49054NS

Sterile Code	Non-Sterile Code	Size/Description	Sterile Case Count	Non-Sterile Case Count
49043	49053NS	Clear Instrument Pouch, 2 Compartments, 30x18 cm	40	200
49044	49054NS	Clear Instrument Pouch, 3 Compartments, 47x25 cm	40	200
49045	49055NS	Large Clear Instrument Pouch, 33x38 cm	200	336



49045 / 49055NS

