


DENTAL PRODUCTS CATALOGUE



GUESS WHO
CAME IN
WITH YOUR
PATIENTS! —

YOU DESERVE THE RIGHT PROTECTION —

CONTENTS

7	Facial Protection & Protective Eyewear
15	Exam Gloves
19	Sterilisation Wrap
23	Other Products

ABOUT HALYARD...

Over the course of our 150-year history, Halyard has built a reputation for moving care forward through innovation and improvement. **Today, healthcare professionals in over 90 countries use Halyard** gloves, facial protection, protection apparel, surgical gowns and drapes, sterilisation solutions and custom procedure trays to prevent infections and help to ensure the safety of their patients – and themselves.

Now, as part of the Global Products Division of Owens & Minor, a healthcare-oriented company with a long tradition of industry leadership, we are even better positioned to **offer best-in-class products and support.**


Dental professionals are increasingly challenged to provide infection control for their patients, while protecting themselves. As one of the world's leading suppliers of infection control products, you can count on Halyard for products you use every day in your practice including exam gloves, face masks, eyewear, sterilisation wrap and dental surgical drapes.

This catalogue details the wide array of products available to dental professionals.

For more information, visit
www.halyardhealth.co.uk



A Partnership with Halyard ensures:



- **Reliability:** reduce disruptions across the entire supply chain
- **Standardisation & Efficiency:** long-term cost savings in an increasingly complicated industry
- **Knowledge & Expertise:** empowering you to prevent infections across the continuum of care

WHY CHOOSE HALYARD AS YOUR DENTAL PARTNER IN INFECTION PREVENTION?



FACIAL PROTECTION & PROTECTIVE EYEWEAR



FACIAL PROTECTION



STANDARD PROCEDURE MASKS

- Standard masks for routine patient exams – no fluid exposure expected
- Pleat-style with earloops
- Meets EN 14683:2019+AC:2019 Type II standards for:
 - Bacterial Filtration Efficiency (BFE \geq 98%)
 - Differential Pressure ($<$ 40 Pa/cm²)
 - Microbial Cleanliness (\leq 30 cfu/g)

Item Code	Description	Type	Items/Box	Boxes/Case
47080	Standard Procedure Mask, Blue	I	50	10
47117	Standard Procedure Mask, Yellow	I	50	10
47295	TEDDY BEAR* Procedure Mask, Pink and Blue Teddy Bear Print	II	50	10
62356	THE LITE ONE* Procedure Mask, Blue	II	50	10

FACIAL PROTECTION



SPECIALTY SURGICAL FACE MASKS

- Standard masks for routine patient exams – no fluid exposure expected
- Pleat-style with ties
- Meets EN 14683:2019+AC:2019 Type II standards for:
 - Bacterial Filtration Efficiency (BFE \geq 98%)
 - Differential Pressure ($<$ 40 Pa/cm²)
 - Microbial Cleanliness (\leq 30 cfu/g)

Item Code	Description	Type	Items/Box	Boxes/Case
49310	Fog-Free ULTI-MASK [®] Surgical Mask, Foam Band, White	II	50	6
49214	Fog-Free Surgical Mask, Foam Band, Blue	II	50	6
49215	Anti-Fog Surgical Mask, DERMA-TOUCH [®] Tape, Green	II	50	6
49235	Anti-Fog Surgical Mask, Film Strip, Green	II	50	6
48390	SO SOFT [®] Surgical Mask, White	II	50	6

FACIAL PROTECTION



FLUIDSHIELD* LEVEL 1 PROCEDURE MASKS

- Meets EN 14683:2019+AC:2019 Type II standards for:
 - Bacterial Filtration Efficiency (BFE \geq 98%)
 - Microbial cleanliness (\leq 30 cfu/g)
- Meets ASTM F2100-11 Level 1 standards for fluid resistance, filtration and breathability
- Assured fluid protection for routine patient care when low level of exposure to fluids expected
- Pleat-style with earloops with inner layer of SO SOFT* fabric for comfort

Item Code	Description	Type	Items/Box	Boxes/Case
25869	FLUIDSHIELD* Level 1 Procedure Mask, Blue	-	50	10
25867	FLUIDSHIELD* 1 Procedure Mask with SO SOFT* Lining and SO SOFT* Earloops, Yellow	-	50	10



FLUIDSHIELD* LEVEL 1 SURGICAL MASKS

- Meets EN 14683:2019+AC:2019 Type II standards for:
 - Bacterial Filtration Efficiency (BFE \geq 98%)
 - Microbial cleanliness (\leq 30 cfu/g)
- Meets ASTM F2100-11 Level 1 standards for fluid resistance, filtration and breathability
- Assured fluid protection for routine patient care when low level of exposure to fluids expected
- Pleat-style with ties with inner layer of SO SOFT* fabric for comfort

Item Code	Description	Type	Items/Box	Boxes/Case
28806	FLUIDSHIELD* Level 1 Fog-Free Surgical Mask, Green	-	50	6

FACIAL PROTECTION



FLUIDSHIELD* LEVEL 2 PROCEDURE MASKS

- Meets EN 14683:2019+AC:2019 Type IIR standards for:
 - Bacterial Filtration Efficiency (BFE \geq 98%)
 - Differential Pressure ($<$ 60 Pa/cm²)
 - Splash Resistance Pressure (\geq 16,0 kPa)
 - Microbial Cleanliness (\leq 30 cfu/g)
- For protection from moderate risk of fluid splash and spray
- Meets ASTM F2100-11 Level 2 standards for fluid resistance (120 mm Hg), filtration and breathability
- Pleat-style with earloops

Item Code	Description	Type	Items/Box	Boxes/Case
62115	FLUIDSHIELD* 2 Splash Resistant Procedure mask, Blue with Orange Diamond Print	IIR	50	10
62116	FLUIDSHIELD* 2 Fog-Free Splash Resistant Procedure Mask, Wraparound Splashguard Visor, Blue	IIR	25	4



FLUIDSHIELD* LEVEL 2 SURGICAL MASKS

- Meets EN 14683:2019+AC:2019 Type IIR standards for:
 - Bacterial Filtration Efficiency (BFE \geq 98%)
 - Differential Pressure ($<$ 60 Pa/cm²)
 - Splash Resistance Pressure (\geq 16,0 kPa)
 - Microbial cleanliness (\leq 30 cfu/g)
- For protection from moderate risk of fluid splash and spray
- Meets ASTM F2100-11 Level 2 standards for fluid resistance (120 mm Hg), filtration and breathability
- Pleat-style with ties

Item Code	Description	Type	Items/Box	Boxes/Case
62113	FLUIDSHIELD* 2 Fog-Free, Splash Resistant Surgical Mask, Blue with Orange Diamond Print	IIR	50	6
62114	FLUIDSHIELD* 2 Fog-Free, Splash Resistant Surgical mask, Wraparound Splashguard Visor, Blue with Orange Diamond Print	IIR	25	4

FACIAL PROTECTION



FLUIDSHIELD* LEVEL 3 PROCEDURE MASKS

- Meets EN 14683:2019+AC:2019 Type IIR requirements for:
 - Bacterial Filtration Efficiency (BFE \geq 98%)
 - Differential Pressure ($<$ 60 Pa/cm²)
 - Splash Resistance Pressure (\geq 16,0 kPa)
 - Microbial Cleanliness (\leq 30 cfu/g)
- Protects from fluid splattering and splashing
- Meets ASTM F2100-11 Level 3 standards for fluid resistance (160 mm Hg), filtration and breathability
- Pleat-style with earloops with inner layer of SO SOFT* fabric for comfort

Item Code	Description	Type	Items/Box	Boxes/Case
47107	SO SOFT* FLUIDSHIELD* Fog-Free Procedure Mask, Foam Band, Orange	IIR	40	10
28797	SO SOFT* FLUIDSHIELD* Fog-Free Procedure Mask, Foam Band, Orange	IIR	40	10
47137	SO SOFT* FLUIDSHIELD* Fog-Free Procedure Mask w/SPLASHGUARD* Visor, Foam Band, Orange	IIR	25	4
47147	SO SOFT* FLUIDSHIELD* Fog-Free Procedure Mask w/Wraparound SPLASHGUARD* Visor, Foam Band, Orange	IIR	25	4



FLUIDSHIELD* LEVEL 3 SPECIALTY PROCEDURE MASKS

- Meets EN 14683:2019 + AC:2019 Type IIR requirements for:
 - Bacterial Filtration Efficiency (BFE \geq 98%)
 - Differential Pressure ($<$ 60 Pa/cm²)
 - Splash Resistance Pressure (\geq 16,0 kPa)
 - Microbial Cleanliness (\leq 30 cfu/g)
- Protects from fluid splattering and splashing
- Meets ASTM F2100-11 Level 3 standards for fluid resistance (160 mm Hg), filtration and breathability
- Pleat-style with earloops with inner layer of SO SOFT* fabric for comfort

Item Code	Description	Type	Items/Box	Boxes/Case
47297	FLUIDSHIELD* TEDDY BEAR* Procedure Mask, Pink and Blue Teddy Bear Print	IIR	40	10
47298	FLUIDSHIELD* Fog-Free TEDDY BEAR* Procedure Mask w/Wraparound SPLASHGUARD* Visor, Foam Band, Pink and Blue Teddy Bear Print	IIR	25	4

FACIAL PROTECTION



FLUIDSHIELD* LEVEL 3 SURGICAL MASKS

- Meets EN 14683:2019 + AC:2019 Type IIR requirements for:
 - Bacterial Filtration Efficiency (BFE \geq 98%)
 - Differential Pressure ($<$ 60 Pa/cm²)
 - Splash Resistance Pressure (\geq 16,0 kPa)
 - Microbial Cleanliness (\leq 30 cfu/g)
- Protects from fluid splattering and splashing
- Meets ASTM F2100-11 Level 3 standards for fluid resistance (160 mm Hg), filtration and breathability
- Pleat-style with ties with inner layer of SO SOFT* fabric for comfort



Item Code	Description	Type	Items/Box	Boxes/Case
48207	SO SOFT* FLUIDSHIELD* Fog-Free Surgical Mask, Foam Band, Orange	IIR	50	6
48208	SO SOFT* FLUIDSHIELD* Surgical Mask, Orange	IIR	50	6
48237	SO SOFT* FLUIDSHIELD* Fog-Free Surgical Mask w/SPLASHGUARD* Visor, Foam Band, Orange	IIR	25	4
48247	SO SOFT* FLUIDSHIELD* Fog-Free Surgical Mask w/Wraparound SPLASHGUARD* Anti-Glare Visor, Foam Band, Orange	IIR	25	4



FLUIDSHIELD* LEVEL 3 SPECIALTY SURGICAL MASK

- Meets EN 14683:2019 + AC:2019 Type IIR requirements for:
 - Bacterial Filtration Efficiency (BFE \geq 98%)
 - Differential Pressure ($<$ 60 Pa/cm²)
 - Splash Resistance Pressure (\geq 16,0 kPa)
 - Microbial Cleanliness (\leq 30 cfu/g)
- Protects from fluid splattering and splashing
- Meets ASTM F2100-11 Level 3 standards for fluid resistance (160 mm Hg), filtration and breathability
- Pleat-style with ties with inner layer of SO SOFT* fabric for comfort

Item Code	Description	Type	Items/Box	Boxes/Case
48297	SO SOFT* FLUIDSHIELD* TEDDY BEAR* Surgical Mask, Pink and Blue Teddy Bear Print	IIR	50	6

FACIAL PROTECTION



CHILD FACE MASK

- Teddy Bear Child's Mask with SO-SOFT* Earloops
- Soft, lightweight materials for cool comfort
- Comfortable knitted earloops
- Formable nose wire for a custom fit on small faces

Item Code	Description	Items/Box	Boxes/Case
47359	Child Face Mask	75	10



SAFEVIEW* GLASSES – DISPOSABLE EYEWEAR

- Protection from splash, splatter or spray
- Multiple colours

Item Code	Description	Items/Box	Boxes/Case
SV50A	Assembled Glasses	10	5
SV100F	Replacement Frames	10	10
SV250L	Replacement Lenses	25	10



EXAM GLOVES



EXAM GLOVES



STERLING* NITRILE EXAM GLOVES

- For procedures with low-to-moderate risk of exposure to fluid
- Flexible with high tactile sensitivity
- Non-sterile • Ambidextrous • Textured fingertips • Beaded cuffs
- Provides nitrile protection while avoiding latex
- Powder-free and natural rubber latex-free • Comply to EN 455, EN 374 & EN 21420
- Acceptable Quality Level (AQL) < 1,0 • PPE Cat III & food handling approved
- Average length > 240 mm

Item Code	Description	Size	Items/Box	Boxes/Case
13938	STERLING* NITRILE Exam Gloves	X-Small	200	10
13940	STERLING* NITRILE Exam Gloves	Small	200	10
13941	STERLING* NITRILE Exam Gloves	Medium	200	10
13942	STERLING* NITRILE Exam Gloves	Large	200	10
13943	STERLING* NITRILE Exam Gloves	X-Large	170	10



STERLING* NITRILE-XTRA* EXAM GLOVES

- For more protection against fluid and splashing
- Flexible with high tactile sensitivity • Nitrile protection while avoiding latex
- Non-sterile • Ambidextrous • Textured fingertips
- Extended, beaded cuff (average length 310 mm) for better coverage over the wrist
- Powder-free and natural rubber latex-free • Comply to EN 455, EN 374 & EN 21420
- Acceptable Quality Level (AQL) < 1,0 • PPE Cat III & food handling approved

Item Code	Description	Size	Items/Box	Boxes/Case
44286	STERLING* NITRILE-XTRA* Exam Gloves	X-Small	100	10
44287	STERLING* NITRILE-XTRA* Exam Gloves	Small	100	10
44288	STERLING* NITRILE-XTRA* Exam Gloves	Medium	100	10
44289	STERLING* NITRILE-XTRA* Exam Gloves	Large	100	10
44290	STERLING* NITRILE-XTRA* Exam Gloves	X-Large	100	10

EXAM GLOVES



PURPLE NITRILE* EXAM GLOVES

- Proven barrier protection from bacteria, bloodborne viruses and chemicals during procedures where fluid exposure risk is moderate to high
- Available in both sterile and non-sterile versions
- Average length > 240 mm
- Powder-free and natural rubber latex-free • Comply to EN 455, EN 374 & EN 21420
- Acceptable Quality Level (AQL) < 1,0 PPE Cat III

Item Code	Description	Size	Items/Box	Boxes/Case
52000M	PURPLE NITRILE* Exam Gloves	X-Small	100	10
52001M	PURPLE NITRILE* Exam Gloves	Small	100	10
52002M	PURPLE NITRILE* Exam Gloves	Medium	100	10
52003M	PURPLE NITRILE* Exam Gloves	Large	100	10
52004M	PURPLE NITRILE* Exam Gloves	X-Large	90	10



PURPLE NITRILE-XTRA* EXAM GLOVES

- Proven barrier protection from bacteria, bloodborne viruses and chemicals during procedures where fluid exposure risk is moderate to high
- Extended cuff (average length 310 mm) for better coverage over the wrist
- Powder-free and natural rubber latex-free • Comply to EN 455, EN 374 & EN 21420
- Acceptable Quality Level (AQL) < 1,0 PPE Cat III

Item Code	Description	Size	Items/Box	Boxes/Case
50600M	PURPLE NITRILE-XTRA* Exam Gloves	X-Small	50	10
50601M	PURPLE NITRILE-XTRA* Exam Gloves	Small	50	10
50602M	PURPLE NITRILE-XTRA* Exam Gloves	Medium	50	10
50603M	PURPLE NITRILE-XTRA* Exam Gloves	Large	50	10
50604M	PURPLE NITRILE-XTRA* Exam Gloves	X-Large	50	10



PURPLE NITRILE* STERILE PAIRS EXAM GLOVES

Item Code	Description	Size	Items/Box	Boxes/Case
52201M	PURPLE NITRILE* Sterile Exam Gloves	Small	50	4
52202M	PURPLE NITRILE* Sterile Exam Gloves	Medium	50	4
52203M	PURPLE NITRILE* Sterile Exam Gloves	Large	50	4

EXAM GLOVES




SKYBREEZE* NITRILE EXAM GLOVES

- Suitable glove choice for a wide range of tasks:
- Tested with various chemotherapy drugs and several commonly used chemicals
- Skin friendly formulations and ISO 10993-10 compliant gloves offer excellent comfort
- Anatomical thin design for enhanced tactile sensitivity; Textured fingertips for better grip
- Viral barrier protection, non sterile
- Powder-free and Natural rubber latex-free formulation reduces the likelihood of skin irritation and Type I allergies
- SMARTPULL* dispenser facilitates glove dispensing and reduces glove waste by 38.2%

Item Code	Description	Size	Items/Box	Boxes/Case
47373	SKYBREEZE* Powder-free Nitrile Exam Gloves	X-Small	200	10
47374	SKYBREEZE* Powder-free Nitrile Exam Gloves	Small	200	10
47375	SKYBREEZE* Powder-free Nitrile Exam Gloves	Medium	200	10
47376	SKYBREEZE* Powder-free Nitrile Exam Gloves	Large	200	10
47377	SKYBREEZE* Powder-free Nitrile Exam Gloves	X-Large	180	10

STERILISATION WRAP




STERILISATION WRAP

HALYARD* SEQUENTIAL STERILISATION WRAP H100







- Sequential Wrap
- Low-Lint SMS Fabric
- Steam and EO Compatible
- Compatible with Low Temperature Gas Plasma Sterilisation
- Compliant with EN ISO 11607 and EN 868-2



Item Code	Size	Quantity
10712	30 cm x 30 cm	1000 sheets/case
10715	38 cm x 38 cm	1000 sheets/case
10718	45 cm x 45 cm	1000 sheets/case
10720	50 cm x 50 cm	1000 sheets/case
10724	60 cm x 60 cm	500 sheets/case
10730	76 cm x 76 cm	300 sheets/case
10736	91 cm x 91 cm	300 sheets/case
10740	101 cm x 101 cm	250 sheets/case
10745	114 cm x 114 cm	250 sheets/case
10748	121 cm x 121 cm	250 sheets/case
10754	137 cm x 137 cm	100 sheets/case
10772	137 cm x 182 cm	100 sheets/case


STERILISATION WRAP



HALYARD* QUICK CHECK* STERILISATION WRAP H100

- Simultaneous (Non-Sequential) Wrap
- One Application Equivalent to 2 Sheets of Wrap: Dual Colouring
- Low-Lint SMS Fabric
- Steam and EO Compatible
- Compatible with Low Temperature Gas Plasma Sterilisation
- Compliant with EN ISO 11607 and EN 868-2

Item Code	Size	Quantity
34156	45 x 45 cm	480 applications/case
34192	50 x 50 cm	480 applications/case
34178	60 x 60 cm	240 applications/case
34176	76 x 76 cm	144 applications/case
34194	91 x 91 cm	144 applications/case
34200	101 x 101 cm	120 applications/case
34161	121 x 121 cm	120 applications/case
34172	137 x 137 cm	96 applications/case



OTHER PRODUCTS



APPAREL

PROFESSIONAL JACKETS



- Protective 3-layer SMS fabric
- Protective in situations with potential for moderate fluid contact
- Knit collar and cuffs
- 5 snaps and 2 reinforced pockets

Item Code	Description	Size	Items/Case
10067	Professional Jacket, Blue	Medium	24
10078	Professional Jacket, Blue	Large	24
10087	Professional Jacket, Blue	X-large	24

SCRUBSUITS



SCRUB SHIRTS

- Protective 3-Layer SMS Fabric
- Potential for Light Fluid Contact
- Scoop Neck
- Reinforced Pockets
- Sealed Seams

Item Code	Description	Size	Items/Case
69701	Scrub Shirt, Blue	Medium	48
69702	Scrub Shirt, Blue	Large	48
69703	Scrub Shirt, Blue	X-large	48



SCRUB PANTS

- Protective 3-Layer SMS Fabric
- Potential for Light Fluid Contact
- Drawstring Waist
- Side and Rear Pockets
- Reinforced In-Seam

Item Code	Description	Size	Items/Case
69711	Scrub Pants, Blue	Medium	48
69712	Scrub Pants, Blue	Large	48
69713	Scrub Pants, Blue	X-large	48

APPAREL

PROTECTIVE GOWNS



PROCEDURE GOWNS WITH KNIT CUFFS

- Protective 3-Layer SMS Fabric
- Potential for Light to Moderate Fluid Contact
- Knit Cuff
- VELCRO® Neck Closure
- Tie Waist
- Full Back


Item Code	Description	Size	Items/Bundle	Items/Case
69025	Procedure Gown w/Knit Cuffs, Blue	Universal	10	6
69028	Procedure Gown w/Knit Cuffs, Blue	X-Large	10	6



IMPERVIOUS COMFORT GOWNS

- Impervious Film Laminate
- Potential for Soaking and Heavy Fluid Contact in Front
- Perforated Open Back

Item Code	Description	Size	Items/Bundle	Items/Case
69600	Impervious Comfort Gown, Knit Cuffs, Blue	Universal	10	10
69601	Impervious Comfort Gown, Knit Cuffs, Blue	XX-Large	10	10
69602	Impervious Comfort Gown, Thumbhooks, Blue	Universal	10	10
69603	Impervious Comfort Gown, Thumbhooks, Blue	XX-Large	10	10



IMPERVIOUS FILM GOWNS

- Potential for Soaking and Heavy Fluid Contact in Front
- Thumbhooks
- Perforated Open Back

Item Code	Description	Size	Items/Bundle	Items/Case
69316	Impervious Film Gown, Blue	X-Large	12	5
69490	Impervious Film Gown, Blue	Universal	15	5

ICE PACKS



BILATERAL FACIAL ICE PACKS

- Consistent cold therapy
- 3-layer construction
- Minimises condensation
- Soft outer cover
- Single ice chamber for easy filling
- Conforms to upper face and eyes

Item Code	Description	Size	Items/Box	Items/Case
33101	Bilateral Facial Ice Pack, 1 strap	12.5 x 30 cm	12	2

GOWNS AND DRAPES



EVOLUTION* 4 SURGICAL GOWNS

Item Code	Description	Size	Items/Case
90002	EVOLUTION* 4 Surgical Gown, Sterile	Small	36
90012	EVOLUTION* 4 Surgical Gown, Sterile	Large	36
90042	EVOLUTION* 4 Surgical Gown, Sterile	X-Large	34

STERILE DRAPES

Item Code	Description	Size	Items/Case
89543	HALYARD* BASICS* Dental Drape	120 x 150 cm	15

**For more information,
please contact your O&M Halyard Customer Care**

Team:

UK & IE (toll-free)

Tel: 0800 085 82 04 - Fax: 0800 066 91 06

customerservice.uk.ie@hyh.com

South Africa (toll-free)

Tel: 0801 111 332 - Fax: 0801 111 335

customerservice.za@hyh.com

Other countries

Tel: +32 2 585 51 14 - Fax: +32 2 897 00 95

customerservice.export@hyh.com

* Registered Trademark or Trademark of O&M Halyard or its affiliates. The colour GREY is a Registered Trademark or Trademark of O&M Halyard or its affiliates. The colour PURPLE is a Registered Trademark or Trademark of O&M Halyard or its affiliates. The colour ORANGE is a Registered Trademark or Trademark of O&M Halyard or its affiliates.
©2019. All rights reserved. HC135-04-UK