

OUR EXPERTISE TO GUARD AGAINST INFECTIONS IN THE OR

**Our
Commitment
to you**

**Better
Infection
Protection**



SOLUTIONS FOR SURGICAL SITE INFECTIONS

O&M Halyard offers a comprehensive range of surgical solutions designed to help reduce surgical site infections ranging from Surgical Drapes & Gowns, Face Masks to Sterilisation Wrap.



DEPEND ON HALYARD EXPERTISE TO HELP GUARD AGAINST INFECTIONS IN THE OR

When you choose HALYARD* Surgical Gowns, Drapes & Packs, you can be assured that you are selecting products that use innovative fabric technology and design to provide excellent strike-through resistance along with comfort.

O&M Halyard is committed to providing quality products that help reduce the toll infections take on patients, clinicians and facilities, assisting you in delivering optimal patient care.


Because we are as serious as you are about preventing infection in the OR, we have expanded our surgical gown and drape portfolio into multiple performance levels, to bring you a wide choice of protection without compromise on quality.

The products referenced in this catalogue do not contain natural rubber latex, nor is this material intentionally added during the manufacturing process.

ADVANCED PROTECTION THROUGH SUPERIOR SMS

Halyard offers a extensive portfolio of Surgical Gowns, Drapes and Sterilisation Wrap based upon the HALYARD* SMS technology (Spunbond, Meltblown, Spunbond). SMS is made up of three layers of polypropylene:

- The Spunbond outer layers provide exceptional strength and durability
- The Meltblown middle layer provides an effective barrier to bacteria and liquids
- The material has a low lint level, helping reduce the risk of contamination in the OR



A WIDE PORTFOLIO OF GOWNS AND DRAPES DESIGNED TO OFFER THE RIGHT PROTECTION FOR DIFFERENT PROCEDURES

Cost-effective solution for standard performance & protection

- For short duration, low fluid-intensive procedures
- For low risk, day & general surgery (minor procedures)

Basic Surgical Gowns

- HALYARD BASICS* Standard gowns
- Standard SMS fabric

HALYARD BASICS* Surgical Drapes & Packs

- Flexible, easy-to-use designs
- Reinforcement

Best-value for enhanced performance & protection

- For low-to-moderate fluid-intensive procedures
- For low-to-moderate risk, general surgery

Standard Surgical Gowns

- EVOLUTION* 4 Standard gowns
- ULTRA Standard gowns
- EVOLUTION* 4 & ULTRA enhanced SMS fabric

Standard Surgical Drapes & Packs

- Wide selection for hundreds of different procedures
- CONTROL* Plus Reinforcement

Cost-effective for high performance & protection

- For more fluid-intensive procedures
- For moderate-to-high risk, specialised surgery

Reinforced Surgical Gowns

- ULTRA fabric-reinforced gowns

Reinforced Surgical Drapes & Packs

- Including ORTHOARTS* and CVARTS* drapes
- CONTROL* Plus Reinforcement

Solution for high performance & perfection

- For long duration, high fluid-intensive procedures
- For complex, high risk specialised surgery

Protective Surgical Gowns

- AeroChrome*
- Cool Shield Core Technology, Aero Vent back
- Soft fabric S/F/SMS

HALYARD* CPTs ARE PROFESSIONALLY TAILORED TO HELP FIT YOUR NEEDS¹

PRECISE

From fully customised trays to excellent value trays with functional equivalent components that helps meet your specifications. Based on your vision, our team can also help with tray standardisation.

PROFESSIONAL

Assembled by CPT professionals at specialised CPT facilities, ensuring logical tray set-up and aseptic presentation in the OR.



PREFERRED

HALYARD* Drapes and Gowns deliver an excellent balance of properties to keep your staff and patients protected in the OR:

- Fluid resistance
- Flame resistance
- Lint and abrasion resistance

¹. Not available for the UK market.

TABLE OF CONTENTS

CONTENTS

ARCROYAL® CUSTOM PROCEDURE TRAYS	10
CUSTOM PROCEDURE TRAYS	10
GENERAL SURGERY	15
GENERAL SURGERY DRAPES	28
FENESTRATED DRAPES	28
U-DRAPES	32
BAR DRAPES	35
BOTTOM DRAPES	36
SIDE DRAPES	38
TOP DRAPES	39
ABSORBENT UNDERLEG DRAPES	40
SINGLE DRAPES	40
UTILITY DRAPES	42
LAPAROTOMY DRAPES	46
BREAST/CHEST DRAPES	48
LAPAROSCOPY DRAPES	49
ORTHOPAEDICS	50
ORTHOPAEDIC PACKS	50
ORTHOPAEDIC DRAPES	62
ORTHOARTS® SPLIT / BAR DRAPES	62
IMPERVIOUS SPLIT DRAPES	62
EXTREMITY DRAPES	63
LOWER EXTREMITY DRAPES	63
UPPER EXTREMITY DRAPES	65
SHOULDER DRAPES	65
HAND DRAPES	66
HIP DRAPES	66
VERTICAL ISOLATION DRAPES	68
KNEE ARTHROSCOPY DRAPES	68
UNIVERSAL SPINE DRAPES	69
U-DRAPES	69

GYNAECOLOGY & OBSTETRICS	70
GYNAECOLOGY & OBSTETRICS PACKS	70
GYNAECOLOGY & OBSTETRICS DRAPES	76
GYNAECOLOGY/UROLOGY DRAPES	76
C-SECTION DRAPES	76
LITHOTOMY DRAPES	78
UROLOGY	79
UROLOGY PACKS	79
UROLOGY PACKS & DRAPES	80
UROLOGY DRAPES	81
OPHTHALMOLOGY	82
OPHTHALMOLOGY PACKS	82
OPHTHALMOLOGY DRAPES	83
ENT	86
ENT PACKS	86
ENT DRAPES	88
ENT DRAPES	88
HEAD DRAPES	88
DENTAL DRAPES	88
NEUROLOGY	89
CRANIOTOMY DRAPES	89
CARDIOVASCULAR & ANGIOGRAPHY	90
CARDIOVASCULAR & ANGIOGRAPHY PACKS	90
CARDIOVASCULAR DRAPES	93
CARDIOVASCULAR & ANGIOGRAPHY DRAPES	94
ANGIOGRAPHY DRAPES	97

CATHETER PLACEMENT	103
CATHETER PLACEMENT DRAPES	103
PAEDIATRICS	104
PAEDIATRICS PACKS	104
PAEDIATRICS DRAPES	108
FENESTRATED DRAPES	108
BREAST/CHEST DRAPES	109
ANGIOGRAPHY DRAPES	109
BABY BLANKETS	109
EQUIPMENT COVERS	110
MAYO STAND COVERS	110
STOOL COVERS	110
BACK TABLE COVERS	111
ACCESSORIES	112
UNDER BUTTOCKS DRAPES	112
CLEAR DRAPES	112
PERINEAL COVER DRAPES	113
IMPERVIOUS STOCKINETTES	113
ARMBOARD COVERS	114
LEGGINGS	114
TOWELS	114
TAPE STRIPS	115
SLEEVES	115
POUCHES	115
FOOT COVERS	116
PITCHER	116
ROOM TURNOVER KITS	117

HALYARD BASICS*	118
HALYARD BASICS* SURGICAL GOWNS	118
EVOLUTION* 4	119
EVOLUTION* 4 SURGICAL GOWNS	119
ULTRA	121
ULTRA SURGICAL GOWNS	121
VALUESELECT*	123
VALUESELECT* SURGICAL GOWNS	123
AERO CHROME*	124
AERO CHROME* SURGICAL GOWNS	124
OPERATING ROOM SAFETY PRODUCTS	125
SAFETY PRODUCTS	125
LIGHT HANDLE COVERS	125
NEEDLE COUNTERS	126
INSTRUMENT PADS	127
EQUIPMENT COVERS	128
SAFETY ACCESSORIES	131
ENDOSCOPY, LAPAROSCOPY & ELECTROSURGICAL PRODUCTS	133
TEMPERATURE MANAGEMENT	134
HOLLOWARE MULTIPURPOSE BOWLS, BASINS AND TRAYS	135
INDEX BY CODE	138

ARCROYAL® CUSTOM PROCEDURE TRAYS

CUSTOM PROCEDURE TRAYS

Owens & Minor's ArcRoyal® Custom procedure trays include all single-use medical devices needed for a specific surgical procedure.

As a brand neutral supplier, we work with all the leading manufacturers of CPT components across the globe to ensure we have the broadest range of components and the highest quality components available, whatever the procedure. You, the customer, can choose every single component contained within an ArcRoyal® CPT.



WHY ARCROYAL® CPTS ARE DIFFERENT:

WE DELIVER GENUINE SERVICE

Our Promise to our Customers:

- No Substitution of Branded Products - No Surprises in the Operating Room!
- Efficient Revision Turnaround - Reductions in Non-Use and Wastage
- World Class Quality & IT Support - Quality you can rely on time & time again

Service is in our DNA

WE DELIVER REAL CHOICE

Customer Choice:

- We have the broadest range of components
- We offer multiple drapes and gown solutions
- Our global sourcing team is dedicated to locating high quality value options
- We are a brand neutral supplier

Your Tray. Your Way

WE DELIVER ABSOLUTE QUALITY

Quality First & Always:

- Proprietary Systems: VidViz - Visual Bill of Materials
- Highest Standards of Cleanroom Control - Full gowning procedures in place
- Meticulous Attention to Detail - Components individually approved for inclusion

We build quality into everything we do

WE DELIVER ON TIME, EVERY TIME

Our Culture:

Owens & Minor teammates go the extra mile for the customer, doing everything in our power to ensure the product is made right, packed right, shipped right and delivered right - and is on time, every time through our network of distribution centres.

Customer satisfaction is paramount



CUSTOM PROCEDURE TRAYS

COMPONENT CHOICE:

With an ArcRoyal® Custom procedure tray you have greater choice, as you decide exactly on your component mix, including the brand of drapes and gowns you prefer.



FIELD OF APPLICATION

At Owens & Minor, we offer both surgical and minor Custom Procedure Trays. These span the following areas:

Robotics, Orthopaedics, Cardiology/Vascular, Cath Lab, Labour and Delivery, Laparoscopic, Gastric, Trauma, General Surgery, Neurology, Plastics, ENT, Angiography, Ophthalmology, and Gynaecology.



ARCROYAL® CPT'S - MANUFACTURED TO HIGHEST QUALITY STANDARDS

Owens & Minor is committed to providing our customers with quality products and services in a timely and cost effective manner. To achieve this, we constantly strive to measure, analyse and improve quality in all aspects of our business, through a comprehensive Quality Management System developed to comply with the requirements of:

- ISO/EN 13485:2016 standards
- Annex V and Annex VII of the Medical Device Directive (93/42/EEC)
- All the relevant regulations of EMEA countries into which we supply our CPTs

Each individual component included in our CPT's is CE-marked and subject to a rigorous quality and testing procedure.

We have an extensive system of traceability in place, with every pack and every individual component therein traceable back to a specific vendor batch.

Within our Class 8 cleanroom facility we have developed a proprietary visual Bill of Materials IT system which ensures that all components are included in each pack and assembled exactly according to customer specification and preference.

GUARANTEED STERILITY MAINTENANCE OF ARCROYAL® CPT'S

ArcRoyal® CPTs are sterilised after they have been assembled, boxed and prepared for final shipping. Sterility is maintained by specialised Tyvek packaging according to the packaging standard ISO 11607-1:2020. Our sterilisation system is compliant with ISO 11135-1:2014, and all chambers used for product sterilisation across various facilities have been validated.

Owens & Minor controls all transport to and from the steriliser, with all steriliser and freight suppliers subject to our strict validation processes.



CUSTOM PROCEDURE TRAYS

PROPERTIES

- ArcRoyal® Custom procedure trays contain a selection of customised single-use medical devices required for a specific surgical procedure.
- All trays are assembled exactly to clinician and hospital specifications and open sequentially as required.
- ArcRoyal® Custom procedure trays are presented in surgically compliant packaging which is compact and storage-friendly.
- All ArcRoyal® Custom procedure trays are sterile and ready for immediate use.



BENEFITS

- Reduced set-up time: fewer products to open, freeing up theatre staff for other duties.
- More efficient product ordering: fewer products to order, resulting in more streamlined ordering and stock control processes.
- Waste Reduction: less packaging to dispose.
- Cost Saving: few component costs per procedure to track, resulting in more efficient cost control.
- Greater Safety: Presented in an easy sterilised format to reduce the risk of contamination.
- Quality Certification: All ArcRoyal® Custom Procedure Trays are ISO/EN 13485:2003 and ISO 9001:2008 certified, and fully compliant with European Medical Device Directives.
- Greater Traceability: Every component used in an ArcRoyal Custom Procedure Tray can be traced back to a specific vendor batch.

GENERAL SURGERY


GENERAL SURGERY PACKS

UNIVERSAL PACK, SMALL

CODE	DESCRIPTION		ITEMS/CASE
88762	Universal Pack, Small, incl. Mayo Stand Cover		11
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 96 x 96 cm
89603		79603	1 Mayo Stand Cover 80 x 145 cm
42222		42428	1 Back Table Cover, Reinforced 152 x 191 cm
89053		79053	1 Tape Strip 10 x 55 cm
89701		79720	1 Absorbent Towel 38 x 56 cm
-	-	-	1 Suture Bag
-	-	-	2 Side Drapes 90 x 75 cm: • Positioning Tape, 2.5 cm width • CONTROL [®] Fabric Reinforcement
-		79214	1 Bottom Drape 149 x 180 cm: • Positioning Tape, 2.5 cm width • CONTROL [®] Fabric Reinforcement
-		79215	1 Top Drape 149 x 180 cm: • Positioning Tape, 2.5 cm width • CONTROL [®] Fabric Reinforcement


UNIVERSAL PACK, SMALL

CODE	DESCRIPTION		ITEMS/CASE
88763	Universal Pack, Small, excl. Mayo Stand Cover		11
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 96 x 96 cm
42222		42428	1 Back Table Cover, Reinforced 152 x 191 cm
89053		79053	1 Tape Strip 10 x 55 cm
89701		79720	1 Absorbent Towel 38 x 56 cm
-	-	-	1 Suture Bag
-	-	-	2 Side Drapes 90 x 75 cm: • Positioning Tape, 2.5 cm width • CONTROL [®] Fabric Reinforcement
-		79214	1 Bottom Drape 149 x 180 cm: • Positioning Tape, 2.5 cm width • CONTROL [®] Fabric Reinforcement
-		79215	1 Top Drape 149 x 180 cm: • Positioning Tape, 2.5 cm width • CONTROL [®] Fabric Reinforcement



All codes in **bold** indicate a sterile product.

GENERAL SURGERY PACKS



UNIVERSAL PACK, LARGE

CODE	DESCRIPTION		ITEMS/CASE
88764	Universal Pack, Large, incl. Mayo Stand Cover		10
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 96 x 96 cm
89603		79603	1 Mayo Stand Cover 80 x 145 cm
42222		42428	1 Back Table Cover, Reinforced 152 x 191 cm
89053		79053	1 Tape Strip 10 x 55 cm
89701		79720	1 Absorbent Towel 38 x 56 cm
	-	-	1 Suture Bag
	-	79216	2 Side Drapes 100 x 90 cm: • Positioning Tape • CONTROL* Fabric Reinforcement
	-	79217	1 Bottom Drape 175 x 180 cm: • Positioning Tape • CONTROL* Fabric Reinforcement
	-	79218	1 Top Drape 149 x 241 cm: • Positioning Tape • CONTROL* Fabric Reinforcement


UNIVERSAL PACK, LARGE

CODE	DESCRIPTION		ITEMS/CASE
88765	Universal Pack, Large, excl. Mayo Stand Cover		10
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 96 x 96 cm
42222		42428	1 Back Table Cover, Reinforced 152 x 191 cm
89053		79053	1 Tape Strip 10 x 55 cm
89701		79720	1 Absorbent Towel 38 x 56 cm
	-	-	1 Suture Bag
	-	79216	2 Side Drapes 100 x 90 cm: • Positioning Tape • CONTROL* Fabric Reinforcement
	-	79217	1 Bottom Drape 175 x 180 cm: • Positioning Tape • CONTROL* Fabric Reinforcement
	-	79218	1 Top Drape 149 x 241 cm: • Positioning Tape • CONTROL* Fabric Reinforcement

All codes in **bold** indicate a sterile product.

UNIVERSAL PACK PLUS

CODE	DESCRIPTION		ITEMS/CASE
88760	Universal Pack Plus, incl. Mayo Stand Cover		8
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
42222	42428	-	1 Back Table Cover, Reinforced 152 x 191 cm
89603	79603	-	1 Mayo Stand Cover 80 x 145 cm
89053	79053	-	1 Tape Strip 10 x 55 cm
89701	79720	-	1 Absorbent Towel 38 x 56 cm
-	-	-	1 Suture Bag
-	-	76365	2 Side Drapes 91 x 112 cm: • Positioning Tape • CONTROL [®] Plus Fabric Reinforcement
-	-	78961	1 Bottom Drape 193 x 239 cm: • Positioning Tape • CONTROL [®] Plus Fabric Reinforcement • Tube Holders
-	-	76363	1 Bar Drape 157 x 269 cm: • Armboard Covers • Positioning Tape • CONTROL [®] Plus Fabric Reinforcement • Tube Holders




UNIVERSAL PACK PLUS

CODE	DESCRIPTION		ITEMS/CASE
88767	Universal Pack Plus, excl. Mayo Stand Cover		9
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
42222	42428	-	1 Back Table Cover, Reinforced 152 x 191 cm
89053	79053	-	1 Tape Strip 10 x 55 cm
89701	79720	-	1 Absorbent Towel 38 x 56 cm
-	-	-	1 Suture Bag
-	-	76365	2 Side Drapes 91 x 112 cm: • Positioning Tape • CONTROL [®] Plus Fabric Reinforcement
-	-	78961	1 Bottom Drape 193 x 239 cm: • Positioning Tape • CONTROL [®] Plus Fabric Reinforcement • Tube Holders
-	-	76363	1 Bar Drape 157 x 269 cm: • Armboard Covers • Positioning Tape • CONTROL [®] Plus Fabric Reinforcement • Tube Holders

All codes in **bold** indicate a sterile product.

GENERAL SURGERY PACKS



UNIVERSAL PACK

CODE	DESCRIPTION		ITEMS/CASE
88100	Universal Pack, incl. Mayo Stand Cover		12
	STERILE	NON-STERILE	CONTENTS
42224		42451	1 Back Table Cover, Reinforced 152 x 229 cm
89603		79603	1 Mayo Stand Cover 80 x 145 cm
89053		79053	1 Tape Strip 10 x 55 cm
89701		79720	2 Absorbent Towels 38 x 56 cm
-		79012	2 Side Drapes 75 x 100 cm: • Positioning Tape • CONTROL[®] Plus Fabric Reinforcement
-		79011	1 Bottom Drape 173 x 180 cm: • Positioning Tape • CONTROL[®] Plus Fabric Reinforcement
-		79010	1 Top Drape 148 x 240 cm: • Positioning Tape • CONTROL[®] Plus Fabric Reinforcement

UNIVERSAL PACK


CODE	DESCRIPTION		ITEMS/CASE
88103	Universal Pack, excl. Mayo Stand Cover		12
	STERILE	NON-STERILE	CONTENTS
-		-	1 Sterilisation Wrap 96 x 96 cm
42224		42451	1 Back Table Cover, Reinforced 152 x 229 cm
89053		79053	1 Tape Strip 10 x 55 cm
89701		79720	2 Absorbent Towels 38 x 56 cm
-		79012	2 Side Drapes 75 x 100 cm: • Positioning Tape • CONTROL[®] Plus Fabric Reinforcement
-		79011	1 Bottom Drape 173 x 180 cm: • Positioning Tape • CONTROL[®] Plus Fabric Reinforcement
-		79010	1 Top Drape 148 x 240 cm: • Positioning Tape • CONTROL[®] Plus Fabric Reinforcement

All codes in **bold** indicate a sterile product.

GENERAL SURGERY PACKS


UNIVERSAL PACK IV

CODE	DESCRIPTION		ITEMS/CASE
7772	Universal Pack IV, incl. Mayo Stand Cover		8
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
-	-	42451	1 Back Table Cover, Reinforced 152 x 229 cm
89603	-	79603	1 Mayo Stand Cover 80 x 145 cm
89701	-	79720	1 Absorbent Towel 38 x 56 cm
-	-	-	1 Suture Bag
-	-	76365	2 Side Drapes 91 x 112 cm: • Positioning Tape • CONTROL* Plus Fabric Reinforcement
-	-	76366	1 Bottom Drape 157 x 178 cm: • Positioning Tape • CONTROL* Plus Fabric Reinforcement • Tube Holders
-	-	76367	1 Top Drape 157 x 241 cm: • Positioning Tape • CONTROL* Plus Fabric Reinforcement



All codes in **bold** indicate a sterile product.

GENERAL SURGERY PACKS



HALYARD BASICS* UNIVERSAL PACK I

CODE	DESCRIPTION		ITEMS/CASE
88900	HALYARD BASICS* Universal Pack I, incl. Mayo Stand Cover		5
	STERILE	NON-STERILE	CONTENTS
88673		78669	1 Back Table Cover, Zone Reinforced 152 x 191 cm
47348		-	1 Mayo Stand Cover, Zone Reinforced 80 x 144 cm
-		79790	1 Top Drape 275 x 150 cm
-		79787	1 Bottom Drape 195 x 195 cm
-		-	2 Side Drapes 90 x 100 cm
-		-	1 Tape Strip 10 x 55 cm
-		-	2 Absorbent Towels 33 x 39 cm

HALYARD BASICS* UNIVERSAL PACK II

CODE	DESCRIPTION		ITEMS/CASE
88901	HALYARD BASICS* Universal Pack II, excl. Mayo Stand Cover		5
	STERILE	NON-STERILE	CONTENTS
88673		78669	1 Back Table Cover, Zone Reinforced 152 x 191 cm
-		79790	1 Top Drape 275 x 150 cm
-		79787	1 Bottom Drape 195 x 195 cm
-		-	2 Side Drapes 90 x 100 cm
-		-	1 Tape Strip 10 x 55 cm
-		-	2 Absorbent Towels 33 x 39 cm

UNIVERSAL PACK

CODE	DESCRIPTION		ITEMS/CASE
K887603	Universal Pack, incl. Mayo Stand Cover		54
	STERILE	NON-STERILE	CONTENTS
-		-	1 Mayo Stand Cover, 80 x 145 cm
-		-	1 Back Table Cover, 152 x 190 cm
-		-	2 Side Drapes 91 x 115 cm • Positioning Tape • CONTROL* Fabric Reinforcement 40.5 x 61 cm
-		-	2 Surgical Drapes 193 x 239 cm • Positioning Tape • CONTROL* Fabric Reinforcement 61 x 81 cm

All codes in **bold** indicate a sterile product.

UNIVERSAL PACK

CODE	DESCRIPTION	ITEMS/CASE
K887621	Universal Pack, Bi-laminated Fabric, Sterile	6
SIZE	DESCRIPTION	
152 x 254 cm	Top Drape -Positioning Tape	
90 x 90 cm	Outer Wrap	
152 x 191 cm	Back Table Cover	
76 x 89 cm	2 Side Drape -Positioning Tape	
196 x 196 cm	Bottom Drape -Positioning Tape	
10 x 51 cm	ISO-SEAL___Drape Tape	
43 x 51 cm	Absorbent Towel	

UNIVERSAL PACK WITH MAYO STAND COVER

CODE	DESCRIPTION	ITEMS/CASE
K887625	Universal Pack, incl. Mayo Stand Cover, Bi-laminated, Sterile	6
SIZE	DESCRIPTION	
152 x 254 cm	Top Drape -Positioning Tape	
75 x 127 cm	Mayo Stand Cover	
90 x 90 cm	Outer Wrap	
152 x 191 cm	Back Table Cover	
76 x 89 cm	2 Side Drape -Positioning Tape	
196 x 196 cm	Bottom Drape -Positioning Tape	
10 x 51 cm	ISO-SEAL Drape Tape	
43 x 51 cm	Absorbent Towel	


UNIVERSAL PACK WITH LARGE TOP DRAPE

CODE	DESCRIPTION	ITEMS/CASE
K887623	Universal Pack, Top Drape L, Bi-laminated Fabric, Sterile	6
SIZE	DESCRIPTION	
176 x 300 cm	Top Drape -Positioning Tape	
90 x 90 cm	Outer Wrap	
152 x 191 cm	Back Table Cover	
76 x 89 cm	2 Side Drape -Positioning Tape	
196 x 196 cm	Bottom Drape -Positioning Tape	
10 x 51 cm	ISO-SEAL Drape Tape	
43 x 51 cm	Absorbent Towel	



All codes in **bold** indicate a sterile product.

GENERAL SURGERY PACKS



UNIVERSAL PACK WITH MAYO STAND COVER AND LARGE TOP DRAPE

CODE	DESCRIPTION	ITEMS/CASE
K887626	Universal Pack, incl. Mayo Stand Cover, Top Drape L, Bi-laminated Fabric, Sterile	6
SIZE	DESCRIPTION	
175 x 300 cm	Top Drape -Positioning Tape	
80 x 127 cm	Mayo Stand Cover	
90 x 90 cm	Outer Wrap	
152 x 191 cm	Back Table Cover	
76 x 89 cm	2 Side Drape -Positioning Tape	
196 x 196 cm	Bottom Drape -Positioning Tape	
10 x 51 cm	ISO-SEAL Drape Tape	
43 x 51 cm	Absorbent Towel	

BASIC PACK VIII

CODE	DESCRIPTION	ITEMS/CASE
88113	Basic Pack VIII	24
STERILE	NON-STERILE	CONTENTS
42224	42451	1 Back Table Cover, Reinforced 152 x 229 cm
89603	79603	1 Mayo Stand Cover 80 x 145 cm
-	-	1 Suture Bag

U-BAR* PACK VI

CODE	DESCRIPTION	ITEMS/CASE
77794	U-BAR* Pack VI	9
STERILE	NON-STERILE	CONTENTS
-	-	1 Sterilisation Wrap 96 x 96 cm
-	42451	1 Back Table Cover, Reinforced 152 x 229 cm
89603	79603	1 Mayo Stand Cover 80 x 145 cm
89701	79720	1 Absorbent Towel 38 x 56 cm
-	-	1 Suture Bag
-	76361	1 U-Drape 193 x 305 cm: • Split 10 x 101 cm • Positioning Tape • CONTROL* Plus Fabric Reinforcement • Tube Holders
-	76363	1 Bar Drape 157 x 269 cm: • Armboard Covers • Positioning Tape • CONTROL* Plus Fabric Reinforcement • Tube Holders


All codes in **bold** indicate a sterile product.

U-BAR* PACK

CODE	DESCRIPTION		ITEMS/CASE
77853	U-BAR* Pack		10
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
-	-	42451	1 Back Table Cover, Reinforced 152 x 229 cm
89701	-	79720	1 Absorbent Towel 38 x 56 cm
-	-	-	1 Suture Bag
-	-	76361	1 U-Drape 193 x 305 cm: • Split 10 x 101 cm • Positioning Tape • CONTROL* Plus Fabric Reinforcement • Tube Holders
-	-	76363	1 Bar Drape 157 x 269 cm: • Armboard Covers • Positioning Tape • CONTROL* Plus Fabric Reinforcement • Tube Holders
89053	-	79053	1 Tape Strip 10 x 55 cm


U-BAR* PACK VII

CODE	DESCRIPTION		ITEMS/CASE
77769	U-BAR* Pack VII		10
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
-	-	42451	1 Back Table Cover, Reinforced 152 x 229 cm
89603	-	79603	1 Mayo Stand Cover 80 x 145 cm
89701	-	79720	1 Absorbent Towel 38 x 56 cm
-	-	-	1 Suture Bag
-	-	76364	1 Breast/Chest U-Drape 193 x 305 cm: • Split 30 x 76 cm • Positioning Tape • CONTROL* Plus Fabric Reinforcement • Tube Holders
-	-	76363	1 Bar Drape 157 x 269 cm: • Armboard Covers • Positioning Tape • CONTROL* Plus Fabric Reinforcement • Tube Holders



All codes in **bold** indicate a sterile product.

GENERAL SURGERY PACKS



HALYARD BASICS* U-BAR* PACK I

CODE	DESCRIPTION		ITEMS/CASE
88902	HALYARD BASICS* U-BAR* Pack I		10
	STERILE	NON-STERILE	CONTENTS
88673		78669	1 Back Table Cover, Zone Reinforced 152 x 191 cm
47348		-	1 Mayo Stand Cover, Zone Reinforced 80 x 144 cm
-		79790	1 Top Drape 275 x 150 cm
-		79532	1 U-Drape 203 x 315 cm: • Split 10 x 112 cm
-		-	1 Tape Strip 10 x 55 cm
-		-	2 Absorbent Towels 33 x 39 cm


HALYARD BASICS* U-BAR* PACK III

CODE	DESCRIPTION		ITEMS/CASE
88904	HALYARD BASICS* U-BAR* Pack III		10
	STERILE	NON-STERILE	CONTENTS
88673		78669	1 Back Table Cover, Zone Reinforced 152 x 191 cm
47348		-	1 Mayo Stand Cover, Zone Reinforced 80 x 144 cm
-		79790	1 Top Drape 275 x 150 cm
-		-	1 Breast and Chest U-Drape 193 x 305 cm: • Split 30 x 76 cm
-		-	1 Tape Strip 10 x 55 cm
-		-	2 Absorbent Towels 33 x 39 cm

All codes in **bold** indicate a sterile product.


VEIN STRIPPING PACK

CODE	DESCRIPTION		ITEMS/CASE
59709	Vein Stripping Pack		5
	STERILE	NON-STERILE	CONTENTS
	42222	42428	1 Back Table Cover, Reinforced 152 x 191 cm
	89603	79603	1 Mayo Stand Cover 80 x 145 cm
	89701	79720	4 Absorbent Towels 38 x 56 cm
	-	-	1 Suture Bag
	47622	47623	1 Large Drape 193 x 139 cm
	-	76356	1 Perineal Drape 187 x 241 cm: • Split 44 x 95 cm • Positioning Tape • CONTROL* Plus Fabric Reinforcement
	76291	76492	1 Drape with Tape 157 x 241 cm



All codes in **bold** indicate a sterile product.

GENERAL SURGERY PACKS



SMALL PROCEDURE PACK II

CODE	DESCRIPTION		ITEMS/CASE
7717	Small Procedure Pack II		10
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
42222	-	42428	1 Back Table Cover, Reinforced 152 x 191 cm
-	-	-	1 Fenestrated Drape 211 x 241 cm: • Circular Fenestration 12 cm, centered • CONTROL[®] Plus Fabric Reinforcement
89701	-	79720	2 Absorbent Towels 38 x 56 cm
89053	-	79053	1 Tape Strip 10 x 55 cm

HALYARD BASICS* SMALL PROCEDURE PACK II

CODE	DESCRIPTION		ITEMS/CASE
89686	HALYARD BASICS* Small Procedure Pack II		10
	STERILE	NON-STERILE	CONTENTS
88673	-	78669	1 Back Table Cover, Zone Reinforced 152 x 191 cm
-	-	79745	1 Fenestrated Drape 211 x 241 cm: • Circular Fenestration 12 cm, centered • Absorbent Reinforcement
-	-	-	1 Tape Strip 10 x 55 cm
-	-	-	2 Absorbent Towels 33 x 39 cm

LAPAROSCOPY PACK

CODE	DESCRIPTION		ITEMS/CASE
76216	Laparoscopy pack		6
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 63 x 63 cm
42224	-	42451	1 Back Table Cover, Reinforced 152 x 229 cm
89603	-	79603	1 Mayo Stand Cover 80 x 145 cm
90013	-	79973	3 Surgical Gowns, Large, X-long
89701	-	79720	7 Absorbent Towels 38 x 56 cm
-	-	-	2 Iso-Seal Tape Strips 12 x 59 cm
-	-	76374	1 Laparoscopy Drape with Attached Leggings 254 x 177 x 279 cm • 2 Pockets 27 x 71 cm • Abdominal Fenestration 30 x 33 cm • Positioning Tape



All codes in **bold** indicate a sterile product.

UNIVERSAL LAPAROSCOPY / PELVISCOPY PACK WITH SIDE POCKETS


CODE	DESCRIPTION		ITEMS/CASE
76263	Universal Laparoscopy / Pelviscopy Pack		8
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
42222	42428	-	1 Back Table Cover, Reinforced 152 x 191 cm
-	-	76377	1 Laparoscopy Drape with Attached Leggings 254 x 279 x 178 cm and side pouches 28 x 58 cm: • Abdominal Fenestration (30 x 33 cm) with 5 cm Tape Width • Perineal Fenestration (22 x 10 cm) with 2.5 cm Tape Width • CONTROL* Plus Fabric Reinforcement • With Flap over both Fenestrations
-	-	-	2 Iso-Seal Tape Strips 12 x 59 cm

LAPAROSCOPY PELVISCOPY PACK V

CODE	DESCRIPTION		ITEMS/CASE
77792	Laparoscopy Pelviscopy Pack V		6
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
42222	42428	-	1 Back Table Cover, Reinforced 152 x 191 cm
89603	79603	-	1 Mayo Stand Cover 80 x 145 cm
-	-	-	1 Suture Bag
-	-	76380	1 Laparoscopy Drape 254 x 312 x 74 cm, with Attached Leggings: • 23 x 28 cm Abdominal Fenestration • 7.6 x 18 cm Trapezoidal Perineal Fenestration with Flap • Positioning Tape on both Fenestrations • Tube Holders • CONTROL* Fabric Reinforcement



GENERAL SURGERY DRAPE



SMALL SURGICAL DRAPE PACK

CODE	DESCRIPTION	ITEMS/CASE
76243	Small Surgical Drape Pack	25
STERILE	NON-STERILE	CONTENTS
-	-	1 Sterilisation Wrap 89 x 89 cm
-	42446	1 Back Table Cover, Reinforced 112 x 112 cm
89363	79293	4 Utility Drapes 90 x 90 cm

FENESTRATED DRAPE

MINOR PROCEDURE FENESTRATED DRAPE, 76 X 76 CM

- Fenestration 3 x 5 cm
- CONTROL* Fabric Reinforcement 12 x 20 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE	ITEMS/DISPENSER	DISPENSERS/CASE
89206	Minor Procedure Fenestrated Drape, Sterile	100	25 x 4	4
79206	Minor Procedure Fenestrated Drape Non-Sterile	250	-	-

MINOR PROCEDURE FENESTRATED DRAPE, 122 X 152 CM

- Fenestration (Diameter 7 cm)
- CONTROL* Fabric Reinforcement 31 x 31 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89207	Minor Procedure Fenestrated Drape, Sterile	40
79207	Minor Procedure Fenestrated Drape, Non-Sterile	60

MINOR PROCEDURE FENESTRATED DRAPE, 122 X 152 CM

- Fenestration 12 x 15 cm
- CONTROL* Fabric Reinforcement 41 x 31 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89208	Minor Procedure Fenestrated Drape, Sterile	40
79208	Minor Procedure Fenestrated Drape, Non-Sterile	60

All codes in **bold** indicate a sterile product.

GENERAL SURGERY DRAPES

MINOR PROCEDURE FENESTRATED DRAPE, 193 X 315 CM

- Fenestration 15 x 15 cm
- CONTROL* Fabric Reinforcement 51 x 67 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89209	Minor Procedure Fenestrated Drape, Sterile	19
79209	Minor Procedure Fenestrated Drape, Non-Sterile	23

MINOR PROCEDURE FENESTRATED DRAPE, 119 X 152 CM

- SMS
- Fenestration 12 x 15 cm
- CONTROL* Plus Fabric Reinforcement 30 x 40 cm

CODE	DESCRIPTION	ITEMS/CASE
K792610	Minor Procedure Fenestrated Drape, Non-Sterile	60

HALYARD BASICS* EPIDURAL DRAPE WITH TAPE, 50 X 75 CM

- Oval Fenestration 8 x 10 cm, centered
- Spunbond Film
- Lateral Tape 50 x 2.5 cm
- Tape around Fenestration
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
14364	HALYARD BASICS* Epidural Drape with Tape, Sterile	40
14365	HALYARD BASICS* Epidural Drape with Tape, Non-Sterile	160

HALYARD BASICS* FENESTRATED DRAPE, 75 X 75 CM


- Oval Fenestration 7 x 10 cm, centered
- Spunbond Film

CODE	DESCRIPTION	ITEMS/CASE
89638	HALYARD BASICS* Fenestrated Drape, Sterile	35

HALYARD BASICS* FENESTRATED DRAPE, 90 X 90 CM

- Circular Fenestration 10 cm, centered
- Spunbond Film

CODE	DESCRIPTION	ITEMS/CASE
89670	HALYARD BASICS* Fenestrated Drape, Sterile	30
79692	HALYARD BASICS* Fenestrated Drape, Non-Sterile	120



All codes in **bold** indicate a sterile product.

GENERAL SURGERY DRAPE

SMALL FENESTRATED DRAPE

- 40GSM Bi-Laminated
- Circular Fenestration 5cm

CODE	SIZE	ITEMS/CASE
K891027	50x60cm	400
K791027	70x60cm	470

HALYARD BASICS* FENESTRATED DRAPE WITH TAPE, 90 X 90 CM

- Circular Fenestration 10 cm, centered
- Spunbond Film
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89671	HALYARD BASICS* Fenestrated Drape with Tape, Sterile	30
79693	HALYARD BASICS* Fenestrated Drape with Tape, Non-Sterile	50

HALYARD BASICS* VARIABLE FENESTRATED DRAPE WITH TAPE, 75 X 100 CM


- Variable Fenestration: half circle on each drape
- Spunbond Film
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89680	HALYARD BASICS* Variable Fenestrated Drape with Tape, Sterile	30
79784	HALYARD BASICS* Variable Fenestrated Drape with Tape, Non-Sterile	120

HALYARD BASICS* FENESTRATED DRAPE, 211 X 241 CM

- Circular Fenestration 12 cm, centered, with Positioning Tape
- SMS
- Absorbent Reinforcement

CODE	DESCRIPTION	ITEMS/CASE
79745	HALYARD BASICS* Fenestrated Drape, Non-Sterile	30



All codes in **bold** indicate a sterile product.

VARIABLE FENESTRATED ABSORBENT DRAPE

- 40GSM Bi-Laminated
- Fenestration 18x2 4x6cm
- Positioning tape

CODE	SIZE	ITEMS/CASE
K892031	45x75cm per half	140
K792031	45x75cm per half	1700

FENESTRATED DRAPE, 68 X 68 CM

- SMS
- Circular Fenestration 18 cm, centered
- CONTROL* Fabric Reinforcement 30 x 30 cm

CODE	DESCRIPTION	ITEMS/CASE
K792015	Fenestrated Drape, Non-Sterile	180

ABSORBENT FENESTRATED DRAPE




- 40GSM Bi-Laminated
- Ovale Fenestration 6.5x8.25cm
- Positioning tape 13cm

CODE	SIZE	ITEMS/CASE
K892036	120x150cm	180
K792036	120x150cm	235
K792035	50x60cm	450
K892034	75x75cm	230
K792034	75x75cm	230

ABSORBENT FENESTRATED DRAPE

- 40GSM Bi-Laminated
- Circular Fenestration 10cm
- Positioning tape 15cm

CODE	SIZE	ITEMS/CASE
K892033	75x90cm	180
K792033	75x90cm	230



GENERAL SURGERY DRAPE

ABSORBENT FENESTRATED DRAPE

- 40GSM Bi-Laminated
- Circular Fenestration 5cm
- Positioning tape 10cm

CODE	SIZE	ITEMS/CASE
K892032	50x60cm	400
K792032	50x60cm	450

U-DRAPE

U-DRAPE, 193 X 305 CM

- Split 10 x 102 cm
- CONTROL* or CONTROL* Plus Fabric Reinforcement 122 x 61 cm
- Tube Holders
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
76361	U-Drape, CONTROL* Plus, Non-Sterile	15
13231	U-Drape, CONTROL*, Non-Sterile - handi-bin	190

U-DRAPE, 193 X 305 CM


- Split 10 x 102 cm
- CONTROL* Plus Fabric Reinforcement 122 x 61 cm
- Tube Holders
- Positioning Tape, width: 2 cm

CODE	DESCRIPTION	ITEMS/CASE
89301	U-Drape, CONTROL* Plus, Sterile	17

U-DRAPE, 200 X 150 CM

- Split 8 x 64 cm
- SMS
- CONTROL* Fabric Reinforcement

CODE	DESCRIPTION	ITEMS/CASE
K790790	U-Drape, Non-Sterile	22



All codes in **bold** indicate a sterile product.

LARGE ORTHO U-DRAPE, 224 X 318 CM

- Split 15 x 137 cm
- CONTROL* or CONTROL* Plus Fabric Reinforcement 81 x 170 cm
- Tube Holders
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89311	Large Ortho U-Drape, CONTROL* Plus, Sterile	11
79363	Large Ortho U-Drape, CONTROL*, Non-Sterile	15
K793630	Large Ortho U-Drape, CONTROL*, Non-Sterile	22
76368	Large Ortho U-Drape, CONTROL* Plus, Non-Sterile	15

HALYARD BASICS* U-DRAPE, 203 X 315 CM

- Split 10 x 112 cm
- SMS
- Absorbent Reinforcement

CODE	DESCRIPTION	ITEMS/CASE
79532	HALYARD BASICS* U-Drape, Non-Sterile	30

SPLIT DRAPE 228 X 260 CM W/ U 20 X 100 CM





- Main Drape Size : 228 x 260 cm
- U-Split Size : 20 x 100 cm
- Reinforcement Size : 81.28 x 121.92 cm
- SMS

CODE	DESCRIPTION	ITEMS/CASE
K793214	Split Drape, 228 X 260 cm. W/split 20 X 100 cm, Non-Sterile	40

TRI LAMINATED FABRIC SPLIT DRAPE, 228 X 260 CM U 20 X 100 CM

- Main Drape Size : 228 x 260 cm
- U-Split Size : 20 x 100 cm
- Tri-Laminated Fabric

CODE	DESCRIPTION	ITEMS/CASE
K793215	Split Drape, 228 X 260 cm. W/split 20 X 100 cm, Non-Sterile	40



GENERAL SURGERY DRAPES

HYBRID SPLIT SHEET, 250 CM X 430 CM

CODE	DESCRIPTION	ITEMS/CASE
K893330	Hybrid Split Sheet, Sterile	10
K793330	Hybrid Split Sheet, Non-Sterile	12

BI LAMINATED FABRIC SPLIT DRAPE, 228 X 256 CM U 20 X 100 CM

- Main Drape Size : 228 x 256 cm
- U-Split Size : 20 cm x 100 cm
- Base material: Bi-Laminated fabric

CODE	DESCRIPTION	ITEMS/CASE
K793216	Split Drape, 228 X 256 cm. W/split 20 X 100 cm, Non-Sterile	40

TRI LAMINATED FABRIC 152 X 183 CM. U 7 X 70 CM

- Main Drape Size : 152 x 183 cm.
- U-Split Size : 7 x 70 cm.
- Reinforcement Size : 81.28 x 121.92 cm
- Base material : SMS

CODE	DESCRIPTION	ITEMS/CASE
K793219	Split Drape, 152 X 183 cm. W/split 7 x 70 cm, Non-Sterile	40

SMS 152 X 183 CM. U 7 X 70 CM

- Tri-Laminated Drape Size : 152 x 183 cm
- U-Split Size surrounded by medical tape : 7 x 70 cm
- Tri-Laminated Blue Fabric

CODE	DESCRIPTION	ITEMS/CASE
K893217	Split Drape, Sterile	40

All codes in **bold** indicate a sterile product.

BAR DRAPES

ORTHO BAR DRAPE, 254 X 152 CM

- CONTROL* Plus Fabric Reinforcement 41 x 61 cm
- Tube Holders
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89381	Ortho Bar Drape, Sterile	18
79360	Ortho Bar Drape, Non-Sterile	28

BAR DRAPE WITH ARMBOARD COVERS, 269 X 157 CM

- Armboard Covers
- CONTROL* or CONTROL* Plus Fabric Reinforcement 41 x 61 cm
- Tube Holders
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79375	Bar Drape, CONTROL*, Non-Sterile	28
76363	Bar Drape, CONTROL* Plus, Non-Sterile	26

SURGICAL DRAPE, 195 X 241 CM





- SMS
- CONTROL* Fabric Reinforcement
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
K791412	Surgical Drape, Non-Sterile	50


BAR DRAPE, 183 X 300 CM

- SMS
- CONTROL* Fabric Reinforcement 20 x 150 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
K793703	Bar Drape, Non-Sterile	40



GENERAL SURGERY DRAPES



BOTTOM DRAPES

HALYARD BASICS* BOTTOM DRAPE, 195 X 195 CM

- SMS

CODE	DESCRIPTION	ITEMS/CASE
79787	HALYARD BASICS* Bottom Drape, Non-Sterile	35

BOTTOM DRAPE, 180 X 149 CM

- CONTROL* Fabric Reinforcement 30 x 89 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79214	Bottom Drape, Non-Sterile	40

BOTTOM DRAPE, 180 X 175 CM

- CONTROL* Fabric Reinforcement 30 x 89 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79217	Bottom Drape, Non-Sterile	35

BOTTOM DRAPE, 157 X 178 CM

- CONTROL* Plus Fabric Reinforcement 81 x 61 cm
- Tube Holders
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
76366	Bottom Drape, Non-Sterile	22

All codes in **bold** indicate a sterile product.

GENERAL SURGERY DRAPES

BOTTOM DRAPE, 173 X 180 CM

- CONTROL* Plus Fabric Reinforcement 63 x 61 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79011	Bottom Drape, Non-Sterile	50

BOTTOM DRAPE, 183 X 230 CM

- CONTROL* Fabric Reinforcement 10.6 x 50.8 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
K795610	Bottom Drape, Non-Sterile	50
K795616	Bottom Drape, Non-Sterile	60

BOTTOM DRAPE WITH ADHESIVE, 220 X 175 CM




- Main Drape size : 220 x 175 cm
- Reinforcement size 30.5 x 101.6 cm
- Adhesive size : 2.5 x 101.6 cm.
- Base material : SMS fabric
- Reinforcement : Control Plus

CODE	DESCRIPTION	ITEMS/CASE
K795613	Bottom Drape W/Adhesive, Non-Sterile	50
K795614	Bottom Drape W/Adhesive, Non-Sterile	50

BOTTOM DRAPE, 196 X 196 CM

- Bi-Laminated Blue Fabric

CODE	DESCRIPTION	ITEMS/CASE
K795615	Bottom Drape, Non-Sterile	50



All codes in **bold** indicate a sterile product.

GENERAL SURGERY DRAPES

SIDE DRAPES

SIDE DRAPE, 90 X 100 CM

- CONTROL* Fabric Reinforcement 99 x 31 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79216	Side Drape, Non-Sterile	70

SIDE DRAPE, 75 X 100 CM

- CONTROL* Plus Fabric Reinforcement 41 x 51 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79012	Side Drape, Non-Sterile	150

SIDE DRAPE, 76 X 91 CM

- SMS
- CONTROL* Fabric Reinforcement
- Positioning Tape 2.5 cm

CODE	DESCRIPTION	ITEMS/CASE
K792130	Side Drape, Non-Sterile	200

SIDE DRAPE, 91 X 112 CM

- CONTROL* Plus Fabric Reinforcement 46 x 112 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
76365	Side Drape, Non-Sterile	54

SIDE DRAPE WITH ADHESIVE, 75 X 89 CM

- Base material : Bi-Laminated

CODE	DESCRIPTION	ITEMS/CASE
K795806	Side Drape, Non-Sterile	50

SIDE DRAPE WITH ADHESIVE, 120 X 100 CM

- SMS
- adhesive size 2,5cm x 120 cm

CODE	DESCRIPTION	ITEMS/CASE
K795805	Side Drape with Adhesive	70

SIDE DRAPE 100 CM X130 CM

CODE	DESCRIPTION	ITEMS/CASE
K795809	Side Drape, Non-Sterile	100

All codes in **bold** indicate a sterile product.

TOP DRAPES

TOP DRAPE, 149 X 241 CM

- CONTROL* Fabric Reinforcement
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79218	Top Drape, Non-Sterile	35

TOP DRAPE, 149 X 180 CM

- CONTROL* Fabric Reinforcement 30 x 89 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79215	Top Drape, Non-Sterile	4

TOP DRAPE, 148 X 240 CM

- CONTROL* Plus Fabric Reinforcement 42 x 61 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79010	Top Drape, Non-Sterile	50

TOP DRAPE 152 CM X 254 CM

CODE	DESCRIPTION	ITEMS/CASE
K795713	Top Drape, Non-Sterile	34

TOP DRAPE, 157 X 241 CM

- CONTROL* Plus Fabric Reinforcement 30 x 61 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
76367	Top Drape, Non-Sterile	33


TOP DRAPE 175 CM X 300 CM

CODE	DESCRIPTION	ITEMS/CASE
K795714	Top Drape, Non-Sterile	34

HALYARD BASICS* TOP DRAPE, 275 X 150 CM

- SMS

CODE	DESCRIPTION	ITEMS/CASE
79790	HALYARD BASICS* Top Drape, Non-Sterile	35



All codes in **bold** indicate a sterile product.

GENERAL SURGERY DRAPE



TOP DRAPE WITH ADHESIVE, 300 X 175 CM

- Main Drape size : 300 x 175 cm
- Reinforcement size : 30.5 x 101.6 cm
- Base material: SMS fabric
- Reinforcement : Control Plus

CODE	DESCRIPTION	ITEMS/CASE
K795711	Top Drape W/Adhesive, Non-Sterile	40

ABSORBENT UNDERLEG DRAPE

POLY DRAPE, 158 X 235 CM

- CONTROL* Plus Fabric Reinforcement 157 x 104 cm

CODE	DESCRIPTION	ITEMS/CASE
59705	Poly Drape, Sterile	14

SINGLE DRAPE

SMALL ABSORBENT DRAPE

- 40 GSM Bi-laminated

CODE	SIZE	ITEMS/CASE
K891022	45x50cm	400
K791022	45x50cm	470
K891025	75x90cm	440
K791025	75x90cm	450
K791026	112x145cm	235
K891023	45x75cm	400
K791023	45x75cm	460
K891024	75x75cm	400
K791024	75x75cm	450

DRAPE, 90 X 90 CM

- SMS

CODE	DESCRIPTION	ITEMS/CASE
89353	Drape, Sterile	120
76323	Drape, Non-Sterile	200

HALYARD BASICS* DRAPE, 90 X 90 CM

- Spunbond Film

CODE	DESCRIPTION	ITEMS/CASE
89678	HALYARD BASICS* Drape, Sterile	30
79782	HALYARD BASICS* Drape, Non-Sterile	90

All codes in **bold** indicate a sterile product.

GENERAL SURGERY DRAPES

HALF DRAPE, 147 X 101 CM

- SMS

CODE	DESCRIPTION	ITEMS/CASE
47614	Half Drape, Sterile	60
47616	Half Drape, Non-Sterile	80

HALYARD BASICS* HALF DRAPE, 112 X 150 CM

- Spunbond Film

CODE	DESCRIPTION	ITEMS/CASE
89672	HALYARD BASICS* Half Drape, Sterile	40
79694	HALYARD BASICS* Half Drape, Non-Sterile	50

LARGE DRAPE, 139 X 193 CM

- SMS

CODE	DESCRIPTION	ITEMS/CASE
47622	Large Drape, Sterile	32
47623	Large Drape, Non-Sterile	32

HALYARD BASICS* LARGE DRAPE, 150 X 193 CM

- Spunbond Film

CODE	DESCRIPTION	ITEMS/CASE
89674	HALYARD BASICS* Large Drape, Sterile	30
79696	HALYARD BASICS* Large Drape, Non-Sterile	40

MEDIUM DRAPE, 185 X 101 CM






- SMS

CODE	DESCRIPTION	ITEMS/CASE
47618	Medium Drape, Non-Sterile	64

X-LARGE DRAPE, 193 X 243 CM


- SMS

CODE	DESCRIPTION	ITEMS/CASE
47631	X-Large Drape, Non-Sterile	24



All codes in **bold** indicate a sterile product.

GENERAL SURGERY DRAPE



ABSORBENT DRAPE, 100 X 152 CM

- Bi-Laminated Drape Size : 100 x 152 cm

CODE	DESCRIPTION	ITEMS/CASE
K791020	Absorbent Drape, Non-Sterile	110

UTILITY DRAPE

SMALL ABSORBENT DRAPE WITH TAPE

- 40 GSM Bi-laminated
- Full length tape on one side

CODE	SIZE	TAPE LENGTH	ITEMS/CASE
K897322	75x75cm	75cm	200
K797322	75x75cm	75cm	300
K897320	75x90cm	90cm	440
K797320	75x90cm	90cm	440
K897321	45x50cm	45cm	460
K797321	45x50cm	45cm	460
K797316	50x50cm	50 cm	780

UTILITY DRAPE WITH TAPE, 76 X 76 CM, PER PAIR

- SMS

CODE	DESCRIPTION	ITEMS/CASE
89313	Utility Drape with Tape, Sterile	180 (2 x 90)
79313	Utility Drape with Tape, Non-Sterile	200

HALYARD BASICS* UTILITY DRAPE WITH TAPE, 75 X 75 CM

- Spunbond Film

CODE	DESCRIPTION	ITEMS/CASE
89677	HALYARD BASICS* Utility Drape with Tape, Sterile	35
79699	HALYARD BASICS* Utility Drape with Tape, Non-Sterile	60

UTILITY DRAPE WITH TAPE, 90 X 90 CM, PER PAIR

- SMS

CODE	DESCRIPTION	ITEMS/CASE
89293	Utility Drape with Tape, Sterile	160 (2 x 80)
79293	Utility Drape with Tape, Non-Sterile	170

All codes in **bold** indicate a sterile product.

GENERAL SURGERY DRAPES

UTILITY DRAPE WITH TAPE, 90 X 90 CM

- SMS

CODE	DESCRIPTION	ITEMS/CASE
89363	Utility Drape with Tape, Sterile	100
79293	Utility Drape with Tape, Non-Sterile	170

HALYARD BASICS* UTILITY DRAPE WITH TAPE, 90 X 90 CM

- Spunbond Film

CODE	DESCRIPTION	ITEMS/CASE
89679	HALYARD BASICS* Utility Drape with Tape, Sterile	30
79783	HALYARD BASICS* Utility Drape with Tape, Non-Sterile	90

STANDARD ABSORBENT DRAPE, 91 X 91 CM

CODE	DESCRIPTION	ITEMS/CASE
89343	Standard Absorbent Drape, Sterile	50
79343	Standard Absorbent Drape, Non-Sterile	70

UTILITY DRAPE, 157 X 241 CM






- SMS
- Full Length Tape on the Long Side

CODE	DESCRIPTION	ITEMS/CASE
76291	Utility Drape, Sterile	30
76492	Utility Drape, Non-Sterile	32

HALYARD BASICS* UTILITY DRAPE WITH TAPE, 150 X 240 CM

- SMS
- Full Length Tape on the Long Side

CODE	DESCRIPTION	ITEMS/CASE
89676	HALYARD BASICS* Utility Drape with Tape, Sterile	30



All codes in **bold** indicate a sterile product.

GENERAL SURGERY DRAPE

PLAIN DRAPE 150 CM X 200 CM

CODE	DESCRIPTION	ITEMS/CASE
K791021	Plain Drape, Non-Sterile	60

UTILITY DRAPE, 90 X 157 CM

- SMS
- Full Length Tape on the Long Side

CODE	DESCRIPTION	ITEMS/CASE
76293	Utility Drape, Sterile	60
78965	Utility Drape, Non-Sterile	120

UTILITY DRAPE WITH TAPE, 64 X 38 CM

- SMS
- Full Length Tape on the Long Side

CODE	DESCRIPTION	ITEMS/CASE
79731	Utility Drape with Tape, Non-Sterile	550

UTILITY DRAPE WITH TAPE, 75 X 90 CM

CODE	DESCRIPTION	ITEMS/CASE
K797310	Utility Drape with Tape, Non-Sterile	200





STANDARD ABSORBENT DRAPE, 50 X 50CM

- Full Length Tape on One Side

CODE	DESCRIPTION	ITEMS/CASE
K897312	Utility Drape with Tape, Sterile	100

XXL SURGICAL DRAPE 195 CM X 241 CM

CODE	DESCRIPTION	ITEMS/CASE
K791413	XXL Surgical Drape, Non-Sterile	50



All codes in **bold** indicate a sterile product.

FENESTRATED DRAPE WITH TAPE 90X190

- Main Drape size : 190 x 90 cm.
- Fenestration size : 10.2 cm (Diamètre)
- Base material: Tri-Laminated Fabric (Hydrophilic)

CODE	DESCRIPTION	ITEMS/CASE
K796710	Fenestrated Drape with Tape, Non-sterile	50

BI LAMINATED FABRIC DRAPE 75X90 W/ ADHESIVE


- Main Drape size : 76 x 90 cm
- Base material : Bi-Laminated Fabric (Hydrophilic)

CODE	DESCRIPTION	ITEMS/CASE
K891013	Drape with Adhesive, Sterile	100

BI LAMINATED FABRIC DRAPE 75X90 WITHOUT ADHESIVE

- Main Drape size : 76 x 90 cm
- Base material : Bi-Laminated Fabric (Hydrophilic)

CODE	DESCRIPTION	ITEMS/CASE
K891012	Drape Plain, Sterile	150



All codes in **bold** indicate a sterile product.

GENERAL SURGERY DRAPE

LAPAROTOMY DRAPE

LAPAROTOMY DRAPE WITH INCISE, 254 X 305 X 183 CM


- Fenestration 20 x 33 cm
- CONTROL* Plus Fabric Reinforcement 61 x 96 cm
- Tube Holders
- Armboard Covers
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89282	Laparotomy Drape with Incise, Sterile	15
79282	Laparotomy Drape with Incise, Non-Sterile	15

LAPAROTOMY DRAPE WITH ATTACHED LEGGINGS, 254 X 312 X 74 CM

- Abdominal Fenestration 23 x 28 cm
- Trapezoidal Perineal Fenestration with Flap 7.6 x 18 cm
- CONTROL* Plus Fabric Reinforcement 76 x 38 cm
- Hook and Loop Tube Holders
- Positioning Tape on both Fenestrations

CODE	DESCRIPTION	ITEMS/CASE
76380	Laparotomy Drape with Attached Leggings, Non-Sterile	8



All codes in **bold** indicate a sterile product.

LAPAROTOMY DRAPE WITH ATTACHED LEGGINGS, 255 X 317 X 284 CM


- SMS
- Abdominal Fenestration 33 x 37 cm, with Incise
- Perineal Fenestration 11 x 16 cm
- CONTROL* Reinforcement
- Velcro and pouches with 2 compartments

CODE	DESCRIPTION	ITEMS/CASE
K892245	Laparotomy Drape with Attached Leggings, Sterile	12
K792245	Laparotomy Drape with Attached Leggings, Non-Sterile	8

LAPAROTOMY DRAPE WITH ATTACHED LEGGINGS, 255 X 317 X 284 CM

- SMS
- Abdominal Fenestration 33 x 37 cm
- Perineal Fenestration 11 x 16 cm
- CONTROL* Reinforcement
- Velcro and pouches with 2 compartments

CODE	DESCRIPTION	ITEMS/CASE
K892246	Laparotomy Drape with Attached Leggings, Sterile	12
K792246	Laparotomy Drape with Attached Leggings, Non-Sterile	12



All codes in **bold** indicate a sterile product.

GENERAL SURGERY DRAPES

BREAST/CHEST DRAPES

BREAST DRAPE, 254 X 315 X 193 CM



- Fenestration 15 x 15 cm
- CONTROL* Fabric Reinforcement 51 x 71 cm
- 2 Velcros - Aperture is adhesive and it has 2,5 cm wide tapes
- Armboard Covers

CODE	DESCRIPTION	ITEMS/CASE
89236	Breast Drape, Sterile	10
79236	Breast Drape, Non-Sterile	20

BREAST & CHEST U-DRAPE, 193 X 305 CM

- CONTROL* Plus Fabric Reinforcement 102 x 61 cm
- Tube Holders
- Positioning Tape
- 30 x 76 cm Split

CODE	DESCRIPTION	ITEMS/CASE
76364	Breast & Chest U-Drape, Non-Sterile	19




All codes in **bold** indicate a sterile product.

LAPAROSCOPY DRAPE

LAPAROSCOPY DRAPE WITH ATTACHED LEGGINGS, 254 X 279 X 178 CM WITH SIDE POUCHES 28 X 58 CM

- Abdominal Fenestration 30 x 33 cm
- Perineal Fenestration 22 x 10 cm
- With Flap over both Fenestrations
- CONTROL* Plus Reinforcement 96 x 152 cm


CODE	DESCRIPTION	ITEMS/CASE
76377	Laparoscopy Drape with Attached Leggings, Non-Sterile	10



LAPAROSCOPY/PELVISCOPY DRAPE, 264 X 89 X 264 CM

- SMS
- CONTROL* Plus Fabric Reinforcement 38 x 76 cm
- Trapezoidal Abdominal Fenestration 18 x 23 x L 18 cm
- Perineal Fenestration 7 x 15 cm
- With Flap over both Fenestrations
- Positioning Tape around Abdominal Fenestration


CODE	DESCRIPTION	ITEMS/CASE
79222	Laparoscopy/Pelviscopy Drape, Non-Sterile	22



LAPAROSCOPY/CHOLECYSTECTOMY DRAPE, 295 X 193 X 304 CM


- SMS
- CONTROL* Plus Fabric Reinforcement 61 x 71 cm
- Fenestration 30 x 33 cm
- Instrument Pouch 48 x 33 cm
- Armboard Covers
- Tube Holders
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79223	Laparoscopy/Cholecystectomy Drape, Non-Sterile	13



ORTHOPAEDICS

ORTHOPAEDIC PACKS



UPPER EXTREMITY PACK

CODE	DESCRIPTION		ITEMS/CASE
77781	Upper Extremity Pack		9
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 89 x 89 cm
42222		42428	1 Back Table Cover, Reinforced 152 x 191 cm
89053		79053	1 Tape Strip 10 x 55 cm
89271		79273	1 Extremity Drape 193 x 305 cm: • Conformable Circular Fenestration 3 cm • CONTROL* Fabric Reinforcement • Tube Holders

HALYARD BASICS* EXTREMITY PACK II

CODE	DESCRIPTION		ITEMS/CASE
89684	HALYARD BASICS* Extremity Pack II		10
	STERILE	NON-STERILE	CONTENTS
88673		78669	1 Back Table Cover, Zone Reinforced 152 x 191 cm
	-	-	1 Extremity Drape 228 x 300 cm: • Oval Fenestration 8 x 6 cm, centered • Spunbond Film
89672		79694	1 Half Drape 112 x 150 cm
	-	-	1 Tape Strip 10 x 55 cm

All codes in **bold** indicate a sterile product.

HAND PACK I

CODE	DESCRIPTION		ITEMS/CASE
76226	Hand Pack I		7
	STERILE	NON-STERILE	CONTENTS
	42345	-	1 Back Table Cover, Reinforced 140 x 190 cm
	89601	79601	1 Mayo Stand Cover 58 x 137 cm
	90012	79970	1 Surgical Gown Set-in-Sleeves Large
	90043	79983	1 Surgical Gown Raglan Sleeves XL, X-Long
	89053	79053	2 Tape Strips 10 x 55 cm
	47622	47623	1 Large Drape 193 x 139 cm
	-	-	1 Hand Drape 193 x 305 cm • Attached Sleeve • Tube Holders • CONTROL* Plus Fabric Reinforcement

HAND PACK


CODE	DESCRIPTION		ITEMS/CASE
77783	Hand Pack		7
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 89 x 89 cm
	42222	42428	1 Back Table Cover, Reinforced 152 x 191 cm
	89053	79053	1 Tape Strip 10 x 55 cm
	77450	13227 (handi-bin)	1 Hand Drape 290 x 361 x 157 cm: • Conformable Circular Fenestration 3 cm • CONTROL* Plus Fabric Reinforcement • Tube Holders • Gusseted Hand Table Cover 157 x 147 cm • Positioning Tape

HALYARD BASICS* HAND PACK


CODE	DESCRIPTION		ITEMS/CASE
88910	HALYARD BASICS* Hand Pack		10
	STERILE	NON-STERILE	CONTENTS
	88673	78669	1 Back Table Cover, Zone Reinforced 152 x 191 cm
	47348	-	1 Mayo Stand Cover, Zone Reinforced 80 x 144 cm
	-	-	1 Hand Drape 290 x 361 cm: • Conformable Circular Fenestration 3 cm, centered • Absorbent Reinforcement
	89672	79694	1 Half Drape 112 x 150 cm

All codes in **bold** indicate a sterile product.

See Back Cover for Trademark Information



ORTHOPAEDIC PACKS



HALYARD BASICS* HAND PACK II

CODE	DESCRIPTION		ITEMS/CASE
89687	HALYARD BASICS* Hand Pack II		10
	STERILE	NON-STERILE	CONTENTS
88673		78669	1 Back Table Cover, Zone Reinforced 152 x 191 cm
-	-	-	1 Hand Drape 290 x 361 cm: • Conformable Circular Fenestration 3 cm, centered • Absorbent Reinforcement
89672		79694	1 Half Drape 112 x 150 cm

HALYARD BASICS* HAND & FOOT PACK

CODE	DESCRIPTION		ITEMS/CASE
88909	HALYARD BASICS* Hand & Foot Pack		10
	STERILE	NON-STERILE	CONTENTS
88673		78669	1 Back Table Cover, Zone Reinforced 152 x 191 cm
47348		-	1 Mayo Stand Cover, Zone Reinforced 80 x 144 cm
89673		-	1 Hand & Foot Drape 180 x 300 cm: • Conformable Circular Fenestration 3 cm, centered • Spunbond Film
89672		79694	1 Half Drape 112 x 150 cm
-	-	-	1 Tape Strip 10 x 55 cm

HALYARD BASICS* HAND & FOOT PACK II

CODE	DESCRIPTION		ITEMS/CASE
89685	HALYARD BASICS* Hand & Foot Pack II		10
	STERILE	NON-STERILE	CONTENTS
88673		78669	1 Back Table Cover, Zone Reinforced 152 x 191 cm
89673		-	1 Hand & Foot Drape 180 x 300 cm: • Conformable Circular Fenestration 3 cm, centered • Spunbond Film
89672		79694	1 Half Drape 112 x 150 cm
-	-	-	1 Tape Strip 10 x 55 cm

All codes in **bold** indicate a sterile product.

ARTHROSCOPY PACK II


CODE	DESCRIPTION	ITEMS/CASE
77881	Arthroscopy Pack II	10
STERILE	NON-STERILE	CONTENTS
-	-	1 Sterilisation Wrap 89 x 89 cm
42222	42428	1 Back Table Cover, Reinforced 152 x 191 cm
47614	47616	1 Half Drape 101 x 147 cm
89279	79279	1 Universal Lower Extremity Drape 229 x 315 cm: • Conformable Circular Fenestration 6 cm • CONTROL* Plus Fabric Reinforcement • Tube Holders

HALYARD BASICS* ARTHROSCOPIC CRUCIATE LIGAMENT (ACL) PACK

CODE	DESCRIPTION	ITEMS/CASE
88907	HALYARD BASICS* Arthroscopic Cruciate Ligament (ACL) Pack	5
STERILE	NON-STERILE	CONTENTS
88673	78669	1 Back Table Cover, Zone Reinforced 152 x 191 cm
47348	-	1 Mayo Stand Cover, Zone Reinforced 80 x 144 cm
-	-	1 Knee Arthroscopy Drape 229 x 315 cm
-	-	1 Stockinette Medium 25 x 122 cm
-	-	1 Tape Strip 10 x 55 cm
-	-	2 Absorbent Towels 33 x 39 cm


HALYARD BASICS* SHOULDER ARTHROSCOPY PACK

CODE	DESCRIPTION	ITEMS/CASE
88908	HALYARD BASICS* Shoulder Arthroscopy Pack	5
STERILE	NON-STERILE	CONTENTS
88673	78669	1 Back Table Cover, Zone Reinforced 152 x 191 cm
47348	-	1 Mayo Stand Cover, Zone Reinforced 80 x 144 cm
-	-	1 Shoulder Arthroscopy Drape 254 x 157 cm: • Fluid Collection Pouch • 10 x 53 cm Split
-	-	1 Shoulder Split Drape 254 x 168 cm: • 10 x 53 cm Split
-	-	1 Stockinette Medium 25 x 122 cm
-	-	1 Tape Strip 10 x 55 cm



All codes in **bold** indicate a sterile product.

ORTHOPAEDIC PACKS



LOWER EXTREMITY PACK I

CODE	DESCRIPTION		ITEMS/CASE
77790	Lower Extremity Pack I		13
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
42222	42428	-	1 Back Table Cover, Reinforced 152 x 191 cm
89053	79053	-	1 Tape Strip 10 x 55 cm
-	-	13232 (handi-bin)	1 Extremity Drape 264 x 193 x 315cm: • Conformable Circular Fenestration 6 cm • Armboard Covers • CONTROL* Plus Fabric Reinforcement • Tube Holders

LOWER EXTREMITY PACK II

CODE	DESCRIPTION		ITEMS/CASE
77798	Lower Extremity Pack II		9
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
42222	42428	-	1 Back Table Cover, Reinforced 152 x 191 cm
89603	79603	-	1 Mayo Stand Cover, 80 x 145 cm
89053	79053	-	1 Tape Strip 10 x 55 cm
-	-	13232 (handi-bin)	1 Extremity Drape 264 x 193 x 315cm: • Conformable Circular Fenestration 6 cm • Armboard Covers • CONTROL* Plus Fabric Reinforcement • Tube Holders

LOWER EXTREMITY PACK

CODE	DESCRIPTION		ITEMS/CASE
59706	Lower Extremity Pack		6
	STERILE	NON-STERILE	CONTENTS
-	-	42451	1 Back Table Cover, Reinforced 152 x 229 cm
89603	79603	-	1 Mayo Stand Cover 80 x 145 cm
-	-	-	1 Suture Bag
89701	79720	-	1 Absorbent Towel 38 x 56 cm
-	-	-	1 Large Drape 152 x 193 cm
89279	79279	-	1 Universal Extremity Drape • Conformable Circular Fenestration 6 cm • CONTROL* Plus Fabric Reinforcement • Tube Holders


All codes in **bold** indicate a sterile product.

LOWER EXTREMITY PACK

CODE	DESCRIPTION		ITEMS/CASE
76246	Lower Extremity Pack		5
	STERILE	NON-STERILE	CONTENTS
	42312	42313	1 Back Table Cover, Reinforced 191 x 201 cm
	89603	79603	1 Mayo Stand Cover 80 x 145 cm
	47622	47623	1 Large Drape 139 x 193 cm
	89053	79053	2 Tape Strips 10 x 55 cm
	89278	-	1 ORTHOARTS* Lower Extremity Drape 290 x 224 x 320 cm • Conformable Circular Fenestration 6 cm • Armboard Covers • CONTROL* Plus Fabric Reinforcement • Tube Holders • Heavy Duty Fabric


ORTHOARTS* PACK VII

CODE	DESCRIPTION		ITEMS/CASE
88203	ORTHOARTS* Pack VII		3
	STERILE	NON-STERILE	CONTENTS
	42224	42451	1 Back Table Cover, Reinforced 152 x 229 cm
	89603	79603	1 Mayo Stand Cover 80 x 145 cm
	89811	-	1 Impervious Stockinette, Medium 23 x 112 cm
	89331	79330	1 Impervious Split Drape 152 x 178 cm: • 10 x 53 cm Split • Positioning Tape
	89278	13233 (handi-bin)	1 ORTHOARTS* Lower Extremity Drape 290 x 320 x 224 cm: • Conformable Circular Fenestration 6 cm • Armboard Covers • CONTROL* Plus Fabric Reinforcement • Tube Holders • Heavy Duty Fabric
	-	-	2 Iso-Seal Tape Strips 12 x 59 cm



All codes in **bold** indicate a sterile product.

ORTHOPAEDIC PACKS



KNEE PACK I

CODE	DESCRIPTION		ITEMS/CASE
48756	Knee Pack I		9
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
42222		42428	1 Back Table Cover, Reinforced 152 x 191 cm
89603		79603	1 Mayo Stand Cover 80 x 145 cm
89053		79053	2 Tape Strips 10 x 55 cm
-		13232 (handi-bin)	1 Extremity Drape 264 x 193 x 315cm: • Conformable Circular Fenestration 6 cm • Armboard Covers • CONTROL[®] Plus Fabric Reinforcement • Tube Holders


KNEE PACK II

CODE	DESCRIPTION		ITEMS/CASE
48757	Knee Pack II		7
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
42222		42428	1 Back Table Cover, Reinforced 152 x 191 cm
47614		47616	1 Half Drape 101 x 147 cm
89053		79053	2 Tape Strips 10 x 55 cm
-		13232 (handi-bin)	1 Extremity Drape 264 x 193 x 315cm: • Conformable Circular Fenestration 6 cm • Armboard Covers • CONTROL[®] Plus Fabric Reinforcement • Tube Holders

All codes in **bold** indicate a sterile product.


ORTHOARTS* ARTHROSCOPY PACK II

CODE	DESCRIPTION		ITEMS/CASE
88213	ORTHOARTS* Arthroscopy Pack II		4
	STERILE	NON-STERILE	CONTENTS
42224		42451	1 Back Table Cover, Reinforced 152 x 229 cm
89603		79603	1 Mayo Stand Cover 80 x 145 cm
-	-	-	1 Suture Bag
-	-	-	4 Iso-Seal Tape Strips 12 x 59 cm
89821		-	1 Impervious Stockinette, Large 30 x 122 cm
89331		79330	1 Impervious Split Drape 152 x 178 cm: • 10 x 53 cm Split • Positioning Tape
89267		79267	1 ORTHOARTS* Knee Arthroscopy Drape 290 x 295 x 224 cm: • Heavy Duty Fabric • Fluid Collection Pouch with Reflux Shields and Suction Port • Armboard Covers • CONTROL* Plus Fabric Reinforcement • Tube Holders




HIP PACK

CODE	DESCRIPTION		ITEMS/CASE
48758	Hip Pack		4
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 96 x 96 cm
42224		42451	1 Back Table Cover, Reinforced 152 x 229 cm
89603		79603	1 Mayo Stand Cover 80 x 145 cm
89053		79053	1 Tape Strip 10 x 55 cm
89331		79330	1 Impervious Split Drape 152 x 178 cm: • 10 x 53 cm Split • Positioning Tape
-	-	76368	1 Ortho U-Drape 224 x 318 cm: • 15 x 137 cm Split • Positioning Tape • CONTROL* Plus Fabric Reinforcement
-	-	76363	1 Bar Drape 157 x 269 cm: • Armboard Covers • Positioning Tape • CONTROL* Plus Fabric Reinforcement • Tube Holders



All codes in **bold** indicate a sterile product.

ORTHOPAEDIC PACKS



HIP PACK

CODE	DESCRIPTION		ITEMS/CASE
76249	Hip Pack		5
	STERILE	NON-STERILE	CONTENTS
42312		42313	1 Back Table Cover, Reinforced 191 x 201 cm
89603		79603	1 Mayo Stand Cover 80 x 145 cm
-		-	1 Drape 119 x 152 cm • CONTROL ⁺ Plus Fabric Reinforcement
-		-	1 ORTHOARTS ⁺ U-Drape 224 x 274 cm • 12 x 101 cm Split • CONTROL ⁺ Plus Fabric Reinforcement • Positioning Tape • Heavy Duty Fabric
89348		79348	1 ORTHOARTS ⁺ Bar Drape 183 x 284 cm • Heavy Duty Fabric • CONTROL ⁺ Plus Fabric Reinforcement • Positioning Tape • Tube Holders

HIP PACK I

CODE	DESCRIPTION		ITEMS/CASE
48755	Hip Pack I		6
	STERILE	NON-STERILE	CONTENTS
42345		-	2 Back Table Covers, Reinforced 140 x 190 cm
89603		79603	1 Mayo Stand Cover 80 x 145 cm
89701		79720	1 Absorbent Towel 38 x 56 cm
89311		79363	1 Ortho U-Drape 224 x 318 cm: • 15 x 137 cm Split • Positioning Tape • CONTROL ⁺ Plus Fabric Reinforcement
89381		79360	1 Ortho Bar Drape 152 x 254 cm: • Positioning Tape • CONTROL ⁺ Fabric Reinforcement • Tube Holders
-		-	1 Fluid Collection Pouch with Tape 25 x 41 cm
89053		79053	1 Tape Strip 10 x 55 cm
-		-	1 Diathermy Pouch 29 x 28 cm, Heavy Duty, Double Tape


All codes in **bold** indicate a sterile product.

ORTHOARTS* HIP PACK II

CODE	DESCRIPTION		ITEMS/CASE
88253	ORTHOARTS* Hip Pack II		4
STERILE	NON-STERILE	CONTENTS	
-	42449	1 Back Table Cover, Reinforced 127 x 229 cm	
89603	79603	1 Mayo Stand Cover 80 x 145 cm	
-	-	1 Suture Bag	
89821	-	1 Impervious Stockinette, Large 30 x 122 cm	
-	79333	1 Large Impervious Split Drape 152 x 206 cm: • 11 x 78 cm Split • Positioning Tape	
-	79346	1 ORTHOARTS* Hip Drape 290 x 320 x 224 cm: • Heavy Duty Fabric • Armboard Covers • CONTROL* Plus Fabric Reinforcement • Tube Holders • Conformable Teardrop Fenestration	


ORTHOARTS* HIP PACK

CODE	DESCRIPTION		ITEMS/CASE
88183	ORTHOARTS* Hip Pack		3
STERILE	NON-STERILE	CONTENTS	
42224	42451	1 Back Table Cover, Reinforced 152 x 229 cm	
89603	79603	1 Mayo Stand Cover 80 x 145 cm	
-	-	1 Suture Bag	
89821	-	1 Impervious Stockinette, Large 30 x 122 cm	
-	79333	1 Large Impervious Split Drape 152 x 206 cm: • 11 x 78 cm Split • Positioning Tape	
89345	79345	1 ORTHOARTS* Hip Drape with Pockets 290 x 320 x 224 cm: • Heavy Duty Fabric • Armboard Covers • CONTROL* Plus Fabric Reinforcement • Tube Holders • Pockets on both sides • Conformable Teardrop Fenestration	



All codes in **bold** indicate a sterile product.

ORTHOPAEDIC PACKS



ORTHO MAJOR PACK III

CODE	DESCRIPTION		ITEMS/CASE
88373	ORTHO Major Pack III		7
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 115 x 115 cm
42222		42428	1 Back Table Cover, Reinforced 152 x 191 cm
47614		47616	1 Half Drape 101 x 147 cm
-		79333	1 Large Impervious Split Drape 152 x 206 cm: • 11 x 78 cm Split • Positioning Tape
-		76368	1 Ortho U-Drape 224 x 318 cm: • 15 x 137 cm Split • Positioning Tape • CONTROL* Plus Fabric Reinforcement
-		78961	1 Top Drape 193 x 239 cm: • Positioning Tape • CONTROL* Plus Fabric Reinforcement • Tube Holders


ORTHOARTS* SHOULDER PACK II

CODE	DESCRIPTION		ITEMS/CASE
88243	ORTHOARTS* Shoulder Pack II		5
	STERILE	NON-STERILE	CONTENTS
42224		42451	1 Back Table Cover, Reinforced 152 x 229 cm
89603		79603	1 Mayo Stand Cover 80 x 145 cm
89701		79720	2 Absorbent Towels 38 x 56 cm
-	-	-	1 Suture Bag
-	-	-	2 Iso-Seal Tape Strips 12 x 59 cm
89275		79275	1 ORTHOARTS* Body Split Drape 224 x 295 cm: • Heavy Duty Fabric • CONTROL* Plus Fabric Reinforcement • Positioning Tape • Tube Holders
-	-	-	1 ORTHOARTS* Shoulder Split Drape 152 x 284 cm: • Heavy Duty Fabric • CONTROL* Plus Fabric Reinforcement • 15 x 20 cm Bilateral Split • Positioning Tape • Tube Holders

All codes in **bold** indicate a sterile product.

ORTHOARTS* U-BAR* PACK

CODE	DESCRIPTION		ITEMS/CASE
88233	ORTHOARTS* U-BAR* Pack		4
	STERILE	NON-STERILE	CONTENTS
	42224	42451	1 Back Table Cover, Reinforced 152 x 229 cm
	89603	79603	1 Mayo Stand Cover 80 x 145 cm
	89701	79720	2 Absorbent Towels 38 x 56 cm
	-	-	1 Suture Bag
	89331	79330	1 Impervious Split Drape 152 x 178 cm: • 10 x 53 cm Split • Positioning Tape
	-	13294 (handi-bin)	1 ORTHOARTS* U-Drape 224 x 274 cm: • Heavy Duty Fabric • CONTROL* Fabric Reinforcement • Positioning Tape
	89348	79348	1 ORTHOARTS* Bar Drape 183 x 284 cm: • Heavy Duty Fabric • CONTROL* Fabric Reinforcement • Positioning Tape • Tube Holders



All codes in **bold** indicate a sterile product.

ORTHOPAEDIC DRAPES

ORTHOARTS* SPLIT / BAR DRAPES

ORTHOARTS* BODY SPLIT DRAPE, 224 X 295 CM

- 11 x 69 cm Split
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement 61 x 56 cm
- Tube Holders
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89275	ORTHOARTS* Body Split Drape, Sterile	10
79275	ORTHOARTS* Body Split Drape, Non-Sterile	12

ORTHOARTS* BAR DRAPE, 284 X 183 CM

- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement 41 x 61 cm
- Tube Holders
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89348	ORTHOARTS* Bar Drape, Sterile	12
79348	ORTHOARTS* Bar Drape, Non-Sterile	14

IMPERVIOUS SPLIT DRAPES

IMPERVIOUS SPLIT DRAPE, 152 X 178 CM


- 10 x 53 cm Split
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89331	Impervious Split Drape, Sterile	42
79330	Impervious Split Drape, Non-Sterile	115

HALYARD BASICS* IMPERVIOUS SPLIT DRAPE, 152 X 178 CM

- 10 x 53 cm Split
- Poly
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
78660	HALYARD BASICS* Impervious Split Drape, Non-Sterile	65



All codes in **bold** indicate a sterile product.

LARGE IMPERVIOUS SPLIT DRAPE, 152 X 206 CM

- 11 x 78 cm Split
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79333	Large Impervious Split Drape, Non-Sterile	100

SPLIT DRAPE, 150 X 230 CM

- 10 x 80 cm Split
- CONTROL* Plus Fabric Reinforcement 61 x 140 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
K793210	Split Drape, Non-Sterile	180

EXTREMITY DRAPE

BILATERAL LIMB DRAPE, 193 X 305 CM

- CONTROL* Plus Fabric Reinforcement 102 x 81 cm
- Two Conformable Circular Fenestrations, 6 cm
- Tube Holders


CODE	DESCRIPTION	ITEMS/CASE
89291	Bilateral Limb Drape, Sterile	14

LOWER EXTREMITY DRAPE

UNIVERSAL EXTREMITY DRAPE, 229 X 315 CM


- CONTROL* Plus Fabric Reinforcement 122 x 71 cm
- Conformable Circular Fenestration 6 cm
- Tube Holders

CODE	DESCRIPTION	ITEMS/CASE
89279	Universal Extremity Drape, Sterile	14
79279	Universal Extremity Drape, Non-Sterile	150
K892791	Universal Extremity Drape, Sterile	14



All codes in **bold** indicate a sterile product.

ORTHOPAEDIC DRAPE



LOWER EXTREMITY DRAPE, 224 X 333 CM

- CONTROL* Plus Fabric Reinforcement 81 x 81 cm
- Conformable Circular Fenestration 6 cm
- Tube Holders

CODE	DESCRIPTION	ITEMS/CASE
89276	Lower Extremity Drape, Sterile	14
13228	Lower Extremity Drape, Non-Sterile - handi-bin	192

LOWER EXTREMITY DRAPE, 274 X 315 X 203 CM

- CONTROL* Plus Fabric Reinforcement 51 x 91 x 122 cm
- Circular Fenestration 6 cm

CODE	DESCRIPTION	ITEMS/CASE
13232	Lower Extremity Drape, Non-Sterile	130

LOWER EXTREMITY DRAPE, 289 X 203 X 365 CM

- CONTROL* Plus Fabric Reinforcement 51 x 91 x 122 cm
- Circular Fenestration 6 cm
- SMS Fabric

CODE	DESCRIPTION	ITEMS/CASE
29222	Extra-Long Extremity Drape	13

ORTHOARTS* LOWER EXTREMITY DRAPE, 290 X 320 X 224 CM

- Armboard Covers
- Conformable Circular Fenestration 6 cm
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement 91 x 122 cm
- Tube Holders

CODE	DESCRIPTION	ITEMS/CASE
89278	ORTHOARTS* Lower Extremity Drape, Sterile	7
13233	ORTHOARTS* Lower Extremity Drape, Non-Sterile - handi-bin	84

ORTHOARTS* EXTRA LONG LOWER EXTREMITY DRAPE, 289 X 381 X 223 CM

- Heavy-duty SMS
- Absorbent CONTROL* Plus fabric reinforcement
- Elastic film with 6 cm circular fenestration
- Armboard covers
- Tube holders
- Extra length for coverage where needed

CODE	DESCRIPTION	ITEMS/CASE
29219	ORTHOARTS* Extra Long Lower Extremity Drape, Sterile	8
29220	ORTHOARTS* Extra Long Lower Extremity Drape, Non-Sterile	8

All codes in **bold** indicate a sterile product.

EXTRA LONG LOWER EXTREMITY DRAPE, 290 X 366 X 203 CM

- Heavy-duty SMS
- Absorbent CONTROL* Plus fabric reinforcement
- Elastic film with 6 cm circular fenestration
- Armboard covers
- Tube holders
- Extra length for coverage where needed

CODE	DESCRIPTION	ITEMS/CASE
29223	Extra Long Lower Extremity Drape, Non-Sterile	14

UPPER EXTREMITY DRAPE

UPPER EXTREMITY DRAPE, 193 X 305 CM

- CONTROL* Plus Fabric Reinforcement 61 x 102 cm
- Conformable Circular Fenestration 3 cm
- Tube Holders

CODE	DESCRIPTION	ITEMS/CASE
89271	Upper Extremity Drape, Sterile	14
79273	Upper Extremity Drape, Non-Sterile	20

SHOULDER DRAPE

SHOULDER ARTHROSCOPY DRAPE, 406 X 259 CM

- Attached Fluid Collection Pouch with Suction Port
- Hook and Loop Tube and Line Control
- Conformable Tear Drop Fenestration 6 cm
- Instrument Pad Fabric Reinforcement 53 x 66 cm

CODE	DESCRIPTION	ITEMS/CASE
89070	Shoulder Arthroscopy Drape, Sterile	10
79070	Shoulder Arthroscopy Drape, Non-Sterile	12


ORTHOARTS* SHOULDER ARTHROSCOPY DRAPE, 254 X 152 CM

- Attached Fluid Collection Pouch with Suction Ports
- Heavy-Duty Fabric
- Positioning Tape


CODE	DESCRIPTION	ITEMS/CASE
89492	ORTHOARTS* Shoulder Arthroscopy Drape, Sterile	12
79492	ORTHOARTS* Shoulder Arthroscopy Drape, Non-Sterile	15

All codes in **bold** indicate a sterile product.

See Back Cover for Trademark Information



ORTHOPAEDIC DRAPES



HAND DRAPES

HAND DRAPE, 290 X 361 CM

- CONTROL* Plus Fabric Reinforcement 86 x 41 cm
- Conformable Circular Fenestration 3 cm
- Gusseted Hand Table Cover 157 x 147 cm
- Tube Holders
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
77450	Hand Drape, Sterile	11
13227	Hand Drape, Non-Sterile - handi-bin	108

HALYARD BASICS* HAND & FOOT DRAPE, 180 X 300 CM

- Conformable Circular Fenestration 3 cm, centered
- Spunbond Film

CODE	DESCRIPTION	ITEMS/CASE
89673	HALYARD BASICS* Hand & Foot Drape, Sterile	20
79695	HALYARD BASICS* Hand & Foot Drape, Non-Sterile	30

HIP DRAPES

ORTHOARTS* HIP DRAPE, 290 X 320 X 224 CM

- Armboard Covers
- Conformable Teardrop Fenestration 21 x 17 cm
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement 152 x 61 cm
- Tube Holders

CODE	DESCRIPTION	ITEMS/CASE
79346	ORTHOARTS* Hip Drape, Non-Sterile	7

ORTHOARTS* HIP DRAPE WITH POCKETS, 290 X 320 X 224 CM

- Pockets on Both Sides
- Armboard Covers
- Conformable Teardrop Fenestration 21 x 17 cm
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement 61 x 152 cm
- Tube Holders

CODE	DESCRIPTION	ITEMS/CASE
89345	ORTHOARTS* Hip Drape with Pockets, Sterile	6
79345	ORTHOARTS* Hip Drape with Pockets, Non-Sterile	6

All codes in **bold** indicate a sterile product.

ANTERIOR HIP DRAPE WITH LEGGINGS, 401 X 231 CM

- Main Drape size : 401 x 231 cm
- Reinforcement size : 137 x 122 cm
- Fenestration size : 15 x 28 cm each one
- Basis weight: Thickness of base material: 78 GSM Tri-lamination Fabric
- High performance fabric as defined by EN 13795:2011
- Test results as per EN13795 available upon request

CODE	DESCRIPTION	ITEMS/CASE
K893462	Anterior Hip drape with leggings, Sterile	12
K793462	Anterior Hip drape with leggings, Non-Sterile	14

ANTERIOR HIP DRAPE WITH INCISE FILM AND FENESTRATION 15X20 CM

- Main drape size: 315 x 400 cm
- Reinforcement size: 122 x 178 cm
- Fenestration size: 15 x 21 cm (x2)
- Leggings: 80 x 122 cm (x2)

CODE	DESCRIPTION	ITEMS/CASE
K893463	Anterior Hip Drape with Incise Film and Fenestration, Sterile	12
K793463	Anterior Hip Drape with Incise Film and Fenestration, Non-Sterile	12

ANTERIOR HIP DRAPE WITH FENESTRATION 15X20 CM




- Main drape size: 315 x 400 cm
- Reinforcement size: 122 x 178 cm
- Fenestration size: 15 x 21 cm (x2)
- Leggings: 80 x 122 cm (x2)

CODE	DESCRIPTION	ITEMS/CASE
K893464	Anterior Hip Drape with Fenestration, Sterile	12
K793464	Anterior Hip Drape with Fenestration, Non-Sterile	12

ANTERIOR HIP DRAPE WITH INCISE FILM AND FENESTRATION 20X28 CM

- Main drape size: 315 x 400 cm
- Reinforcement size: 122 x 178 cm
- Fenestration size: 15 x 28 cm (x2)
- Leggings: 80 x 122 cm (x2)

CODE	DESCRIPTION	ITEMS/CASE
K893467	Anterior Hip Drape with Incise Film and Fenestration, Sterile	6
K793467	Anterior Hip Drape with Incise Film and Fenestration, Non-Sterile	6



ORTHOPAEDIC DRAPES

ANTERIOR HIP DRAPE WITH FENESTRATION 20X28 CM

- Main drape size: 315 x 400 cm
- Reinforcement size: 122 x 178 cm
- Fenestration size: 15 x 28 cm (x2)
- Leggings: 80 x 122 cm (x2)

CODE	DESCRIPTION	ITEMS/CASE
K893468	Anterior Hip Drape with Fenestration, Sterile	6
K793468	Anterior Hip Drape with Fenestration, Non-Sterile	6

VERTICAL ISOLATION DRAPES

VERTICAL ISOLATION DRAPE, 320 X 229 CM

- Clear Film
- Adhesive Incise Area 41 x 76 cm
- Fluid Collection Pouch with Suction Port
- Tube Holders
- With Tape Strip for Attachment to Horizontal Rod
- Attachable Instrument Pouch

CODE	DESCRIPTION	ITEMS/CASE
77723	Vertical Isolation Drape, Sterile	14
76375	Vertical Isolation Drape, Non-Sterile	15

VERTICAL ISOLATION DRAPE, 320 X 279 CM

- Clear Film
- Adhesive Incise Area
- Fluid Collection Pouch with Suction Port
- Tube Holders
- With Tape Strip for Attachment to Horizontal Rod


CODE	DESCRIPTION	ITEMS/CASE
33415	Vertical Isolation Drape, Sterile	9
33412	Vertical Isolation Drape, Non-Sterile	11

KNEE ARTHROSCOPY DRAPES

KNEE ARTHROSCOPY DRAPE, 229 X 315 CM

- Attached Fluid Collection Pouch with Suction Port
- Conformable Circular Fenestration 6 cm
- Tube Holders

CODE	DESCRIPTION	ITEMS/CASE
89265	Knee Arthroscopy Drape, Sterile	10
79265	Knee Arthroscopy Drape, Non-Sterile	150



All codes in **bold** indicate a sterile product.

ORTHOARTS* KNEE ARTHROSCOPY DRAPE, 290 X 295 X 224 CM

- Armboard Covers
- Fluid Collection Pouch with Reflux Shields and Suction Port
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement 91 x 121 cm
- Tube Holders

CODE	DESCRIPTION	ITEMS/CASE
89267	ORTHOARTS* Knee Arthroscopy Drape, Sterile	7
79267	ORTHOARTS* Knee Arthroscopy Drape, Non-Sterile	7

UNIVERSAL SPINE DRAPES

UNIVERSAL SPINE DRAPE, 284 X 350 X 251 CM



- SMS
- Absorbent CONTROL* Plus fabric reinforcement
- 10 x 30 cm fenestration
- Armboard covers
- Hook-n-loop tube holders
- Clear instrument pouch

CODE	DESCRIPTION	ITEMS/CASE
29221	Universal Spine Drape, Sterile	8
29224	Universal Spine Drape, Non-Sterile	9

U-DRAPES

ORTHOPAEDIC U-DRAPE

CODE	DESCRIPTION	ITEMS/CASE
K793632	Orthopaedic U-drape, Non-Sterile	22



GYNAECOLOGY & OBSTETRICS

GYNAECOLOGY & OBSTETRICS PACKS

C-SECTION PACK I

CODE	DESCRIPTION		ITEMS/CASE
77454	C-Section Pack I		9
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 96 x 96 cm
42222		42428	1 Back Table Cover, Reinforced 152 x 191 cm
89701		79720	1 Absorbent Towel 38 x 56 cm
-		76394	1 Baby Blanket 80 x 130 cm: • Absorbent
89529		79529	1 C-Section Fluid Collection Drape 254 x 305 x 183 cm: • Full Adhesive Incise Area • Fluid Collection Pouch with Suction Port and Malleable Band • Armboard Covers • Tube Holders

C-SECTION PACK II


CODE	DESCRIPTION		ITEMS/CASE
77455	C-Section Pack II		8
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 96 x 96 cm
42222		42428	1 Back Table Cover, Reinforced 152 x 191 cm
89701		79720	1 Absorbent Towel 38 x 56 cm
-		76394	1 Baby Blanket 80 x 130 cm: • Absorbent
89528		79528	1 C-Section Fluid Collection Drape 254 x 305 x 183 cm: • Triangular Fenestrated Adhesive Incise Area • Fluid Collection Pouch with Suction Port and Malleable Band • Armboard Covers • Tube Holders

All codes in **bold** indicate a sterile product.

GYNAECOLOGY & OBSTETRICS PACKS

C-SECTION PACK III

CODE	DESCRIPTION		ITEMS/CASE
88354	C-Section Pack III		11
	STERILE	NON-STERILE	CONTENTS
	42222	42428	1 Back Table Cover, Reinforced 152 x 191 cm
	89603	79603	1 Mayo Stand Cover 80 x 145 cm
	89701	79720	4 Absorbent Towels 38 x 56 cm
		76353	1 C-Section Drape with Attached Leggings 295 x 318 x 224 cm: • Armboard Covers • Attached Leggings • 36 x 33 cm Adhesive Incise Area with Triangular Fenestration • Fluid Collection Pouch with Suction Port
	-	76394	1 Baby Blanket 80 x 130 cm: • Absorbent
	89053	79053	2 Tape Strips 10 x 55 cm




C-SECTION PACK IV

CODE	DESCRIPTION		ITEMS/CASE
77491	C-Section Pack IV		10
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 96 x 96 cm
	42222	42428	1 Back Table Cover, Reinforced 152 x 191 cm
	89603	79603	1 Mayo Stand Cover 80 x 145 cm
	89701	79720	1 Absorbent Towel 38 x 56 cm
	-	76394	1 Baby Blanket 80 x 130 cm: • Absorbent
	89529	79529	1 C-Section Fluid Collection Drape 254 x 305 x 183 cm: • Full Adhesive Incise Area • Fluid Collection Pouch with Suction Port and Malleable Band • Armboard Covers • Tube Holders

All codes in **bold** indicate a sterile product.


GYNAECOLOGY & OBSTETRICS PACKS

C-SECTION PACK V



CODE	DESCRIPTION		ITEMS/CASE
77492	C-Section Pack V		8
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 96 x 96 cm
42222		42428	1 Back Table Cover, Reinforced 152 x 191 cm
89603		79603	1 Mayo Stand Cover 80 x 145 cm
89701		79720	1 Absorbent Towel 38 x 56 cm
	-	76394	1 Baby Blanket 80 x 130 cm: • Absorbent
89528		79528	1 C-Section Fluid Collection Drape 254 x 305 x 183 cm: • Triangular Fenestrated Adhesive Incise Area • Fluid Collection Pouch with Suction Port and Malleable Band • Armboard Covers • Tube Holders

GYNAECOLOGY PACK



CODE	DESCRIPTION		ITEMS/CASE
77865	Gynaecology Pack		10
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 96 x 96 cm
42345		42365	1 Back Table Cover, Reinforced 140 x 190 cm
89701		79720	2 Absorbent Towels 38 x 56 cm
	-	79411	1 Under Buttocks Drape, Zone Reinforced, Medium, 97 x 112 cm
89408		79408	2 Leggings 79 x 122 cm: • With 15 cm Cuff
	-	-	1 Coeloscopy Sheet 193 x 238 cm

All codes in **bold** indicate a sterile product.


GYNAECOLOGY & OBSTETRICS PACKS

GYNAECOLOGY/UROLOGY PACK I

CODE	DESCRIPTION		ITEMS/CASE
77888	Gynaecology/Urology Pack I		12
STERILE	NON-STERILE	CONTENTS	
-	-	1 Sterilisation Wrap 96 x 96 cm	
42222	42428	1 Back Table Cover, Reinforced 152 x 191 cm	
89701	79720	1 Absorbent Towel 38 x 56 cm	
-	76888	1 Gynaecology/Urology Drape 254 x 287 cm: • Adjustable Fenestration with Flap • Fluid Collection Pouch with Screen and Suction Port • With Attached Leggings • Perineal Incise Film • CONTROL[®] Fabric Reinforcement • Positioning Tape	


GYNAECOLOGY/UROLOGY PACK II

CODE	DESCRIPTION		ITEMS/CASE
77890	Gynaecology/Urology Pack II		11
STERILE	NON-STERILE	CONTENTS	
-	-	1 Sterilisation Wrap 96 x 96 cm	
42222	42428	1 Back Table Cover, Reinforced 152 x 191 cm	
89603	79603	1 Mayo Stand Cover 80 x 145 cm	
89701	79720	1 Absorbent Towel 38 x 56 cm	
-	76888	1 Gynaecology/Urology Drape 254 x 287 cm: • Adjustable Fenestration with Flap • Fluid Collection Pouch with Screen and Suction Port • With Attached Leggings • Perineal Incise Film • CONTROL[®] Fabric Reinforcement • Positioning Tape	



All codes in **bold** indicate a sterile product.

GYNAECOLOGY & OBSTETRICS PACKS



LITHOTOMY PACK I

CODE	DESCRIPTION		ITEMS/CASE
88511	Lithotomy Pack I		10
	STERILE	NON-STERILE	CONTENTS
42217			1 Back Table Cover, Reinforced 112 x 229 cm
89701		79720	1 Absorbent Towel 38 x 56 cm
-	-	-	1 Lithotomy Drape with Attached Leggings, 221 x 254 x 76 cm: • 15 x 7 cm Perineal Fenestration • Positioning Tape • CONTROL* Fabric Reinforcement • Tube Holders

LITHOTOMY PACK II

CODE	DESCRIPTION		ITEMS/CASE
88521	Lithotomy Pack II		10
	STERILE	NON-STERILE	CONTENTS
42217			1 Back Table Cover, Reinforced 112 x 229 cm
89701		79720	1 Absorbent Towel 38 x 56 cm
-		79731	1 Utility Drape with Tape, 38 x 64 cm
-	-	-	1 Under Buttocks Drape, Zone Reinforced, Medium 102 x 117 cm
89408		79408	2 Leggings 79 x 122 cm: • With 15 cm Cuff
47614		47616	1 Half Drape 101 x 147 cm

LITHOTOMY PACK IV


CODE	DESCRIPTION		ITEMS/CASE
77779	Lithotomy Pack IV		12
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
42222		42428	1 Back Table Cover, Reinforced 152 x 191 cm
89411		-	1 Under Buttocks Drape, Zone Reinforced, Medium 112 x 75 cm
-		76384	1 Lithotomy Drape with Attached Leggings, 254 x 272 x 76 cm: • Tape around Fenestration • CONTROL* Plus Fabric Reinforcement • Tube Holders

All codes in **bold** indicate a sterile product.

GYNAECOLOGY & OBSTETRICS PACKS


LITHOTOMY PACK V

CODE	DESCRIPTION		ITEMS/CASE
88551	Lithotomy Pack V		12
	STERILE	NON-STERILE	CONTENTS
	-	42302	1 Back Table Cover, Reinforced 112 x 224 cm
89701		79720	1 Absorbent Towel 38 x 56 cm
	-	79731	1 Utility Drape with Tape, 38 x 64 cm
89408		79408	2 Leggings 79 x 122 cm: • With 15 cm Cuff
	-	79551	1 Lithotomy T-Drape 152 x 79 x 191 cm: • 7 x 15 cm Fenestration • Positioning Tape • CONTROL [®] Fabric Reinforcement • Tube Holders



HALYARD BASICS^{*} DELIVERY PACK

CODE	DESCRIPTION		ITEMS/CASE
88911	HALYARD BASICS [*] Delivery Pack		5
	STERILE	NON-STERILE	CONTENTS
89632		79741	1 Back Table Cover, Zone Reinforced 100 x 150 cm
89677		79699	1 Utility Drape 75 x 75 cm
89584		-	1 Under Buttocks Drape 102 x 112 cm: • With Graduated Fluid Collection Pouch • Spunbond
	-	-	1 Baby Blanket



All codes in **bold** indicate a sterile product.

GYNAECOLOGY & OBSTETRICS DRAPES

GYNAECOLOGY/UROLOGY DRAPES

GYNAECOLOGY/UROLOGY DRAPE, 254 X 287 CM

- Perineal Incise Film
- Fluid Collection Pouch with Screen and Suction Port
- CONTROL* Fabric Reinforcement 39 x 72 cm
- Adjustable Fenestration with Flap
- Attached Leggings
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
76888	Gynaecology/Urology Drape, Non-Sterile	22

C-SECTION DRAPES

C-SECTION FLUID COLLECTION DRAPE, 254 X 305 X 183 CM

- Triangular Fenestrated Adhesive Incise Area
- Fluid Collection Pouch with Suction Port and Malleable Band
- Tube Holders
- Armboard Covers

CODE	DESCRIPTION	ITEMS/CASE
89528	C-Section Fluid Collection Drape with Fenestration, Sterile	8
79528	C-Section Fluid Collection Drape with Fenestration, Non-Sterile	150

C-SECTION FLUID COLLECTION DRAPE, 254 X 305 X 183 CM

- Full Adhesive Incise Area
- Fluid Collection Pouch with Suction Port and Malleable Band
- Tube Holders
- Armboard Covers

CODE	DESCRIPTION	ITEMS/CASE
89529	C-Section Fluid Collection Drape, Sterile	8
79529	C-Section Fluid Collection Drape, Non-Sterile	150

All codes in **bold** indicate a sterile product.

GYNAECOLOGY & OBSTETRICS DRAPES

C-SECTION DRAPE WITH ATTACHED LEGGINGS, 295 X 318 X 224 CM

- Armboard Covers
- Attached Leggings 36 x 33 cm
- Adhesive Incise Area with Triangular Fenestration
- Fluid Collection Pouch with Suction Port

CODE	DESCRIPTION	ITEMS/CASE
76353	C-Section Drape with Attached Leggings, Non-Sterile	112

GENTLE C-SECTION DRAPE WITH FLUID COLLECTION POUCH, 246 X 297 CM




- SMS
- CONTROL* Fabric Reinforcement 43 x 56 cm
- Fluid Collection Pouch
- Full Incise
- Clear window fenestration 39.4 x 49.5 cm
- Fenestration 43.2 x 38.1 cm
- Two Flaps 76.2 x 63.5 cm
- Tube Holders
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
K897213	Gentle C-Section Drape with Fluid Collection Pouch, Sterile	15
K797213	Gentle C-Section Drape with Fluid Collection Pouch, Sterile	20

GENTLE C-SECTION DRAPE WITH FLUID COLLECTION POUCH, 246 X 297 CM


- SMS
- CONTROL* Fabric Reinforcement 43 x 56 cm
- Fluid Collection Pouch
- 22.9 x 17.8 cm Incise Triangle Fenestration
- Window Fenestration 39.4 x 49.5 cm
- Two Flaps 76.2 x 63.5 cm
- Tube Holders

CODE	DESCRIPTION	ITEMS/CASE
K897215	Gentle C-Section Drape with Fluid Collection Pouch, Sterile	15
K797215	Gentle C-Section Drape with Fluid Collection Pouch, Non-Sterile	18



All codes in **bold** indicate a sterile product.

GYNAECOLOGY & OBSTETRICS DRAPES



LITHOTOMY DRAPES

LITHOTOMY DRAPE WITH ATTACHED LEGGINGS, 254 X 272 X 76 CM

- CONTROL* Plus Fabric Reinforcement 61 x 27 cm
- Tape around Fenestration
- Tube Holders

CODE	DESCRIPTION	ITEMS/CASE
76384	Lithotomy Drape with Attached Leggings, Non-Sterile	11

All codes in **bold** indicate a sterile product.

UROLOGY


UROLOGY PACKS

CYSTOSCOPY PACK I

CODE	DESCRIPTION		ITEMS/CASE
88611	Cystoscopy Pack I		12
STERILE	NON-STERILE	CONTENTS	
-	42447	1 Back Table Cover, Reinforced 112 x 198 cm	
89701	79720	1 Absorbent Towel 38 x 56 cm	
89408	79408	2 Leggings 79 x 122 cm: • With 15 cm Cuff	
-	79540	1 Cystoscopy T-Drape 152 x 76 x 213 cm: • Perineal Fenestration 7 x 15 cm • 22 x 22 cm Screen • Positioning Tape • CONTROL [®] Fabric Reinforcement • Tube Holders	


CYSTOSCOPY PACK II

CODE	DESCRIPTION		ITEMS/CASE
88621	Cystoscopy Pack II		12
STERILE	NON-STERILE	CONTENTS	
42222	42428	1 Back Table Cover, Reinforced 152 x 191 cm	
89701	79720	1 Absorbent Towel 38 x 56 cm	
-	79545	1 Cystoscopy Drape 254 x 152 x 221 cm: • With Attached Leggings • Perineal Fenestration 7 x 15 cm • 22 x 22 cm Screen • Positioning Tape • CONTROL [®] Fabric Reinforcement • Tube Holders	



All codes in **bold** indicate a sterile product.

UROLOGY PACKS & DRAPES



HALYARD BASICS* CYSTOSCOPY PACK

CODE	DESCRIPTION		ITEMS/CASE
89681	HALYARD BASICS* Cystoscopy Pack		10
	STERILE	NON-STERILE	CONTENTS
	88673	78669	1 Back Table Cover, Zone Reinforced 152 x 191 cm
	-	-	2 Leggings 79 x 122 cm
	89679	79783	1 Utility Drape with Tape 90 x 90 cm: • Spunbond Film
	-	79742	1 Cystoscopy Drape, 75 x 185 cm: • Diamond Fenestration 5.5 x 12.5 cm • Spunbond Film

T.U.R. PACK

CODE	DESCRIPTION		ITEMS/CASE
77885	T.U.R. Pack		12
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 96 x 96 cm
	42222	42428	1 Back Table Cover, Reinforced 152 x 191 cm
	89701	79720	1 Absorbent Towel 38 x 56 cm
	-	76885	1 T.U.R. Drape 175 x 238 x 170 cm: • With Attached Leggings • Fluid Collection Pouch with Screen, Suction Port and Tie • With Finger Cot

T.U.R. DRAPE MESH & LEGGINGS 220 X 255 X 190 CM

- With Attached Leggings
- For Manufacturing/Processing/Repacking
- With Drape, Delnet, Finger Cot, Reinforcement

CODE	DESCRIPTION	ITEMS/CASE
K768550	T.U.R. Drape w/Mesh & Builtin Leggings, Non-Sterile	13

All codes in **bold** indicate a sterile product.

CYSTOSCOPY DRAPE, 254 X 221 X 152 CM

- CONTROL* Plus Fabric Reinforcement 40 x 22 cm & Poly Reinforcement 40 x 101 cm
- Attached Leggings
- Fenestration 7 x 15 cm
- Screen 22 x 22 cm

CODE	DESCRIPTION	ITEMS/CASE
79545	Cystoscopy Drape, Non-Sterile	20

HALYARD BASICS* CYSTOSCOPY DRAPE, 75 X 185 CM

- Diamond Fenestration 5.5 x 12.5 cm
- Spunbond Film

CODE	DESCRIPTION	ITEMS/CASE
79742	HALYARD BASICS* Cystoscopy Drape, Non-Sterile	120

T.U.R. DRAPE, 175 X 238 X 170 CM


- With Attached Leggings
- Fluid Collection Pouch with Screen, Suction Port and Tie
- With Finger Cot

CODE	DESCRIPTION	ITEMS/CASE
76885	T.U.R. Drape, Non-Sterile	18

GYN/UROLOGY DRAPE


- Attached Fluid Collection Pouch
- With Suction Port/Screen
- Attached Leggings
- Positioning Tape
- CONTROL* Plus Fabric Reinforcement 36 x 43 cm
- Fenestration 8 x 15 cm
- Tube Holder/2 holes

CODE	DESCRIPTION	ITEMS/CASE
89204	Gyn/Urology Drape, Sterile	10



OPHTHALMOLOGY

OPHTHALMOLOGY PACKS



EYE PACK

CODE	DESCRIPTION		ITEMS/CASE
77744	Eye Pack		10
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 89 x 89 cm
	-	-	1 Back Table Cover, 140 x 160 cm, with Tape
	-	-	1 Eye Drape 157 x 279 cm: • With Fluid Collection Pouch

OPHTHALMIC PACK

CODE	DESCRIPTION		ITEMS/CASE
77760	Ophthalmic Pack		8
	STERILE	NON-STERILE	CONTENTS
	42312	42313	1 Back Table Cover, Reinforced 191 x 201 cm
	47614	47616	1 Half Drape 101 x 147 cm
	90013	79973	2 Surgical Gowns Large, X-long
	89701	79720	2 Absorbent Towels 38 x 56 cm
	89304	79304	1 Ophthalmic Drape 132 x 142 cm • With Fluid Collection Pouch • 13 x 9 cm Rectangular Fenestration, with Incise • Applicable for both Left and Right Eye
	76214	76300	2 Armboard Covers 61 x 72 cm
	89945	79945	1 Diathermy Pouch 38 x 41 cm

All codes in **bold** indicate a sterile product.

OPHTHALMIC DRAPE, 132 X 142 CM

- With One Single Fluid Collection Pouch
- 13 x 9 cm Rectangular Fenestration, with Incise
- Applicable for both Left and Right Eye

CODE	DESCRIPTION	ITEMS/CASE
89304	Ophthalmic Drape, Sterile	40
79304	Ophthalmic Drape, Non-Sterile	100

OPHTHALMIC DRAPE, 132 X 142 CM




- With One Single Fluid Collection Pouch
- 8 cm Oval Fenestration, with Apertured Incise 9 x 12.7 cm
- Applicable for both Left and Right Eye

CODE	DESCRIPTION	ITEMS/CASE
89305	Ophthalmic Drape, Oval Fenestration, Sterile	40
79305	Ophthalmic Drape, Oval Fenestration, Non-Sterile	100


OPHTHALMIC DRAPE, 140 X 127 CM

- Supadrape Cellulose based Material with Clear Plastic Reinforcement around Fenestration
- 11.5 x 8.8 cm Rectangular Fenestration, with Incise
- With One Single Fluid Collection Pouch
- Applicable for both Left and Right Eye

CODE	DESCRIPTION	ITEMS/CASE
K793043	Ophthalmic Drape, Non-Sterile	110



OPHTHALMOLOGY DRAPES



OPHTHALMIC DRAPE, 176 X 193 CM

- SMS
- With Fluid Collection Pouch
- Fenestration 7 x 17 cm

CODE	DESCRIPTION	ITEMS/CASE
K793050	Ophthalmic Drape, Non-Sterile	120

OCUARTS* OPHTHALMIC DRAPE, 183 X 292 CM

- 3 x 5 cm Oval Fenestration
- 15 x 20 cm Incise
- Fluid Collection Pouch around Fenestration, with Foam Band
- Breather Bar
- Tube Holders

CODE	DESCRIPTION	ITEMS/CASE
79356	OCUARTS* Ophthalmic Drape, Non-Sterile	20

OPHTHALMIC LASIK DRAPE 38 X 20.3CM

- Bi-Laminated Drape Size: 20 x 38 cm
- Fenestration Size with incise: 2 x 6 cm

CODE	DESCRIPTION	ITEMS/CASE
K793053	Ophthalmic Lasik Drape , Non-Sterile	1200

EYE DRAPE, 140 X 193 CM

- With One Single Fluid Collection Pouch
- 7,6 x 10 cm Rectangular Fenestration, with Incise 18 x 13 cm
- Applicable for both Left and Right Eye

CODE	DESCRIPTION	ITEMS/CASE
77873	Eye Drape, Sterile	25
76321	Eye Drape, Non-Sterile	50

All codes in **bold** indicate a sterile product.

EYE DRAPE, 140 X 193 CM

- With Double Fluid Collection Pouches
- 4,5 x 7 cm Oval Fenestration, without Incise
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
77707	Eye Drape, Sterile	30
76346	Eye Drape, Non-Sterile	36

EYE DRAPE, 176 X 193 CM




- With Double Fluid Collection Pouches
- 4,5 x 7 cm Oval Fenestration, with Incise 9,5 x 12 cm

CODE	DESCRIPTION	ITEMS/CASE
59710	Eye Drape, Sterile	28
76393	Eye Drape, Non-Sterile	50

HALYARD BASICS* BILATERAL EYE FENESTRATED DRAPE WITH TAPE, 120 X 150 CM


- Horizontal Rounded Rectangular Fenestration 7 x 16 cm, centered
- Spunbond Film
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79744	HALYARD BASICS* Bilateral Eye Fenestrated Drape with Tape, Non-Sterile	80



ENT


ENT PACKS



ENT PACK III

CODE	DESCRIPTION		ITEMS/CASE
88223	ENT Pack III		6
	STERILE	NON-STERILE	CONTENTS
42222		42428	1 Back Table Cover, Reinforced, 152 x 191 cm
89603		79603	1 Mayo Stand Cover 80 x 145 cm
-		-	1 Suture Bag
-		79430	1 Head Drape 91 x 127 cm: • Turban • Positioning Tape
89352		79351	1 ENT Split Drape 193 x 315 cm: • With 7 x 71 cm Split • Positioning Tape • CONTROL* Fabric Reinforcement • Tube Holders • 2 SMS tube holders

ENT PACK IV



CODE	DESCRIPTION		ITEMS/CASE
77823	ENT Pack IV		10
	STERILE	NON-STERILE	CONTENTS
-		-	1 Sterilisation Wrap 89 x 89 cm
42222		42428	1 Back Table Cover, Reinforced, 152 x 191 cm
-		-	2 Iso-Seal Tape Strips 12 x 59 cm
-		-	1 ENT Drape with Incise, 193 x 224 cm: • Incise 15 x 15 cm • With 11 x 11 cm Fenestration • CONTROL* Plus Fabric Reinforcement



All codes in **bold** indicate a sterile product.

ENT PACK V – HEAD AND NECK





CODE	DESCRIPTION		ITEMS/CASE
76209	ENT Pack V – Head and Neck		7
STERILE	NON-STERILE	CONTENTS	
-	-	1 Sterilisation Wrap 89 x 89 cm	
42224	42451	1 Back Table Cover, Reinforced, 152 x 229 cm	
42345	42365	1 Back Table Cover, Reinforced, 140 x 190 cm	
-	79430	1 Head Drape 91 x 127 cm: • Turban • Positioning Tape	
89352	79351	1 ENT Split Drape 193 x 315 cm: • With 7 x 71 cm Split • Positioning Tape • CONTROL® Fabric Reinforcement • Tube Holders	

HEAD PACK

CODE	DESCRIPTION		ITEMS/CASE
76234	Head Pack		20
STERILE	NON-STERILE	CONTENTS	
-	-	1 Sterilisation Wrap 89 x 89 cm	
-	-	1 Three-Quarter Drape 152 x 193 cm	
89353	76323	1 Utility Drape 90 x 90 cm	
-	-	1 Head Drape, 61 x 84 cm: • 5 cm tape over full length	



ENT DRAPES



ENT DRAPES

ENT SPLIT DRAPE, 193 X 315 CM

- CONTROL* Fabric Reinforcement 58 x 61 cm
- 7 x 71 cm Split
- Tube Holders
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89352	ENT Split Drape, Sterile	18
79351	ENT Split Drape, Non-Sterile	25

ENT DRAPE, 193 X 229 CM

- With 5 x 6 cm Fenestration
- CONTROL* Plus Fabric Reinforcement 84 x 81 cm
- Positioning Tape
- Tube Holders

CODE	DESCRIPTION	ITEMS/CASE
89010	ENT Drape, Sterile	18

HEAD DRAPES

HEAD DRAPE, 127 X 91 CM

- Turban
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79430	Head Drape, Non-Sterile	60

DENTAL DRAPES

HALYARD BASICS* DENTAL FENESTRATED DRAPE, 120 X 150 CM

- Tear Drop Fenestration 10 x 15 cm, 40 cm from top
- SMS
- Absorbent Reinforcement

CODE	DESCRIPTION	ITEMS/CASE
89543	HALYARD BASICS* Dental Fenestrated Drape, Sterile	15
79543	HALYARD BASICS* Dental Fenestrated Drape, Non-Sterile	68

All codes in **bold** indicate a sterile product.

NEUROLOGY

CRANIOTOMY DRAPES

CRANIOTOMY DRAPE, 310 X 340 X 188 CM


- CONTROL* Plus Fabric Reinforcement 91 x 61 cm
- 20 x 30 cm Oval Fenestration
- Clear Gusseted Anesthesia Panels
- Screen and Suction Port
- 58 x 114 cm Pouch

CODE	DESCRIPTION	ITEMS/CASE
89062	Craniotomy Drape, Sterile	9
79075	Craniotomy Drape, Non-Sterile	15

CRANIOTOMY DRAPE, 310 X 340 X 188 CM

- CONTROL* Plus Fabric Reinforcement 91 x 61 cm
- 15 x 23 cm Oval Fenestration
- 58 x 92 cm Pouch with Screen and Suction Port
- Gusseted Anaesthesia Panels

CODE	DESCRIPTION	ITEMS/CASE
77466	Craniotomy Drape, Sterile	12
76396	Craniotomy Drape, Non-Sterile	14



All codes in **bold** indicate a sterile product.

CARDIOVASCULAR & ANGIOGRAPHY

CARDIOVASCULAR & ANGIOGRAPHY PACKS

CARDIO PACK II

CODE	DESCRIPTION		ITEMS/CASE
77748	Cardio Pack II		3
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 115 x 115 cm
	42312	42313	1 Back Table Cover, Reinforced, 191 x 201 cm
	89603	79603	1 Mayo Stand Cover, 80 x 145 cm
	-	76340	2 Foot Covers
	89701	79720	3 Absorbent Towels 38 x 56 cm
	89053	79053	2 Tape Strips 10 x 55 cm
	-	-	1 Suture Bag: • Double Tape
	89313	79313	1 Utility Drape with Tape, 76 x 76 cm
	-	76356	1 Perineal Drape 187 x 241 cm: • Split 44 x 95 cm • Positioning Tape • CONTROL [®] Plus Fabric Reinforcement
	86357	76357	1 CV Incise Sheet 264 x 363 x 193 cm: • With Pouches • CONTROL [®] Plus Fabric Reinforcement • Tube Holders

CARDIO SPLIT PACK


CODE	DESCRIPTION		ITEMS/CASE
77749	Cardio Split Pack		3
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 115 x 115 cm
	42224	42451	1 Back Table Cover, Reinforced, 152 x 229 cm
	89053	79053	2 Tape Strips 10 x 55 cm
	-	-	2 Suture Bags: • Double Tape
	42345	-	1 Back Table Cover, Reinforced 140 x 190 cm
	89313	79313	1 Utility Drape with Tape, 76 x 76 cm
	-	76356	1 Perineal Drape 187 x 241 cm: • Split 44 x 95 cm • Positioning Tape • CONTROL [®] Plus Fabric Reinforcement
	-	76378	1 CV Split Drape 320 x 224 x 295 cm: • Positioning Tape • CONTROL [®] Plus Fabric Reinforcement
	-	78961	1 Bottom Drape 193 x 239 cm: • Positioning Tape • CONTROL [®] Plus Fabric Reinforcement • Tube Holders

All codes in **bold** indicate a sterile product.

CARDIOVASCULAR & ANGIOGRAPHY PACKS


DILATATION PACK

CODE	DESCRIPTION		ITEMS/CASE
76248	Dilatation Pack		8
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
42345	-	-	1 Back Table Cover, Reinforced, 140 x 190 cm
-	-	-	1 X-Ray Protection Cover
89701		79720	1 Absorbent Towel 38 x 56 cm
89048		79109	1 Femoral Angiography Drape 211 x 330 cm: • 68 x 330 cm Window • Two 12 cm Diameter Fenestrations • CONTROL⁺ Plus Fabric Reinforcement 81 x 122 cm • Positioning Tape




CARDIOVASCULAR PACK III

CODE	DESCRIPTION		ITEMS/CASE
77892	Cardiovascular Pack III		3
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 115 x 115 cm
42224		42451	1 Back Table Cover, Reinforced, 152 x 229 cm
89053		79053	2 Tape Strips 10 x 55 cm
89701		79720	2 Absorbent Towels 38 x 56 cm
89603		79603	2 Mayo Stand Covers 80 x 145 cm
-	-	-	1 Perineal Cover 15 x 38 cm
89945		79945	2 Diathermy Pouches 38 x 41 cm
-	-	-	1 CV Split Drape 295 x 224 x 356 cm • 17 x 233 cm Split • CONTROL⁺ Plus Fabric Reinforcement • Armboard Covers
89381		79360	1 Ortho Bar Drape 152 x 254 cm
47614		47616	2 Half Drapes 147 x 101 cm



All codes in **bold** indicate a sterile product.

CARDIOVASCULAR & ANGIOGRAPHY PACKS



PACEMAKER PACK

CODE	DESCRIPTION		ITEMS/CASE
77870	Pacemaker Pack		8
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
-	-	-	1 Back Table Cover, Reinforced 140 x 190 cm
-	-	-	1 Bottom Drape 193 x 239 cm • CONTROL* Plus Fabric Reinforcement • Positioning Tape
-	-	-	1 Bar Drape 112 x 244 cm • CONTROL* Plus Fabric Reinforcement • Positioning Tape
-	-	-	1 Utility Drape with Tape 90 x 90 cm
-	-	-	1 Standard Absorbent Drape 91 x 91 cm
-	-	-	2 Absorbent Towels 38 x 56 cm
-	-	-	1 X-Ray Cover 71 cm

All codes in **bold** indicate a sterile product.

CARDIOVASCULAR DRAPE

CVARTS* HEAVY-DUTY SPLIT DRAPE, 274 X 315 X 203 CM

- Armboard Covers
- Split 15 x 213 cm
- Longitudinal Line-Control Flaps
- Heavy-Duty Fabric
- CONTROL* Plus Fabric Reinforcement 58 x 203 cm around Split
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89452	CVARTS* Heavy-Duty Split Drape, Sterile	6

CVARTS* BILATERAL UNDERLEG PERINEAL DRAPE, 213 X 203 CM




- Impervious Plastic
- Impervious Perineal Cover
- Umbilicus Cover
- CONTROL* Plus Fabric Reinforcement 122 x 76 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79453	CVARTS* Bilateral Underleg Perineal Drape, Non-Sterile	18





CVARTS* BILATERAL UNDERLEG PERINEAL DRAPE, 200 X 277 CM

- CONTROL* Fabric Reinforcement 127 x 76 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
K894530	CVARTS* Bilateral Underleg Perineal Drape, Sterile	18
K794530	CVARTS* Bilateral Underleg Perineal Drape, Non-Sterile	18



CARDIOVASCULAR & ANGIOGRAPHY DRAPES



CV INCISE DRAPE, 264 X 363 X 193 CM

- With Pouches
- CONTROL* Plus Fabric Reinforcement 122 x 53 cm
- Tube Holders

CODE	DESCRIPTION	ITEMS/CASE
86357	CV Incise Drape, Sterile	12
76357	CV Incise Drape, Non-Sterile	11
K763570	CV Incise Drape, Non-Sterile	12

CV SPLIT DRAPE, 224 X 320 X 295 CM

- Armboard Covers
- CONTROL* Plus Fabric Reinforcement 107 x 216 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
76378	CV Split Drape, Non-Sterile	9

FENESTRATED DRAPE, 68 X 68 CM

- SMS
- Circular Fenestration 18 cm, centered
- CONTROL* Fabric Reinforcement 30 x 30 cm

CODE	DESCRIPTION	ITEMS/CASE
K792015	Fenestrated Drape, Non-Sterile	180

PERINEAL DRAPE, 187 X 241 CM

- SMS
- CONTROL* Plus Fabric Reinforcement 99 x 141 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
76356	Perineal Drape, Non-Sterile	17

All codes in **bold** indicate a sterile product.

CARDIOVASCULAR & ANGIOGRAPHY DRAPES

DEFIBRILLATOR DRAPE I, 264 X 193 X 345 CM

- Armboard Covers
- Heavy-Duty fabric
- 30 x 51 cm Fenestration with Incise Film
- Two Side Pouches 30 x 66 cm
- One Pocket 38 x 38 cm
- Two Pockets 12 x 17 cm
- Tube Holders


CODE	DESCRIPTION	ITEMS/CASE
76354	Defibrillator Drape I, Non-Sterile	5

SCATTER RADIATION DRAPES

2 Levels of Protection:


- Level 1 - Attenuation: Greater than 75% at 90kVp and Lead Equivalency: Greater than 0.125mm (Lead Free)
- Level 2 - Attenuation: Greater than 90% at 90kVp and Lead Equivalency: Greater than 0.25mm (Lead Free)

CODE	DESCRIPTION	ITEMS/CASE
48958	HALYARD® Multi-Purpose Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 1), Sterile	45
48957	HALYARD® Multi-Purpose Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 1), Non-Sterile	50
49037	HALYARD® Multi-Purpose Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 2), Sterile	45
49036	HALYARD® Multi-Purpose Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 2), Non-Sterile	50
48986	HALYARD® Multi-Purpose Shield Drape W Fenestration, 30.4 cm x 41.9 cm (Attenuation Level: 1), Sterile	45
48959	HALYARD® Multi-Purpose Shield Drape W Fenestration, 30.4 cm x 41.9 cm (Attenuation Level: 1), Non-Sterile	50
49039	HALYARD® Multi-Purpose Shield Drape W Fenestration, 30.4 cm x 41.9 cm (Attenuation Level: 2), Sterile	45
49038	HALYARD® Multi-Purpose Shield Drape W Fenestration, 30.4 cm x 41.9 cm (Attenuation Level: 2), Non-Sterile	50
48988	HALYARD® Angiography Shield Drape, 33 cm x 40.6 cm (Attenuation Level: 1), Sterile	45
48987	HALYARD® Angiography Shield Drape, 33 cm x 40.6 cm (Attenuation Level: 1), Non-Sterile	50
49072	HALYARD® Angiography Shield Drape, 33 cm x 40.6 cm (Attenuation Level: 2), Sterile	45
49071	HALYARD® Angiography Shield Drape, 33 cm x 40.6 cm (Attenuation Level: 2), Non-Sterile	50



All codes in **bold** indicate a sterile product.

CARDIOVASCULAR & ANGIOGRAPHY DRAPES



48990	HALYARD* Subclavian Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 1), Sterile	45
48989	HALYARD* Subclavian Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 1), Non-Sterile	50
49074	HALYARD* Subclavian Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 2), Sterile	45
49073	HALYARD* Subclavian Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 2), Non-Sterile	50
49016	HALYARD* Jugular Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 1), Sterile	45
49000	HALYARD* Jugular Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 1), Non-Sterile	50
49076	HALYARD* Jugular Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 2), Sterile	45
49075	HALYARD* Jugular Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 2), Non-Sterile	50
49020	HALYARD* Biliary Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 1), Sterile	45
49017	HALYARD* Biliary Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 1), Non-Sterile	50
49078	HALYARD* Biliary Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 2), Sterile	45
49077	HALYARD* Biliary Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 2), Non-Sterile	50
49035	HALYARD* Peripheral Shield Drape, 30.4 cm x 86.3 cm (Attenuation Level: 1), Sterile	45
49021	HALYARD* Peripheral Shield Drape, 30.4 cm x 86.3 cm (Attenuation Level: 1), Non-Sterile	50
49080	HALYARD* Peripheral Shield Drape, 30.4 cm x 86.3 cm (Attenuation Level: 2), Sterile	45
49079	HALYARD* Peripheral Shield Drape, 30.4 cm x 86.3 cm (Attenuation Level: 2), Non-Sterile	50
49070	HALYARD* Biopsy Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 2), Sterile	45
49069	HALYARD* Biopsy Shield Drape, 30.4 cm x 41.9 cm (Attenuation Level: 2), Non-Sterile	50

All codes in **bold** indicate a sterile product.

CARDIOVASCULAR & ANGIOGRAPHY DRAPE

ANGIOGRAPHY DRAPE

FEMORAL ANGIOGRAPHY DRAPE, 193 X 279 CM

- SMS
- Two 12 cm Diameter Fenestrations
- CONTROL* Plus Fabric Reinforcement 81 x 122 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89714	Femoral Angiography Drape, Sterile	14
79714	Femoral Angiography Drape, Non-Sterile	160

HALYARD BASICS* FEMORAL ANGIOGRAPHY DRAPE, 239 X 350 CM, WITH SINGLE WINDOW




- SMS
- Two 12 cm Diameter Fenestrations
- 68 x 315 cm Window

CODE	DESCRIPTION	ITEMS/CASE
79582	HALYARD BASICS* Femoral Angiography Drape with Single Window, Non-Sterile	24

EXTRA-LONG FEMORAL ANGIOGRAPHY DRAPE, 220 X 342 CM, WITH DUAL WINDOWS


- SMS
- Two 12.7 cm Diameter Fenestrations (8.8 cm Diameter Incise)
- Two 69 x 342 cm Windows
- CONTROL* Plus Fabric Reinforcement 81 x 152 cm

CODE	DESCRIPTION	ITEMS/CASE
13296	Extra-Long Femoral Angiography Drape, Non-Sterile	180



All codes in **bold** indicate a sterile product.

CARDIOVASCULAR & ANGIOGRAPHY DRAPES



ACIST CVI® DRAPE, 287 X 363 CM

- SMS
- Enlarged Windows fit to cover the ACIST® Equipment
- Special Opening for passing ACIST® Manifold
- Two Femoral Fenestrations
- Two Radial Fenestrations
- CONTROL® Fabric Reinforcement
- Positioning Tabs on the Windows

CODE	DESCRIPTION	ITEMS/CASE
K797065	ACIST CVi® Drape, Non-Sterile	22

FEMORAL ANGIOGRAPHY DRAPE, 211 X 292 CM, WITH SINGLE WINDOW

- SMS
- Two 12 cm Diameter Fenestrations
- 68 x 292 cm Window
- CONTROL® Plus Fabric Reinforcement 81 x 122 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89715	Femoral Angiography Drape with Single Window, Sterile	16
79715	Femoral Angiography Drape with Single Window, Non-Sterile	200

FEMORAL ANGIOGRAPHY DRAPE, 225 X 292 CM, WITH DUAL WINDOWS

- SMS
- Two 12 cm Diameter Fenestrations
- Dual 68 x 292 cm Windows
- CONTROL® Plus Fabric Reinforcement 81 x 122 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89716	Femoral Angiography Drape with Dual Windows, Sterile	18
79716	Femoral Angiography Drape with Dual Windows, Non-Sterile	200

BRACHIAL ANGIOGRAPHY DRAPE, 112 X 79 CM

- SMS
- 8 cm Circular Fenestration
- CONTROL® Plus Fabric Reinforcement, 27 x 61 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89471	Brachial Angiography Drape, Sterile	50

All codes in **bold** indicate a sterile product.

CARDIOVASCULAR & ANGIOGRAPHY DRAPE

IMPERVIOUS EXTRA-LONG FEMORAL ANGIOGRAPHY DRAPE, 221 X 343 CM, WITH DUAL WINDOWS

- SMS
- Two 12 cm Diameter Fenestrations
- Fabric Reinforcement 343 x 81 cm
- Two Dual Plastic Windows (343 x 70 cm)
- Extra Long Film Based Drape
- Incise with Release Liner

CODE	DESCRIPTION	ITEMS/CASE
79749	Impervious Extra-Long Femoral Angiography Drape, Non-Sterile	144

FEMORAL ANGIOGRAPHY DRAPE, 193 X 330 CM

- SMS
- Two 12 cm Diameter Fenestrations
- CONTROL* Plus Fabric Reinforcement 81 x 122 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
79063	Femoral Angiography Drape, Non-Sterile	160

FEMORAL ANGIOGRAPHY DRAPE, 211 X 330 CM, WITH SINGLE WINDOW





- SMS
- 68 x 330 cm Window
- Two 12 cm Diameter Fenestrations
- CONTROL* Plus Fabric Reinforcement 81 x 122 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89048	Femoral Angiography Drape with Single Window, Sterile	16
79109	Femoral Angiography Drape with Single Window, Non-Sterile	200

PACEMAKER DRAPE WITH FENESTRATED INCISE, 231 X 343 CM


- SMS
- Absorbent CONTROL* Plus fabric reinforcement with pockets 25 cm deep
- Two 7 x 12 cm fenestrations
- Tube Holders
- Dual Clear Film Poly Windows

CODE	DESCRIPTION	ITEMS/CASE
29229	Pacemaker Drape with Fenestrated Incise, Sterile	12
29231	Pacemaker Drape with Fenestrated Incise, Non-Sterile	12



All codes in **bold** indicate a sterile product.

CARDIOVASCULAR & ANGIOGRAPHY DRAPES



PACEMAKER DRAPE, 221 X 343 CM

- SMS
- CONTROL* Plus Fabric Reinforcement 81 x 89 cm
- Two 20 x 15 cm Fenestrations with Incise 25 x 20 cm
- Dual Clear Film Poly Windows 56 x 343 cm

CODE	DESCRIPTION	ITEMS/CASE
89368	Pacemaker Drape, Sterile	11
79368	Pacemaker Drape, Non-Sterile	12

FEMORAL ANGIOGRAPHY/IJ ACCESS COMBINATION DRAPE, 200 X 360 CM

- SMS
- Absorbent CONTROL* Plus fabric reinforcement
- Two 9 cm adhesive circular fenestrations
- Two 18 x 18 cm trapezoidal fenestrations
- Clear Side Poly Windows

CODE	DESCRIPTION	ITEMS/CASE
29225	Femoral Angiography/IJ Access Combination Drape, Sterile	15
29226	Femoral Angiography/IJ Access Combination Drape, Non-Sterile	15

FEMORAL ANGIOGRAPHY DRAPE, 225 X 330 CM, WITH DUAL WINDOWS

- SMS
- Dual 68 x 330 cm Windows
- Two 12 cm Diameter Fenestrations
- CONTROL* Plus Fabric Reinforcement 81 x 122 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
77776	Femoral Angiography Drape with Dual Windows, Sterile	14
76386	Femoral Angiography Drape with Dual Windows, Non-Sterile	200

RADIAL FEMORAL ANGIOGRAPHY DRAPE, 244 X 351 CM

- SMS
- Dual Windows 63 x 350 cm
- CONTROL* Plus Fabric Reinforcement 122 x 154 cm
- Two Oval Radial Fenestrations 12 x 7.6 cm
- One Round Femoral Fenestration, Diameter 12.7 cm

CODE	DESCRIPTION	ITEMS/CASE
76681	Radial Femoral Angiography Drape, Non-Sterile	150

All codes in **bold** indicate a sterile product.

CARDIOVASCULAR & ANGIOGRAPHY DRAPE

PAEDIATRIC FEMORAL ANGIOGRAPHY DRAPE, 193 X 310 CM, WITH SINGLE WINDOW

- SMS
- 64 x 215 cm Window
- Two 7.6 cm Diameter Fenestrations
- CONTROL* Plus Fabric Reinforcement 122 x 61 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
77382	Paediatric Femoral Angiography Drape with Single Window, Sterile	25
76382	Paediatric Femoral Angiography Drape with Single Window, Non-Sterile	14

HALYARD BASICS* RADIAL FEMORAL ANGIOGRAPHY DRAPE, 240 X 380 CM, WITH DUAL WINDOWS




- SMS
- Dual Windows 63 x 380 cm
- Two Oval Radial Fenestrations 12 x 7.6 cm
- Two Round Femoral Fenestrations, Diameter 12.7 cm

CODE	DESCRIPTION	ITEMS/CASE
89583	HALYARD BASICS* Radial Femoral Angiography Drape with Dual Windows, Sterile	12
79583	HALYARD BASICS* Radial Femoral Angiography Drape with Dual Windows, Non-Sterile	14

IMPERVIOUS RADIAL FEMORAL ANGIOGRAPHY DRAPE, WITH DUAL WINDOWS


- SMS
- Dual Windows 72 x 381 cm
- Two Round Fenestrations, Diameter 13.3 cm
- Two Oval Fenestrations 7.6 x 12.7 cm
- CONTROL* Fabric Reinforcement
- Extra CONTROL* Fabric Reinforcement at the Bottom
- Non-Sterile

CODE	SIZE	REINFORCEMENT		ITEMS/CASE
K797067	264 x 381 cm	137 x 180 cm	Extra 120 x 161 cm	20
K797160	234 x 381 cm	157 x 122 cm	Extra 91 x 220 cm	16



All codes in **bold** indicate a sterile product.

CARDIOVASCULAR & ANGIOGRAPHY DRAPES



FEMORAL ANGIOGRAPHY DRAPE 230 CM X 360 CM

CODE	DESCRIPTION	ITEMS/CASE
K797166	Femoral Angiography Drape, Non-Sterile	12

RADIAL FEMORAL ANGIOGRAPHY DRAPE 264X381 CM

CODE	DESCRIPTION	ITEMS/CASE
K797169	Radial Femoral Angiography Drape, Non-Sterile	20

IMPERVIOUS RADIAL FEMORAL ANGIOGRAPHY DRAPE 292 CM X 368 CM

CODE	DESCRIPTION	ITEMS/CASE
K797165	Impervious Radial Femoral Angiography Drape, Non-Sterile	16

IMPERVIOUS RADIAL FEMORAL ANGIOGRAPHY DRAPE, 269 X 370 CM, WITH DUAL WINDOWS

- SMS/Film Laminate
- Two Plastic Dual Windows 73.7 x 274 cm & 127 x 274 cm
- Reinforcement with Fenestrations size 121.9 x 157.5 cm
- Two Oval Fenestrations 11.7 x 8.9 cm
- Two Round Fenestrations, Diameter 12.7 cm

CODE	DESCRIPTION	ITEMS/CASE
K797162	Impervious Radial Femoral Angiography Drape with Dual Windows, Non-Sterile	20

RADIAL FEMORAL ANGIOGRAPHY DRAPE, 259 X 376 CM, WITH DUAL WINDOWS

- SMS
- CONTROL* Plus Fabric Reinforcement 165 x 122 cm & 128 x 122 cm & 30.5 x 122 cm
- Two Oval Radial Fenestrations 8.9 x 12.1 cm
- Two Round Femoral Fenestrations, Diameter 12.7 cm
- Dual Clear Film Poly Windows
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
K797064	Radial Femoral Angiography Drape with Dual Windows, Non-Sterile	20
K795830	Radial Femoral Angiography Drape, 4 holes	14

All codes in **bold** indicate a sterile product.

CATHETER PLACEMENT

CATHETER PLACEMENT DRAPES


CVC DRAPE, SINGLE FENESTRATION, 172 X 279 CM

- SMS
- Absorbent CONTROL® Plus Fabric Reinforcement
- 10 cm Fenestration
- Incise Window

CODE	DESCRIPTION	ITEMS/CASE
29149	CVC Drape, Single Fenestration, Sterile	20
29148	CVC Drape, Single Fenestration, Non-Sterile	27

DRAINAGE/CVC DRAPE 230X290 CM


CODE	DESCRIPTION	ITEMS/CASE
K797167	Drainage/CVC Drape, Sterile	22



All codes in **bold** indicate a sterile product.

PAEDIATRICS

PAEDIATRICS PACKS



UNIVERSAL PACK, SMALL

CODE	DESCRIPTION		ITEMS/CASE
88762	Universal Pack, Small, incl. Mayo Stand Cover		11
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 96 x 96 cm
89603		79603	1 Mayo Stand Cover 80 x 145 cm
42222		42428	1 Back Table Cover, Reinforced 152 x 191 cm
89053		79053	1 Tape Strip 10 x 55 cm
89701		79720	1 Absorbent Towel 38 x 56 cm
	-	-	1 Suture Bag
	-	-	2 Side Drapes 90 x 75 cm: • Positioning Tape, 2.5 cm width • CONTROL* Fabric Reinforcement
	-	79214	1 Bottom Drape 149 x 180 cm: • Positioning Tape, 2.5 cm width • CONTROL* Fabric Reinforcement
	-	79215	1 Top Drape 149 x 180 cm: • Positioning Tape, 2.5 cm width • CONTROL* Fabric Reinforcement


UNIVERSAL PACK, SMALL

CODE	DESCRIPTION		ITEMS/CASE
88763	Universal Pack, Small, excl. Mayo Stand Cover		11
	STERILE	NON-STERILE	CONTENTS
	-	-	1 Sterilisation Wrap 96 x 96 cm
42222		42428	1 Back Table Cover, Reinforced 152 x 191 cm
89053		79053	1 Tape Strip 10 x 55 cm
89701		79720	1 Absorbent Towel 38 x 56 cm
	-	-	2 Side Drapes 90 x 75 cm: • Positioning Tape, 2.5 cm width • CONTROL* Fabric Reinforcement
	-	79214	1 Bottom Drape 149 x 180 cm: • Positioning Tape, 2.5 cm width • CONTROL* Fabric Reinforcement
	-	79215	1 Top Drape 149 x 180 cm: • Positioning Tape, 2.5 cm width • CONTROL* Fabric Reinforcement

All codes in **bold** indicate a sterile product.


UNIVERSAL PACK

CODE	DESCRIPTION		ITEMS/CASE
88100	Universal Pack, incl. Mayo Stand Cover		12
STERILE	NON-STERILE	CONTENTS	
42224	42451	1 Back Table Cover, Reinforced 152 x 229 cm	
89603	79603	1 Mayo Stand Cover 80 x 145 cm	
89053	79053	1 Tape Strip 10 x 55 cm	
89701	79720	2 Absorbent Towels 38 x 56 cm	
-	79012	2 Side Drapes 75 x 100 cm: • Positioning Tape • CONTROL ⁺ Plus Fabric Reinforcement	
-	79011	1 Bottom Drape 173 x 180 cm: • Positioning Tape • CONTROL ⁺ Plus Fabric Reinforcement	
-	79010	1 Top Drape 148 x 240 cm: • Positioning Tape • CONTROL ⁺ Plus Fabric Reinforcement	




UNIVERSAL PACK

CODE	DESCRIPTION		ITEMS/CASE
88103	Universal Pack, excl. Mayo Stand Cover		12
STERILE	NON-STERILE	CONTENTS	
-	-	1 Sterilisation Wrap 96 x 96 cm	
42224	42451	1 Back Table Cover, Reinforced 152 x 229 cm	
89053	79053	1 Tape Strip 10 x 55 cm	
89701	79720	2 Absorbent Towels 38 x 56 cm	
-	79012	2 Side Drapes 75 x 100 cm: • Positioning Tape • CONTROL ⁺ Plus Fabric Reinforcement	
-	79011	1 Bottom Drape 173 x 180 cm: • Positioning Tape • CONTROL ⁺ Plus Fabric Reinforcement	
-	79010	1 Top Drape 148 x 240 cm: • Positioning Tape • CONTROL ⁺ Plus Fabric Reinforcement	



All codes in **bold** indicate a sterile product.

PAEDIATRICS PACKS



U-BAR* PACK VI

CODE	DESCRIPTION		ITEMS/CASE
77794	U-BAR* Pack VI		9
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 96 x 96 cm
-	-	42451	1 Back Table Cover, Reinforced 152 x 229 cm
89603	-	79603	1 Mayo Stand Cover 80 x 145 cm
89701	-	79720	1 Absorbent Towel 38 x 56 cm
-	-	-	1 Suture Bag
-	-	76361	1 U-Drape 193 x 305 cm: • Split 10 x 101 cm • Positioning Tape • CONTROL* Plus Fabric Reinforcement • Tube Holders
-	-	76363	1 Bar Drape 157 x 269 cm: • Armboard Covers • Positioning Tape • CONTROL* Plus Fabric Reinforcement • Tube Holders

U-BAR* PACK

CODE	DESCRIPTION		ITEMS/CASE
77853	U-BAR* Pack		10
	STERILE	NON-STERILE	CONTENTS
-	-	-	1 Sterilisation Wrap 89 x 89 cm
-	-	42451	1 Back Table Cover, Reinforced 152 x 229 cm
89701	-	79720	1 Absorbent Towel 38 x 56 cm
-	-	-	1 Suture Bag
-	-	76361	1 U-Drape 193 x 305 cm: • Split 10 x 101 cm • Positioning Tape • CONTROL* Plus Fabric Reinforcement • Tube Holders
-	-	76363	1 Bar Drape 157 x 269 cm: • Armboard Covers • Positioning Tape • CONTROL* Plus Fabric Reinforcement • Tube Holders
89053	-	79053	1 Tape Strip 10 x 55 cm

All codes in **bold** indicate a sterile product.

SMALL PROCEDURE PACK II

CODE	DESCRIPTION		ITEMS/CASE
77717	Small Procedure Pack II		10
STERILE	NON-STERILE	CONTENTS	
-	-	Sterilisation Wrap 89 x 89 cm	
42222	42428	1 Back Table Cover, Reinforced 152 x 191 cm	
-	-	1 Fenestrated Drape 211 x 241 cm: • Circular Fenestration 12 cm centered • CONTROL* Plus Fabric Reinforcement	
89701	79720	2 Absorbent Towels 38 x 56 cm	
89053	79053	1 Tape Strip 10 x 55 cm	


HALYARD BASICS* BASIC PACK

CODE	DESCRIPTION		ITEMS/CASE
89682	HALYARD BASICS* Basic Pack		10
STERILE	NON-STERILE	CONTENTS	
89632	79741	1 Back Table Cover, Zone Reinforced, 100 x 150 cm: • Spunbond	
89679	79783	4 Utility Drapes with Tape, 90 x 90 cm: • Spunbond Film	

SIDE DRAPE WITH ADHESIVE, 120 X 100 CM

- SMS
- Adhesive size 2,5cm x 120 cm

CODE	DESCRIPTION	ITEMS/CASE
K795805	Side Drape with Adhesive	70



All codes in **bold** indicate a sterile product.

PAEDIATRICS DRAPES

FENESTRATED DRAPES

MINOR PROCEDURE FENESTRATED DRAPE, 76 X 76 CM

- Fenestration 3 x 5 cm
- CONTROL* Fabric Reinforcement 12 x 20 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/ CASE	ITEMS/ DISPENSER	DISPENSERS/ CASE
89206	Minor Procedure Fenestrated Drape, Sterile	100	25 x 4	4
79206	Minor Procedure Fenestrated Drape, Non-Sterile	250 -	-	

MINOR PROCEDURE FENESTRATED DRAPE, 122 X 152 CM


- Fenestration (Diameter 7 cm)
- CONTROL* Fabric Reinforcement 31 x 31 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89207	Minor Procedure Fenestrated Drape, Sterile	40
79207	Minor Procedure Fenestrated Drape, Non-Sterile	60

MINOR PROCEDURE FENESTRATED DRAPE, 122 X 152 CM

- Fenestration 12 x 15 cm
- CONTROL* Fabric Reinforcement 41 x 31 cm
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89208	Minor Procedure Fenestrated Drape, Sterile	40
79208	Minor Procedure Fenestrated Drape, Non-Sterile	60



All codes in **bold** indicate a sterile product.

BREAST/CHEST DRAPES

BREAST DRAPE, 254 X 315 X 193 CM

- Fenestration 15 x 15 cm
- Hook and Loop Tube and Line Control
- CONTROL* Fabric Reinforcement 51 x 71 cm
- Armboard Covers

CODE	DESCRIPTION	ITEMS/CASE
89236	Breast Drape, Sterile	10
79236	Breast Drape, Non-Sterile	20

ANGIOGRAPHY DRAPES

PAEDIATRIC FEMORAL ANGIOGRAPHY DRAPE, 193 X 310 CM, WITH SINGLE WINDOW

- 64 x 215 cm Window
- CONTROL* Plus Fabric Reinforcement 122 x 61 cm
- Two 7.6 cm Diameter Fenestrations
- Positioning Tape



CODE	DESCRIPTION	ITEMS/CASE
77382	Paediatric Femoral Angiography Drape with Single Window, Sterile	25
76382	Paediatric Femoral Angiography Drape with Single Window, Non-Sterile	14

BABY BLANKETS

BABY BLANKET, 80 X 130 CM

- Absorbent

CODE	DESCRIPTION	ITEMS/CASE
76394	Baby Blanket, Non-Sterile	141




EQUIPMENT COVERS

EQUIPMENT COVERS

MAYO STAND COVERS


MAYO STAND COVERS



CODE STERILE	NON-STERILE	DESCRIPTION	STERILE	ITEMS/CASE NON-STERILE
89601	79601	58 x 137 cm	54	150
-	79602	76 x 145 cm	-	120
89603	79603	80 x 145 cm	44	120
-	K796019	80 x 170 cm	-	60
47348	78657	HALYARD BASICS® Mayo Stand Cover, 80 x 144 cm: • Zone-Reinforced • Spunbond	40	90
14007	76358	Reinforced Mayo Stand Cover, 76 x 145 cm: • CONTROL® Plus Fabric Reinforcement 61 x 76 cm	44	68
77883	-	Mayo Cover Reinforcement 60 x 80 cm: • CONTROL® Plus Fabric Reinforcement	50	-
77452	-	Large Mayo Stand Cover, 121 x 182 cm • CONTROL® Plus Fabric Reinforcement 119 x 88 cm • Cuff 18 cm	15	-
-	K796011	Reinforced Mayo Stand Cover, 80 x 145 cm • CONTROL® Fabric Reinforcement	-	60
-	K796012	Reinforced Mayo Stand Cover, 95 x 195 cm • CONTROL® Fabric Reinforcement	-	54

STOOL COVERS

STOOL COVER



CODE STERILE	NON-STERILE	DESCRIPTION	STERILE	ITEMS/CASE NON-STERILE
K103600	K103502	Stool Cover, 73 x 73 cm, Non-Sterile • SMS	140	150

All codes in **bold** indicate a sterile product.



BACK TABLE COVERS

BACK TABLE COVERS, ZONE-REINFORCED

CODE STERILE	NON-STERILE	DESCRIPTION	STERILE	ITEMS/CASE NON-STERILE
89632	79741	HALYARD BASICS® Back Table Cover, Zone-Reinforced, 100 x 150 cm: • Spunbond	50	70
88673	78669	HALYARD BASICS® Back Table Cover, Zone-Reinforced, 152 x 191 cm: • Spunbond	20	50

BACK TABLE COVERS, REINFORCED

CODE STERILE	NON-STERILE	DESCRIPTION	STERILE	ITEMS/CASE NON-STERILE
42217		Back Table Cover, Reinforced, 112 x 229 cm	28	
-	42447	Back Table Cover, Reinforced, 112 x 198 cm	-	800
-	42446	Back Table Cover, Reinforced, 112 x 112 cm	-	1280
-	42228	Back Table Cover, Reinforced, 178 x 279 cm	36	320
-	42449	Back Table Cover, Reinforced, 127 x 229 cm	-	560
42222	42428	Back Table Cover, Reinforced, 152 x 191 cm	28	540
-	42481	Back Table Cover, Reinforced, 152 x 152 cm	-	660
-	42225	Back Table Cover, Reinforced, 152 x 229 cm	-	480
-	42409	Back Table Cover, Reinforced, Fan Fold, 178 x 279 cm	36	200
-	42366	Back Table Cover, Reinforced, Heavy Duty, Fan Fold, 152 x 229 cm	-	200
-	42367	Back Table Cover, Reinforced, Heavy Duty, 178 x 279 cm	-	160
-	42302	Back Table Cover, Reinforced, Fan Fold, 112 x 224 cm	-	60
42345	42365	Back Table Cover, Reinforced, 140 x 190 cm	24	420
-	42480	Back Table Cover, Reinforced, 150 x 150 cm	-	720
42224	42451	Back Table Cover, Reinforced, 152 x 229 cm	36	480
42312	42313	Back Table Cover, Reinforced, 191 x 201 cm	20	32
-	76372	Back Table Cover 190 x 240 cm with CONTROL® Plus Fabric Reinforcement 102 x 240 cm	-	8
42309	-	Back Table Cover, Reinforced, Heavy Duty, 177 x 279 cm	15	-
-	47520	Back Table Cover, Standard, 152 x 191 cm	-	480



All codes in **bold** indicate a sterile product.

ACCESSORIES

ACCESSORIES

UNDER BUTTOCKS DRAPES

UNDER BUTTOCKS DRAPE, PLASTIC, 75 X 112 CM

CODE	DESCRIPTION	ITEMS/CASE
89411	Under Buttocks Drape, Plastic, Sterile	48

UNDER BUTTOCKS DRAPE, ZONE REINFORCED, MEDIUM, 97 X 112 CM

- Film Fabric Reinforcement

CODE	DESCRIPTION	ITEMS/CASE
79411	Under Buttocks Drape, Zone Reinforced, Medium, Non-Sterile	115

UNDER BUTTOCKS DRAPE, PLASTIC REINFORCED, LARGE, 102 X 117 CM

CODE	DESCRIPTION	ITEMS/CASE
89413	Under Buttocks Drape, Plastic Reinforced, Large, Sterile	50

UNDER BUTTOCKS DRAPE, 102 X 112 CM

- With Fluid Collection Pouch, Screen and Suction Port
- Adhesive Pouch Closure

CODE	DESCRIPTION	ITEMS/CASE
89415	Under Buttocks Drape, Sterile	40
79415	Under Buttocks Drape, Non-Sterile	100

HALYARD BASICS* UNDER BUTTOCKS DRAPE, 102 X 112 CM






- With Graduated Fluid Collection Pouch
- Spunbond

CODE	DESCRIPTION	ITEMS/CASE
89584	HALYARD BASICS* Under Buttocks Drape, Sterile	65

CLEAR DRAPES

CLEAR TOWEL DRAPE 45 CM X 30 CM

CODE	DESCRIPTION	ITEMS/CASE
49046	Clear Towel Drape, Sterile	40
49056	Clear Towel Drape, Non-Sterile	200



All codes in **bold** indicate a sterile product.

CLEAR TOWEL DRAPE, 45 CM X 60 CM

CODE	DESCRIPTION	ITEMS/CASE
49047	Clear Towel Drape, Sterile	40
49057	Clear Towel Drape, Non-Sterile	300

CLEAR U-DRAPE WITH APERTURE 120 CM X 130 CM

CODE	DESCRIPTION	ITEMS/CASE
49048	Clear U-Drape with Aperture, Sterile	40
49058	Clear U-Drape with Aperture, Non-Sterile	200

CLEAR U-DRAPE WITH APERTURE 120 CM X 180 CM

CODE	DESCRIPTION	ITEMS/CASE
49049	Clear U-Drape with Aperture, Sterile	40
49059	Clear U-Drape with Aperture, Non-Sterile	200

CLEAR INCISE DRAPE 40 CM X 40 CM

CODE	DESCRIPTION	ITEMS/CASE
49050	Clear Incise Drape, Sterile	80
49060	Clear Incise Drape, Non-Sterile	200

PERINEAL COVER DRAPES

HALYARD BASICS* PERINEAL COVER DRAPE WITH TAPE, 15 X 40 CM



- Spunbond Film
- Positioning Tape

CODE	DESCRIPTION	ITEMS/CASE
89675	HALYARD BASICS* Perineal Cover Drape with Tape, Sterile	100
79697	HALYARD BASICS* Perineal Cover Drape with Tape, Non-Sterile	500

IMPERVIOUS STOCKINETTES

IMPERVIOUS STOCKINETTE, MEDIUM, 23 X 112 CM

CODE	DESCRIPTION	ITEMS/CASE
89811	Impervious Stockinette, Medium, Sterile	36
48759	Impervious Stockinette, Medium, Non-Sterile	90




All codes in **bold** indicate a sterile product.



ACCESSORIES



76214



44591



IMPERVIOUS STOCKINETTE, LARGE, 30 X 122 CM

CODE	DESCRIPTION	ITEMS/CASE
89821	Impervious Stockinette, Large, Sterile	25
48760	Impervious Stockinette, Large, Non-Sterile	56

ARMBOARD COVERS

ARMBOARD COVER

- Cuff 15 cm

CODE	DESCRIPTION	ITEMS/CASE
76214	Armboard Cover, 61 x 72 cm, Sterile	164 (2 covers/pouch)
76300	Armboard Cover, 61 x 72 cm, Non-Sterile	166
44591	Armboard Cover, 35.5 x 80 cm, Non-Sterile	3000 (100 covers/poly bag)

LEGGINGS

LEGGINGS, 79 X 122 CM, WITH 15 CM CUFF

CODE	DESCRIPTION	ITEMS/PACKAGE	ITEMS/CASE
89408	Leggings, Sterile Pair	2	52 (2x26)
79408	Leggings, Non-Sterile	-	70

TOWELS

ABSORBENT TOWELS, 38 X 56 CM

CODE	DESCRIPTION	ITEMS/PACKAGE	ITEMS/CASE
89701	Absorbent Towel, Sterile	2	200 (2x100)
79720	Absorbent Towel, Non-Sterile	-	1.000

OR TOWELS

CODE	DESCRIPTION	SIZE	ITEMS/PACKAGE	ITEMS/CASE
AR-28502	Blue OR towel XR, Sterile	45 x 45 cm	2	80
AR-28503	Blue OR towel XR, Sterile	43 x 69 cm	2	80
AR-28504	Blue OR towel XR, Sterile	70 x 90 cm	1	40
AR-63001NS	Blue OR towel, Non-Sterile	66 x 43 cm	-	100

All codes in **bold** indicate a sterile product.

TAPE STRIPS

TAPE STRIPS, 10 X 55 CM

CODE	DESCRIPTION	ITEMS/PACKAGE	ITEMS/CASE
89053	Tape Strips, Sterile	2	600
79053	Tape Strips, Non-Sterile	-	1.000

SLEEVES

SLEEVE, 15 X 59 CM

CODE	DESCRIPTION	ITEMS/PACKAGE	ITEMS/CASE
89791	Sleeve, Sterile	1	100

POUCHES

DIATHERMY POUCH, 38 X 41 CM

- 2 compartments

CODE	DESCRIPTION	ITEMS/PACKAGE	ITEMS/CASE
89945	Diathermy Pouch, Sterile	1	150
79945	Diathermy Pouch, Non-Sterile	-	220

DIATHERMY POUCH, SMS HEAVY DUTY, 27 X 29 CM

- 2 Tape Strips at the Back


CODE	DESCRIPTION	ITEMS/PACKAGE	ITEMS/CASE
76270	Diathermy Pouch, SMS Heavy Duty, Sterile	1	200

CLEAR U-POUCH 89 CM X 76 CM



CODE	DESCRIPTION	ITEMS/CASE
49041	Clear U-pouch, Sterile	20
49051	Clear U-pouch, Non-Sterile	50

CLEAR IRRIGATION POUCH 50 CM X 60 CM


CODE	DESCRIPTION	ITEMS/CASE
49042	Irrigation Pouch, Clear, Sterile	40
49052	Clear Irrigation Pouch, Non-Sterile	200




10 cm x 55 cm



ACCESSORIES



76334



76340

CLEAR INSTRUMENT POUCH 18 CM X 30 CM

CODE	DESCRIPTION	ITEMS/CASE
49053	Clear Instrument Pouch, Non-Sterile	200

CLEAR INSTRUMENT POUCH 33 CM X 38 CM

CODE	DESCRIPTION	ITEMS/CASE
49045	Clear Instrument Pouch, Sterile	200
49055	Clear Instrument Pouch, Non-Sterile	336

CLEAR INSTRUMENT POUCH 25 CM X 47 CM

CODE	DESCRIPTION	ITEMS/CASE
49044	Clear Instrument Pouch, Sterile	40
49054	Clear Instrument Pouch, Non-Sterile	200

CLEAR INSTRUMENT POUCH 50 CM X 60 CM

CODE	DESCRIPTION	ITEMS/CASE
49043	Clear Instrument Pouch, Sterile	40

FOOT COVERS

FOOT COVERS

CODE	DESCRIPTION	ITEMS/PACKAGE	ITEMS/CASE
76334	Foot Cover, Non-Sterile	1	300
76340	Foot Cover, Non-Sterile	1	400

PITCHER

PITCHER GRADUATED 1200 ML

CODE	DESCRIPTION	ITEMS/CASE
AR-1200	Pitcher graduated, Non-Sterile	70

All codes in **bold** indicate a sterile product.

ROOM TURNOVER KITS

ROOM TURNOVER KIT I

CODE	DESCRIPTION		ITEMS/CASE
68500	Room Turnover Kit I		40
	STERILE	NON-STERILE	CONTENTS
	-	68474	1 Absorbent Table Cover, 101.6 cm x 228.6 cm, Grippy 1 Lift Sheet, 101.6 cm x 152.4 cm

ROOM TURNOVER KIT II









CODE	DESCRIPTION		ITEMS/CASE
68501	Room Turnover Kit II		36
	STERILE	NON-STERILE	CONTENTS
	-	68474	1 Absorbent Table Cover, 101.6 cm x 228.6 cm, Grippy 1 Lift Sheet, 101.6 cm x 152.4 cm 2 Armboard Cover, 33 cm x 76.2 cm with Strap

ROOM TURNOVER KIT III

CODE	DESCRIPTION		ITEMS/CASE
68502	Room Turnover Kit III		36
	STERILE	NON-STERILE	CONTENTS
	-	68474	1 Absorbent Table Cover, 101.6 cm x 228.6 cm, Grippy 1 Lift Sheet, 101.6 cm x 152.4 cm 2 Armboard Cover, 33 cm x 76.2 cm with Strap 1 Head Rest Cover, 35.5 cm x 38.1 cm



GRIPPY ABSORBENT TABLE COVER

CODE	DESCRIPTION	ITEMS/CASE
68474	Absorbent Table Cover, 101.6 cm x 228.6 cm, Grippy Finish, Non-Sterile	20



HALYARD BASICS*

HALYARD BASICS* SURGICAL GOWNS



HALYARD BASICS* FABRIC-REINFORCED SURGICAL GOWNS

- SMS Fabric
- Fabric-Reinforced Chest and Sleeves
- Overlapping Back
- Raglan Sleeves
- Low Linting Fabric
- Colour Coded Neck Binding: Green
- Sterile Gown packaged with 2 Towels

CODE		SIZE	LENGTH*	ITEMS/CASE	
STERILE	NON-STERILE			STERILE	NON-STERILE
13815	-	X-Large, XX-Long	175 cm	32	-

HALYARD BASICS* SURGICAL GOWNS

- SMS Fabric
- Overlapping Back
- Raglan Sleeves
- Low Linting Fabric
- Colour Coded Neck Binding: White
- Sterile Gown packaged with 2 Towels
- Meet EN13795 Standard Performance

CODE	SIZE	LENGTH*	ITEMS/CASE		
STERILE	NON-STERILE		STERILE NON-STERILE		
13337	-	X-Large, XX-Long	175 cm	40	
11783	33916	Large - Towel in Overwrap (handi-bin)	125 cm	40	440
11784	33917	X-Large - Towel in Overwrap (handi-bin)	131 cm	40	400
99294	32109	Large, X-Long - Towel in Overwrap (handi-bin)	145 cm	40	460
99295	32113	X-Large, X-Long - Towel in Overwrap (handi-bin)	150 cm	40	380
-	79284	Large (handi-bin)	125 cm	-	640
-	79294	Large, X-Long (handi-bin)	145 cm	-	560
-	79285	X-Large (handi-bin)	130 cm	-	560
-	79295	X-Large, X-Long (handi-bin)	150 cm	-	560

○ Non-Reinforced Fabric

◐ Fabric-Reinforced

● High Performance Area

All codes in **bold** indicate a sterile product.


EVOLUTION* 4

EVOLUTION* 4 SURGICAL GOWNS

EVOLUTION* 4 FABRIC-REINFORCED SURGICAL GOWNS

- EVOLUTION* 4 SMS Fabric
- Fabric-Reinforced Chest and Sleeves
- Raglan Sleeves
- Low Linting Fabric
- Colour Coded Neck Binding: Green
- Adjustable Neckline - Hook and Loop
- Sterile Gown packaged with 2 Towels
- Meet EN13795 Standard Performance

CODE	SIZE		LENGTH ⁺	ITEMS/CASE	
	STERILE	NON-STERILE		STERILE	NON-STERILE
90112	78911	Large	123 cm	32	36
90113	78913	Large, X-Long	143 cm	30	54
90142	78921	X-Large	127 cm	30	36
90143	78923	X-Large, X-Long	148 cm	32	48



⁺ length back (neck to bottom)



Non-Reinforced Fabric





Fabric-Reinforced



High Performance Area

All codes in **bold** indicate a sterile product.

EVOLUTION* 4 SURGICAL GOWNS



EVOLUTION* 4 STANDARD SURGICAL GOWNS SET-IN SLEEVES

- EVOLUTION* 4 SMS Fabric
- Set-In-Sleeves
- Low Linting Fabric
- Colour Coded Neck Binding: Yellow
- Adjustable Neckline – Hook and Loop
- Sterile Gown packaged with 2 Towels
- Meet EN13795 Standard Performance

CODE STERILE	NON-STERILE	SIZE	LENGTH*	STERILE	ITEMS/CASE NON-STERILE
90016	72978	Large	117 cm	36	640
90046	72988	X-Large	127 cm	32	640
90019	70391	Large, X-Long	132 cm	34	66
90049	70392	X-Large, X-Long	137 cm	34	60
-	79419	Large – Towel in Overwrap	111 cm	-	512
-	79434	X-Large – Towel in Overwrap	121 cm	-	448

EVOLUTION* 4 STANDARD SURGICAL GOWNS RAGLAN SLEEVES

CODE STERILE	NON-STERILE	SIZE	LENGTH*	STERILE	ITEMS/CASE NON-STERILE
90002	-	Small	115 cm	36	-
90012	79970	Large	123 cm	36	40
90042	79980	X-Large	127 cm	34	36
90013	79973	Large, X-Long	143 cm	34	60
90043	79983	X-Large, X-Long	148 cm	34	48

* length back (neck to bottom)



Non-Reinforced Fabric



Fabric-Reinforced



High Performance Area

All codes in **bold** indicate a sterile product.


ULTRA

ULTRA SURGICAL GOWNS

ULTRA FABRIC-REINFORCED SURGICAL GOWNS WITH NECK PROTECTION & THUMB LOOPS

- SMS Fabric
- Extra Wide Fabric Reinforcements
- Raglan Sleeves
- Low Linting Fabric
- Colour Coded Neck Binding: Green
- Adjustable Neckline – Hook and Loop
- Neck Protection
- Thumb Loops
- Sterile Gown packaged with 1 Towel
- Meet EN13795 Standard Performance

CODE STERILE	NON-STERILE	SIZE	LENGTH ⁺	ITEMS/CASE STERILE	NON-STERILE
90363	77337	Large, X-Long	143 cm	30	48
90373	77336	X-Large, X-Long	148 cm	28	36



⁺ length back (neck to bottom)



Non-Reinforced Fabric



Fabric-Reinforced



High Performance Area

All codes in **bold** indicate a sterile product.

ULTRA SURGICAL GOWNS

ULTRA FABRIC-REINFORCED SURGICAL GOWNS

- SMS Fabric
- Extra Wide Fabric Reinforcements
- Raglan Sleeves
- Low Linting Fabric
- Colour Coded Neck Binding: Green
- Adjustable Neckline – Hook and Loop
- Sterile Gown packaged with 1 Towel
- Meet EN13795 Standard Performance

CODE STERILE	NON-STERILE	SIZE	LENGTH ⁺	STERILE	ITEMS/CASE NON-STERILE
	74911	Large	125 cm		32
95221	74921	X-Large	130 cm	28	28
-	74931	XX-Large	137 cm	-	28
95224	77329	X-Large, X-Long	150 cm	28	54

ULTRA SURGICAL GOWNS

- SMS Fabric
- Low Linting Fabric
- Raglan Sleeves
- Colour Coded Neck Binding: Yellow
- Adjustable Neckline - Hook and Loop
- Generous Cut
- Sterile Gown packaged with 2 Towels (Gown **95131** packaged with 1 Towel)
- Meet EN13795 Standard Performance

CODE STERILE	NON-STERILE	SIZE	LENGTH ⁺	STERILE	ITEMS/CASE NON-STERILE
-	47247	Small	117 cm	-	640
95113	-	Large	123 cm	32	-
95123	-	X-Large	128 cm	30	-
95131	71331	XX-Large	135 cm	28	36
95114	71626	Large, X-Long	138 cm	32	60
95124	71627	X-Large, X-Long	144 cm	30	48

ULTRA SURGICAL GOWN, WITH POLYPROPYLENE FABRIC

CODE	DESCRIPTION	SIZE	ITEMS/CASE
71321	ULTRA Surgical Gown with Polypropylene Fabric, Non-Sterile	X-Large	36

⁺ length back (neck to bottom)



Non-Reinforced Fabric



Fabric-Reinforced



High Performance Area

All codes in **bold** indicate a sterile product.

VALUESELECT*

VALUESELECT* SURGICAL GOWNS

STANDARD SURGICAL GOWNS

- SMS Fabric
- Book-fold
- Set-in-sleeves
- Low-lint fabric
- Sterile gown packaged with two towels
- Comfortable
- Lightweight
- Soft cuffs
- Long Velcro adjustable collar for maximum flexibility

CODE STERILE	NON-STERILE	SIZE	LENGTH*	ITEMS/CASE STERILE	ITEMS/CASE NON-STERILE
46665	46669	Medium	108 cm	26	513
46666	46670	Large	117 cm	24	486
46667	46671	X-Large	126 cm	22	486
46668	46672	XX-Large	135 cm	22	486
47285	47281	Large, X-Long	150 cm	20	480
47286	47282	X-Large, X-Long	150 cm	20	480
70433	70432	XX-Large, X-Long	155 cm	20	420

HIGH-PERFORMANCE PE-REINFORCED SURGICAL GOWNS


- SMS Fabric
- Book-fold
- Set-in-sleeves
- Low-lint fabric
- Sterile gown packaged with two towels
- Comfortable
- Lightweight
- Soft cuffs
- Long Velcro adjustable collar for maximum flexibility

CODE STERILE	NON-STERILE	SIZE	LENGTH*	ITEMS/CASE STERILE	ITEMS/CASE NON-STERILE
46673	46677	Medium	108 cm	22	459
46674	46678	Large	117 cm	22	432
46675	46679	X-Large	126 cm	22	432
46676	46680	XX-Large	135 cm	22	478
47287	47283	Large, X-Long	150 cm	20	420
47288	47284	X-Large, X-Long	150 cm	20	360
70335	70434	XX-Large, X-Long	155 cm	18	360

* length back (neck to bottom)

All codes in **bold** indicate a sterile product.

See Back Cover for Trademark Information



AERO CHROME*

AERO CHROME* SURGICAL GOWNS



AERO CHROME* BREATHABLE PERFORMANCE SURGICAL GOWNS

- Red collar to identify a high performance gown
- Breathable back panel maximises airflow
- Distinctive colour makes it easy to choose the AAMI level 4 protection
- Excellent ignition resistance, with the highest rating for protection against ignition with surgical lasers
- Soft, light and breathable fabric keeps staff protected yet cool and comfortable
- Low linting to help reduce the risk of wound contamination and infection
- Meet the High Performance per EN13795-1-2019
- Meets AAMI Level 4 standard for surgical gowns, providing the highest level of fluid and microbial protection in the critical zones
- Sterile gown packaged with 2 towels

CODE		SIZE	LENGTH*	ITEMS/CASE
STERILE	NON-STERILE			STERILE NON-STERILE
48395	48392	Large	125 cm	30 40
48397	48399	Large, X-Long	145 cm	30 36
48396	48393	X-Large	131 cm	28 32
48398	48400	X-Large, X-Long	151 cm	30 32
48402	48401	XX-Large, X-Long	137 cm	28 30
48403	48394	XX-Large	137 cm	28 32

* length back (neck to bottom)



Non-Reinforced Fabric



Fabric-Reinforced



High Performance Area

All codes in **bold** indicate a sterile product.

OPERATING ROOM SAFETY PRODUCTS

SAFETY PRODUCTS

LIGHT HANDLE COVERS

STANDARD LIGHT HANDLE COVER

- Ensures peace of mind for clinicians that the sterile field is protected
- Needs separate adaptor that fits specific operating room lights
- Phthalates-free
- Latex-Free
- PVC-free

CODE	DESCRIPTION	QUANTITY/CASE
AR-LHCB01	1 PACK, Sterile	100
AR-LHCB02	2 PACK, Sterile	120
AR-LHCBNS	Non-Sterile	300

A large range of Light Handle Adaptors for disposable light handle covers is available. Please contact your account manager.

SKIN MARKER

- With flexible 15cm ruler
- Convenient packaging
- Non-toxic
- Non-tattooing
- Latex-free
- Phthalates-free
- PVC-free

CODE	DESCRIPTION	BOX QUANTITY	CASE QUANTITY
AR-26001	Standard Tip, blue barrel, Sterile	25	150
AR-26001NS	Standard Tip w/Ruler, Non-Sterile	-	100
AR-26002	Fine Tip, white barrel, Sterile	25	150
AR-26002NS	Fine Tip w/Ruler, white barrel, Non-Sterile	-	150

ESMARK BANDAGE

- Advanced washing technology results in minimum powder
- Latex-free
- Phthalates-free
- PVC-free

CODE	SIZE	DESCRIPTION	QUANTITY/CASE
AR-35309	7.5 x 275 cm	Sterile	20
AR-35409	10 x 275 cm	Sterile	20
AR-35609	15 x 275 cm	Sterile	20
AR-35612	15 x 365 cm	Sterile	20
AR-35609NS	15cm x 275cm	Non-Sterile	20

All codes in **bold** indicate a sterile product.

SAFETY PRODUCTS

TUBE POSITIONER

- Adhesive Backing
- Easily detachable by gloved personnel
- With integrated radiopaque to ensure x-ray detectability
- Latex-free
- Phthalates-free
- PVC-free
- Sterile

CODE	DESCRIPTION	QUANTITY/CASE
AR-30802	Tube Positioner	50

NEEDLE COUNTERS

STANDARD NEEDLE COUNTER

- Surface magnet for safe storage
- Foam for secure needle parking
- Clearly markings on magnets and foam
- Simple latch opening mechanism
- High contrast numbers
- Latex-free
- Phthalates-free
- PVC-free
- Sterile

CODE	DESCRIPTION	COUNTS	CAPACITY	Q/CASE
AR-250101	Double Magnet	20	40	100
AR-250103	Single Magnet, Clear top	10	20	100
AR-250104	Double Magnet, Horizontal	30	60	100
AR-250401	Foam Block/Magnet	20	40	100
AR-250802	Foam Strip	15	15	100
AR-250504	Foam Adhesive	40	80	50
AR-250103Y	Yellow, Single Magnet, Clear Top, Sterile	10	20	100
AR-250103NS	Red, Single Magnet, Clear Top, Non-Sterile	10	20	100
AR-250104NS	Red, Double Magnet, Horizontal, Non-Sterile	30	60	100
AR-250401NS	Red, Foam Block/Magnet, Non-Sterile	20	40	100
AR-250502NS	Red, Foam Block/Adhesive, Non-Sterile	20	40	100

All codes in **bold** indicate a sterile product.

PREMIUM NEEDLE COUNTER

- Including blade removal feature
- Overlapping box edges
- Automatic passive locking latch
- Three interlocking hinges
- Latex-free
- Phthalates-free
- PVC-free

CODE	DESCRIPTION	COUNTS	CAPACITY	Q/CASE
AR-251104	Double Magnet, Sterile	30	60	100
AR-251108	Double Magnet, Sterile	60	120	108
AR-251401	Foam Block/Magnet, Sterile	20	40	80
AR-251502	Foam Block/Adhesive, Sterile	20	40	80
AR-251701	Magnet Foam/strips, Sterile	15	30	100
AR-251104Y	Yellow, Double Magnet, Blade Remover, Sterile	30	60	100
AR-251104NS	Red, Double Magnet, Blade Remover, Non-Sterile	30	60	100
AR-251104YNS	Yellow, Double Magnet, Blade Remover, Non-Sterile	30	60	100
AR-251108NS	Red, Double Magnet, Blade Remover, Non-Sterile	60	120	100
AR-251401NS	Red, Foam Block/Magnet, Blade Remover, Non-Sterile	20	40	64
AR-251502NS	Red, Foam Block/Adhesive, Blade Remover, Non-Sterile	20	40	100
AR-251504NS	Red, Foam Block/Adhesive, Blade Remover, Non-Sterile	40	80	100

INSTRUMENT PADS

MAGNETIC INSTRUMENT PAD

- Non-slip foam backing
- Special magnetic surface
- Latex-free
- Phthalates-free
- PVC-free

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
AR-25001	Sterile	25 x 40cm	30
AR-25011	Non-slip, Sterile	50cm x 40 cm	30
AR-25001NS	Non-Sterile	25 cm x 40 cm	100

INSTRUMENT PAD, NON-MAGNETIC

CODE	DESCRIPTION	QUANTITY/CASE
AR-25011NS	Non-magnetic Instrument pad, Non-Sterile	200

All codes in **bold** indicate a sterile product.

EQUIPMENT COVERS

C-ARM DRAPE WITH STRAPS

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
51717	C-Arm Drape with Straps, Sterile	Large, 104 cm x 187 cm	40
51718	C-Arm Drape with Straps, Sterile	Mini, 137 cm x 213 cm	25
51817	C-Arm Drape with Straps, Non-Sterile	Large, NS, 104 cm X 187 cm	56
51818	C-Arm Drape with Straps, Non-Sterile	Mini, NS, 137 cm X 213 cm	50

PE SHEET GENERAL SELF-ADHESIVE TAPE ON 4 SIDES

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
51719	Sterile	104 cm x 162 cm	70
51819	Non-Sterile	NS, 104 cm X 162 cm	100

DOME BANDED BAG

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
51720	Dome Banded Bag, Sterile	45 cm Depth	80
51721	Dome Banded Bag, Sterile	55 cm Depth	80
51722	Dome Banded Bag, Sterile	66 cm Depth	80
51723	Dome Banded Bag, Sterile	76 cm Depth	80
51820	Dome Banded Bag, Non-Sterile	45 cm Depth	120
51821	Dome Banded Bag, Non-Sterile	55 cm Depth	120
51822	Dome Banded Bag, Non-Sterile	66 cm Depth	120
51823	Dome Banded Bag, Non-Sterile	76 cm Depth	114

All codes in **bold** indicate a sterile product.

SQUARE BANDED BAG WITH TAPE

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
51724	Square Banded Bag with Tape, Sterile	76 cm X 76 cm	90
51727	Square Banded Bag with Tape, Sterile	76 cm X 91 cm	90
51728	Square Banded Bag with Tape, Sterile	91 cm X 76 cm	86
51725	Square Banded Bag with Tape, Sterile	91 cm X 91 cm	100
51729	Square Banded Bag with Tape, Sterile	91 cm X 104 cm	80
51726	Square Banded Bag with Tape, Sterile	101 cm X 101 cm	80
51824	Square Banded Bag with Tape, Non-Sterile	NS, 76 cm X 76 cm	135
51827	Square Banded Bag with Tape, Non-Sterile	NS, 76 cm X 91 cm	135
51828	Square Banded Bag with Tape, Non-Sterile	NS, 91 cm X 76 cm	100
51825	Square Banded Bag with Tape, Non-Sterile	NS, 91 cm X 91 cm	150
51829	Square Banded Bag with Tape, Non-Sterile	NS, 91 cm X 104 cm	100
51826	Square Banded Bag with Tape, Non-Sterile	NS, 101 cm X 101 cm	120

MICROSCOPE DRAPE

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
51730	Zeiss MD, 65mm, Sterile	116 cm X 304 cm	10
51734	Zeiss MD, 65mm, Sterile	137 cm X 381 cm	16
51731	Zeiss OPML, 48 mm, Sterile	137 cm X 381 cm	10
51732	Leica, 68mm, Sterile	116 cm X 304 cm	10
51733	Leica, 68mm, Sterile	137 cm X 381 cm	16
51830	Zeiss 65mm, Non-Sterile	NS, 116 cm X 304 cm	50
51831	Zeiss 48mm, Non-Sterile	NS, 137 cm X 381 cm	50
51832	Leica 68mm, Non-Sterile	NS, 116 cm X 304 cm	50
51833	Leica, 68mm, Non-Sterile	NS, 132 cm X 381 cm	50
51834	65mm	NS, 137 cm X 381 cm	50

PROBE COVER PE, WITH GEL, ELASTIC AND STRIPS

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
51735	Probe Cover PE, with Gel, Elastic and Strips, Sterile	15 cm X 122 cm	60
51736	Probe Cover PE, with Gel, Elastic and Strips, Sterile	15 cm X 243 cm	60
51835	Probe Cover with Gel, Elastic and Strips, Non-Sterile	NS, 15 cm X 122 cm	100
51836	Probe Cover with Gel, Elastic and Strips, Non-Sterile	NS, 15 cm X 243 cm	100

All codes in **bold** indicate a sterile product.

SAFETY PRODUCTS

ADHESIVE PROBE COVER WITH GEL

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
51737	Adhesive Probe Cover With Gel, Sterile	15 cm X 238 cm	100
51837	Adhesive Probe Cover With Gel, Non-Sterile	NS, 15 cm X 238 cm	100

PULL-UP PROBE COVER WITH GEL AND TAPE

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
51738	Pull-Up Probe Cover With Gel and Tape, Sterile	15 cm X 122 cm	100
51838	Pull-Up Probe Cover With Gel and Tape, Non-Sterile	NS, 15 cm X 122 cm	200

VIDEO DRAPE

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
51739	Video Drape, Sterile	12 cm X 243 cm	60
51839	Video Drape, Non-Sterile	NS, 12 cm X 243cm	100

O.R. TOURNIQUET COVER

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
53600	O.R. Tourniquet Cover, Non-Sterile	Small	30
53605	O.R. Tourniquet Cover, Non-Sterile	Large	30

All codes in **bold** indicate a sterile product.

SAFETY ACCESSORIES

PERSONAL STIRRUP STRAP WITH SLIP RING

CODE	DESCRIPTION	QUANTITY/CASE
52712	Personal Stirrup Strap with Slip Ring, Non-Sterile	40

O.R. TABLE/STRETCHER STRAP

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
54800	O.R. Table/Stretcher Strap, Non-Sterile	8 cm X 76 cm	48

SPONGE STICK

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
AR-000182NS	Orange, Non-Sterile	20 cm	1000
AR-003415NS	Blue, Non-Sterile	15 cm	1000

SPONGE FORCEPS NYLON

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
AR-260175NS	Green, Non-Sterile	19 cm	500

TOWEL CLAMP

CODE	DESCRIPTION	QUANTITY/CASE
AR-30900NS	Towel Clamp, Non-Sterile	500

TEMPORARY SHARPS HOLDER

CODE	DESCRIPTION	QUANTITY/CASE
AR-400456	Sterile	50
AR-400456NS	Non-Sterile	100

All codes in **bold** indicate a sterile product.

SAFETY PRODUCTS

SPONGE COUNTER BAG

CODE	DESCRIPTION	QUANTITY/CASE
AR-25500NS	Sponge Counter Bag, Non-Sterile	250

SCALPEL HOLDER

CODE	DESCRIPTION	QUANTITY/CASE
AR-30920	Sterile	120
AR-30920NS	Non-Sterile	200

CAUTERY TIP CLEANER

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
AR-32200	Cautery Tip Cleaner, Sterile	5 cm x 5 cm	300
AR-32200NS	Cautery Tip Cleaner, Non-Sterile		1000

WASTE COLLECTION DEVICE

CODE	DESCRIPTION	QUANTITY/CASE
AR-77301371	Waste Disposal Device, Sterile	50
AR-77301371NS	Waste Disposal Device, Non-Sterile	100

All codes in **bold** indicate a sterile product.

ENDOSCOPY, LAPAROSCOPY & ELECTROSURGICAL PRODUCTS

ELECTROSURGICAL TIP CLEANER

- Coated Scratch Pad
- Extended adhesive peel-off backing
- Easy removal by gloved staff
- Facilitates quick and easy instrument cleaning
- Latex-free
- Phthalates-free
- PVC-free
- Sterile

CODE	DESCRIPTION	QUANTITY/CASE
AR-32200	Cautery Tip Cleaner, 5 x 5cm	1000

ANTI-FOG SOLUTION KIT

- Bottle capacity: 6ml
- Adhesive foam application pad
- X-ray detectable
- Latex-free
- Phthalates-free
- PVC-free
- Sterile

CODE	DESCRIPTION	QUANTITY/CASE
AR-7304412	Anti-fog Solution Kit, Sterile	48

All codes in **bold** indicate a sterile product.

TEMPERATURE MANAGEMENT

PATIENT WARMING BLANKETS

- Three Layer Warmth
- Shrink-wrapped
- Cost-Effective, electricity free
- For reflective blanket - green side to be placed towards the body
- Latex-free
- Phthalates-free
- PVC-free
- Non-Sterile

CODE	DESCRIPTION	QUANTITY/CASE
AR-4115175	Standard Patient Warming Blanket/Large 130x 190 cm	50

BABY BLANKETS

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
AR-003177	White with Pink/Blue Stripe, Sterile	100cm X 85cm	60
AR-003177NS	Baby Blanket, Pink/Blue, Non-Sterile	-	100

BABY CAPS

CODE	DESCRIPTION	QUANTITY/CASE
AR-1336	Pink/Blue, Stripe, Sterile	50
AR-1336NS	Pink/Blue, Stripe, Non-Sterile	100

All codes in **bold** indicate a sterile product.

HOLLOWARE MULTIPURPOSE BOWLS, BASINS AND TRAYS

PLACENTA BASINS

- Graduated in ml. cc and oz
- Optional lid available
- ETO Sterilisable
- No Latex
- No Phthalates
- Lot Traceable

CODE	DESCRIPTION	QUANTITY/CASE
AR-3000	Basin Placenta, 3000 ml	100
AR-9300	Lid Basin Placenta	100

ESMESIS BASIN

- Flexible Design
- Burr-free Edges
- Smooth round Edges
- Graduated in ml. cc and oz
- ETO Sterilisable
- No Latex
- No Phthalates
- Lot Traceable

CODE	DESCRIPTION	QUANTITY/CASE
AR-0700	Basis Emesis, 700 ml	320

RING BASINS

- Fits all standard ring stands
- AR-6000 graduated in ml. cc and oz
- ETO sterilisable
- No Latex
- No Phthalates
- Lot Traceable

CODE	DESCRIPTION	QUANTITY/CASE
AR-6000	Basin Ring, 6000 ml	60

All codes in **bold** indicate a sterile product.

HOLLOWARE MULTIPURPOSE BOWLS, BASINS AND TRAYS

MULTIPURPOSE TRAYS

- Flexible design
- Burr-free edges preventing particle shedding
- Thick corners preventing cracking or breaking
- 'Clean' Product, no wiping necessary
- ETO sterilisable
- No Latex
- No Phthalates
- Lot Traceable

CODE	DESCRIPTION	SIZE	QUANTITY/CASE
AR-0541	Tray Shallow Small	11 x 13 x 2.5 cm	600
AR-0542	Tray Deep Small	11 x 13 x 5 cm	600
AR-0904	Tray 2 Compartment	24 x 11 x 4.5 cm	300
AR-0910	Tray 3 Compartment	25 x 24 x 4.9 cm	150
AR-0911	Tray J-Wire	27 x 22 x 5 cm	150
AR-0952	Tray Deep Large	24.5 x 14 x 5 cm	300
AR-1115	Tray Platform Small	28 x 38 x 2.5 cm	80
AR-1219	Tray Mayo	43 x 28 x 1.5 cm	40
AR-1416	Tray Platform Flat	38 x 24 x 2.5 cm	50
AR-1622	Tray Mayo Large	52 x 38 x 2 cm	30
AR-0951	Tray Shallow Large	14 x 24 x 2.5 cm	300

BOWLS

- Graduated in ml, cc and oz
- Optional lids available
- ETO sterilisable
- No Latex
- No Phthalates
- Lot Traceable

CODE	DESCRIPTION	QUANTITY/CASE
AR-0250	Sponge Bowl, 250 ml	600
AR-9025	Lid Bowl, 250 ml	300
AR-0500	Sponge Bowl, 500 ml	390
AR-1000	Sponge Bowl, 1000 ml	340
AR-9100	Lid Bowl, 1000 & 1500 ml	160
AR-1500	Sponge Bowl, 1500 ml	320
AR-2500	Sponge Bowl, 2500 ml	160
AR-9250	Lid Bowl, 2500ml & Bowl G/W	120
AR-GW2000CT	Guidewire Bowl, 2500 ml	160

All codes in **bold** indicate a sterile product.

HOLLOWARE MULTIPURPOSE BOWLS, BASINS AND TRAYS

RECTANGULAR BASIN

- ETO Sterilisable
- No Latex
- No Phthalates
- Lot Traceable

CODE	DESCRIPTION	QUANTITY/CASE
AR-8000	Rectangular Basin 8000ml	50

MEDICINE CUPS

- Various colours available
- Tapered design for stacking
- Graduated in cc and oz
- ETO Sterilisable
- No Latex
- No Phthalates
- Lot Traceable

CODE	DESCRIPTION	COLOUR	QUANTITY/CASE
AR-0060	Cup medicine, 60 ml	Clear	1150
AR-0060-BL	Cup medicine, 60 ml	Blue	1150
AR-0060-RD	Cup medicine, 60 ml	Red	1150
AR-0060-GN	Cup medicine, 60 ml	Green	1150

All codes in **bold** indicate a sterile product.

INDEX BY CODE

CODE	PAGE
11783	118
11784	118
13227	51, 66
13228	64
13231	32
13232	54, 56, 64
13233	55, 64
13294	61
13296	97
13337	118
13815	118
14007	110
14364	29
14365	29
29148	103
29149	103
29219	64
29220	64
29221	69
29222	64
29223	65
29224	69
29225	100
29226	100
29229	99
29231	99
32109	118
32113	118
33412	68
33415	68
33916	118
33917	118
42217	74, 111
42222	15, 16, 17, 25, 26, 27, 50, 51, 53, 54, 56, 60, 70, 71, 72, 73, 74, 79, 80, 86, 104, 107, 111
42224	18, 22, 26, 55, 57, 59, 60, 61, 87, 90, 91, 105, 111
42225	111
42228	111
42302	75, 111
42309	111
42312	55, 58, 82, 90, 111
42313	55, 58, 82, 90, 111
42345	51, 58, 72, 87, 90, 91, 111
42365	72, 87, 111
42366	111
42367	111
42409	111

CODE	PAGE
42428	15, 16, 17, 25, 26, 27, 50, 51, 53, 54, 56, 60, 70, 71, 72, 73, 74, 79, 80, 86, 104, 107, 111
42446	28, 111
42447	79, 111
42449	59, 111
42451	18, 19, 22, 23, 26, 54, 55, 57, 59, 60, 61, 87, 90, 91, 105, 106, 111
42480	111
42481	111
44591	114
46665	123
46666	123
46667	123
46668	123
46669	123
46670	123
46671	123
46672	123
46673	123
46674	123
46675	123
46676	123
46677	123
46678	123
46679	123
46680	123
47247	122
47281	123
47282	123
47283	123
47284	123
47285	123
47286	123
47287	123
47288	123
47348	20, 24, 51, 52, 53, 110
47614	41, 53, 56, 60, 74, 82, 91
47616	41, 53, 56, 60, 74, 82, 91
47618	41
47622	25, 41, 51, 55
47623	25, 41, 51, 55
47631	41
48392	124
48393	124
48394	124
48395	124
48396	124

CODE	PAGE
48397	124
48398	124
48399	124
48400	124
48401	124
48402	124
48403	124
48755	58
48756	56
48757	56
48758	57, 117
48759	113
48760	114
48957	95
48958	95
48959	95
48986	95
48987	95
48988	95
48989	96
48990	96
49000	96
49016	96
49017	96
49020	96
49021	96
49035	96
49036	95
49037	95
49038	95
49039	
49041	115
49042	115
49043	116
49044	116
49045	116
49046	112
49047	113
49048	113
49049	113
49050	113
49051	115
49052	115
49053	116
49054	116
49055	116
49056	112
49057	113
49058	113

CODE	PAGE
49059	113
49060	113
49069	96
49070	96
49071	95
49072	95
49073	96
49074	96
49075	96
49076	96
49077	96
49078	96
49079	96
49080	96
51717	128
51718	128
51719	128
51720	128
51721	128
51722	128
51723	128
51724	129
51725	129
51726	129
51727	129
51728	129
51729	129
51730	129
51731	129
51732	129
51733	129
51734	129
51735	129
51736	129
51737	130
51738	130
51739	130
51817	128
51818	128
51819	128
51820	128
51821	128
51822	128
51823	128
51824	129
51825	129
51826	129
51827	129
51828	129

INDEX BY CODE

CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE
51829	129	76291	25, 43, 45	77723	68	79207	28, 108
51830	129	76293	44	77744	82	79208	28, 108
51831	129	76300	82, 114	77748	90	79209	29
51832	129	76321	84	77749	90	79214	15, 36, 104
51833	129	76323	40, 87	77760	82	79215	15, 39, 104
51834	129	76334	116	77769	23	79216	16, 38
51835	129	76340	90, 116	77772	19	79217	16, 36
51836	129	76346	85	77776	100	79218	16, 39
51837	130	76353	71, 77	77779	74	79222	49
51838	130	76354	95	77781	50	79223	49
51839	130	76356	25, 90, 94	77783	51	79236	48, 109
52712	131	76357	90, 94	77790	54	79265	68
53600	130	76358	110	77792	27	79267	57, 69
53605	130	76361	22, 23, 32, 106	77794	22, 106	79273	50, 65
54800	131	76363	17, 22, 23, 35, 57, 106	77798	54	79275	60, 62
59705	40	76364	23, 48	77823	86	79279	53, 54, 63
59706	54	76365	17, 19, 38	77853	23, 106	79282	46
59709	25	76366	19, 36	77865	72	79284	118
59710	85	76367	19, 39	77870	92	79285	118
68474	117	76368	33, 57, 60	77873	84	79293	28, 42, 43
68500	117	76372	111	77881	53	79294	118
68501	117	76374	26	77883	110	79295	118
68502	117	76375	68	77885	80	79304	82, 83
70335	123	76377	27, 49	77888	73	79305	83
70391	120	76378	90, 94	77890	73	79313	42, 90
70392	120	76380	27, 46	77892	91	79330	55, 57, 61, 62
70432	123	76382	101, 109	78657	110	79333	59, 60, 63
70433	123	76384	74, 78	78660	62	79343	43
70434	123	76386	100	78669	20,	79345	59, 66
71321	122	76393	85	24, 26, 50, 51, 52, 53, 80, 111		79346	59, 66
71331	122	76394	70, 71, 72, 109	78911	119	79348	58, 61, 62
71626	122	76396	89	78913	119	79351	86, 87, 88
71627	122	76492	25, 43	78921	119	79356	84
72978	120	76681	100	78923	119	79360	35, 58, 91
72988	120	76885	80, 81	78961	17, 60, 90	79363	33, 58
74911	122	76888	73, 76	78965	44	79368	100
74921	122	77329	122	79010	18, 39, 105	79375	35
74931	122	77336	121	79011	18, 37, 105	79408	72, 74, 75, 79, 114
76209	87	77337	121	79012	18, 38, 105	79411	72, 112
76214	82, 114	77382	101, 109	79053	15, 16, 17, 18, 23,	79415	112
76216	26	77450	51, 66	26, 50, 51, 54, 55, 56, 57, 58, 71,		79419	120
76226	51	77452	110	90, 91, 104, 105, 106, 107, 115		79430	86, 87, 88
76234	87	77454	70	79063	99	79434	120
76243	28	77455	70	79070	65	79453	93
76246	55	77466	89	79075	89	79492	65
76248	91	77491	71	79102	56, 74, 91	79528	70, 72, 76
76249	58	77492	72	79109	91, 99	79529	70, 71, 76
76263	27	77707	85	79132	42	79532	24, 33, 34
76270	115	77717	26, 107	79206	28, 108	79540	79

INDEX BY CODE

CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE
79543	88	88233	61	89291	63	89685	52
79545	79, 81	88243	60	89293	42	89686	26
79551	75	88253	59	89301	32	89687	52
79582	97	88354	71	89304	82, 83	89701	15, 16, 17, 18, 19, 22, 23, 25, 26, 54, 58, 60, 61, 70, 71, 72, 73, 74, 75, 79, 80, 82, 90, 91, 104, 105, 106, 107, 114
79583	101	88373	60	89305	83	89714	97
79601	51, 110	88511	74	89311	33, 58	89715	98
79602	110	88521	74	89313	42, 90	89716	98
79603	15, 16, 17, 18, 19, 22, 23, 25, 26, 27, 54, 55, 56, 57, 58, 59, 60, 61, 71, 72, 73, 86, 90, 91, 104, 105, 106, 110	88551	75	89331	55, 57, 61, 62	89791	115
79692	29	88611	79	89343	43	89811	55, 113
79693	30	88621	79	89345	59, 66	89821	57, 59, 114
79694	41, 50, 51, 52	88673	20, 24, 26, 50, 51, 52, 53, 80, 111	89348	58, 61, 62	89945	82, 91, 115
79695	66	88760	17	89352	86, 87, 88	90002	120
79696	41	88762	15, 104	89353	40, 87	90012	51, 120
79697	113	88763	15, 104	89363	28, 43	90013	26, 82, 120
79699	42, 75	88764	16	89368	100	90016	120
79714	97	88765	16	89381	35, 58, 91	90019	120
79715	98	88767	17	89408	72, 74, 75, 79, 114	90042	120
79716	98	88900	20	89411	74, 112	90043	51, 120
79720	15, 16, 17, 18, 19, 22, 23, 25, 26, 54, 58, 60, 61, 70, 71, 72, 73, 74, 75, 79, 80, 82, 90, 91, 104, 105, 106, 107, 114	88901	20	89413	112	90046	120
79731	44, 74, 75	88902	24	89415	112	90049	120
79741	75, 107, 111	88904	24	89452	93	90112	119
79742	80, 81	88907	53	89471	98	90113	119
79744	85	88908	53	89492	65	90142	119
79745	26, 30	88909	52	89528	70, 72, 76	90143	119
79749	99	88910	51	89529	70, 71, 76	90363	121
79782	40	88911	75	89543	88	90373	121
79783	43, 80, 107	89010	88	89583	101	95113	122
79784	30	89048	91, 99	89584	75, 112	95114	122
79787	20, 36	89053	15, 16, 17, 18, 23, 26, 50, 51, 54, 55, 56, 57, 58, 71, 90, 91, 104, 105, 106, 107, 115	89601	51, 110	95123	122
79790	20, 24, 39	89062	89	89603	15, 16, 17, 18, 19, 22, 23, 25, 26, 27, 54, 55, 56, 57, 58, 59, 60, 61, 71, 72, 73, 86, 90, 91, 104, 105, 106, 110	95124	122
79945	82, 91, 115	89070	65	89632	75, 107, 111	95131	122
79970	51, 120	89101	53, 56, 60, 74, 82	89638	29	95221	122
79973	26, 82, 120	89204	81	89670	29	95224	122
79980	120	89206	28, 108	89671	30	99294	118
79983	51, 120	89207	28, 108	89672	41, 50, 51, 52	99295	118
86357	90, 94	89208	28, 108	89673	52, 66	AR-0060	137
88100	18, 105	89209	29	89674	41	AR-0060-BL	137
88103	18, 105	89236	48, 109	89675	113	AR-0060-GN	137
88113	22	89265	68	89676	43	AR-0060-RD	137
88183	59	89267	57, 69	89677	42, 75	AR-000182NS	131
88203	55	89271	50, 65	89678	40	AR-0250	136
88213	57	89275	60, 62	89679	43, 80, 107	AR-0500	136
88223	86	89276	64	89680	30	AR-0541	136
		89278	55, 64	89681	80	AR-0542	136
		89279	53, 54, 63	89682	107	AR-0700	135
		89282	46	89684	50	AR-0904	136

INDEX BY CODE

CODE	PAGE	CODE	PAGE	CODE	PAGE	CODE	PAGE
AR-0910	136	AR-250103	126	K792032	32	K797316	42
AR-0911	136	AR-250103NS	126	K792033	31	K797320	42
AR-0951	136	AR-250103Y	126	K792034	31	K797321	42
AR-0952	136	AR-250104	126	K792035	31	K797322	42
AR-1000	136	AR-250104NS	126	K792036	31	K887603	20
AR-1115	136	AR-250401	126	K792130	38	K887621	21
AR-1200	116	AR-250401NS	126	K792245	47	K887623	21
AR-1219	136	AR-250502NS	126	K792246	47	K887625	21
AR-1336	134	AR-250504	126	K792610	29	K887626	22
AR-1336NS	134	AR-250802	126	K793043	83	K891022	40
AR-1416	136	AR-251104	127	K793050	84	K891023	40
AR-1500	136	AR-251104NS	127	K793053	84	K891024	40
AR-1622	136	AR-251104Y	127	K793210	63	K891025	40
AR-2500	136	AR-251104YNS	127	K793214	33	K891027	30
AR-3000	135	AR-251108	127	K793330	34	K892031	31
AR-003177	134	AR-251108NS	127	K793463	67	K892032	32
AR-003177NS	134	AR-251401	127	K793464	67	K892033	31
AR-003415NS	131	AR-251401NS	127	K793467	67	K892034	31
AR-6000	135	AR-251502	127	K793468	68	K892036	31
AR-8000	137	AR-251502NS	127	K793630	33	K892245	47
AR-9025	136	AR-251504NS	127	K793632	69	K892246	47
AR-9100	136	AR-251701	127	K793703	35	K892791	63
AR-9250	136	AR-260175NS	131	K794530	93	K893330	34
AR-9300	135	AR-400456	131	K795610	37	K893463	67
AR-25001	127	AR-400456NS	131	K795613	37	K893464	67
AR-25001NS	127	AR-4115175	134	K795614	37	K893467	67
AR-25011	127	AR-7304412	133	K795615	37	K893468	68
AR-25011NS	127	AR-77301371	132	K795616	37	K894530	93
AR-25500NS	132	AR-77301371NS	132	K795711	40	K897213	77
AR-26001	125	AR-GW2000CT	136	K795713	39	K897215	77
AR-26001NS	125	AR-LHCB01	125	K795714	39	K897312	44
AR-26002	125	AR-LHCB02	125	K795805	38, 107	K897320	42
AR-26002NS	125	AR-LHCBNS	125	K795809	38	K897321	42
AR-28502	114	K103502	110	K795830	102	K897322	42
AR-28503	114	K103600	110	K796011	110		
AR-28504	114	K763570	94	K796012	110		
AR-30802	126	K768550	80	K796019	110		
AR-30900NS	131	K790790	32	K796710	45		
AR-30920	132	K791021	44	K797064	102		
AR-30920NS	132	K791022	40	K797065	98		
AR-32200	132, 133	K791023	40	K797067	101		
AR-32200NS	132	K791024	40	K797160	101		
AR-35309	125	K791025	40	K797162	102		
AR-35409	125	K791026	40	K797165	102		
AR-35609	125	K791027	30	K797166	102		
AR-35609NS	125	K791412	35	K797167	103		
AR-35612	125	K791413	44	K797169	102		
AR-63001NS	114	K792015	31, 94	K797215	77		
AR-250101	126	K792031	31	K797310	44		

Halyard, now part of the Owens & Minor family, is focused on advancing health and healthcare by delivering clinically-superior infection prevention and surgical solutions. Halyard sells its recognised sterilisation wrap, facial protection, gloves, protective apparel, surgical drapes and gowns in more than 100 countries, and holds leading market positions across the portfolio. For more information, visit www.halyardhealth.eu.

**Our
Comittement
to you**

**Better
infection
prevention**

For more information,
please contact your O&M Halyard Customer Care Team:

UK & IE (toll-free)

Tel: 0800 085 82 04 - Fax: 0800 066 91 06
customerservice.uk.ie@owens-minor.com

Other countries

Tel: +32 2 585 51 14 - Fax: +32 2 897 00 95
customerservice.export@owens-minor.com



www.halyardhealth.eu

*Registered Trademark or Trademark of O&M Halyard or its affiliates.
©2018. All rights reserved. 008734 HC109-08 - EN GB