

PURPLE NITRILE* MAX

Powder-Free Exam Gloves



THE NEW PURPLE STANDARD FOR PROTECTION

Once again, Halyard Health proves "Purple Means Protection" with the introduction of PURPLE NITRILE* MAX, the newest member of our single-use exam glove family.

Designed specifically to meet the needs of regular glove users who require the extra protection of a glove that reaches to the mid-forearm.



Durable, easy-to-don PURPLE NITRILE* MAX offers these critical safety features:

Longer Length – PURPLE NITRILE* MAX cuff extends a full 400 mm, providing the right level of protection from dangerous microorganisms, all the way to the mid-forearm.

Heavier Weight – At 20 grams, the PURPLE NITRILE* MAX glove is more than twice the weight of the PURPLE NITRILE-XTRA* glove, for added confidence when handling sharps and chemicals.

Thicker Fingertips – Its increased fingertip thickness provides almost double the protection of any other exam glove in the Halyard portfolio, yet still offers excellent tactile sensitivity.

Fully Textured – PURPLE NITRILE* MAX features fully textured fingertips and palm, designed for excellent wet and dry grip.

NEW!

PURPLE NITRILE* MAX POWDER-FREE EXAM GLOVES



With PURPLE NITRILE* MAX, glove users can be sure they're getting the protection they require, eliminating the need for double or triple gloving.

Combined with the savings realised from a disposable exam glove vs. the costs associated with reusable industrial gloves, PURPLE NITRILE* MAX is an excellent value.



PURPLE NITRILE* MAX Powder-Free Exam Gloves

Item	Description	Size	Gloves/Box	Boxes/Case	Gloves/Case
47195	PURPLE NITRILE* MAX Powder-Free Exam Glove	Small	50	8	400
47196	PURPLE NITRILE* MAX Powder-Free Exam Glove	Medium	50	8	400
47197	PURPLE NITRILE* MAX Powder-Free Exam Glove	Large	50	8	400
47198	PURPLE NITRILE* MAX Powder-Free Exam Glove	X-Large	50	8	400

Product Specifications

Gauge Thickness Measurements	MM
Middle Finger:	.28
Palm:	.24
Cuff:	.15
Average Length	415 mm


Physical Properties

Before Aging	
Force at break:	18 N
Ultimate Elongation:	550%

After Aging	
Force at break:	24 N
Ultimate Elongation:	450%

AQL - EN455-1

Acceptance Quality Limit	1.0 AQL
--------------------------	---------



KNOWLEDGE NETWORK* Clinical Education
 Knowledgeable Customer Support
 Expert Sales Force
 Tools & Best Practices
 Clinical Research
 Commitment to Excellence

For more information, please send an email to customerservice.uk.ie@hyh.com or visit www.halyardhealth.co.uk