



HALYARD*

Sterilization Wrap



PROTECT EVERY
INSTRUMENT
WITH CONFIDENCE



HALYARD*

Sterilization Wrap

The world's leading sterilization wrap, now available in sheet sizes and compact packages, designed for individual medical and dental practices

Invented in 1978 as an affordable, high-performance alternative to traditional wrap materials, our sterilization wrap remains the industry standard for sterile processing.

Now, medical and dental practices can take advantage of the same performance demanded by major medical and surgery centers. Independent laboratory testing assures that HALYARD* Wrap delivers superior bacterial filtration efficiency (BFE) from airborne pathogens, maintaining sterility of your instruments, at a competitive price.

EACH SHEET OF HALYARD* WRAP DELIVERS:

- Strength and tear resistance
- Superior bacterial filtration efficiency (BFE) for airborne microbes
- Low linting
- Flame resistance
- Compatibility with steam and Ethylene Oxide (EO) sterilant; not appropriate for dry heat sterilization



**WORLD-CLASS PERFORMANCE
AT A COMPETITIVE PRICE**

**THE RIGHT SIZES FOR
SMALLER PACKAGES**

**COMPACT, 100-COUNT BAGS
SAVE STORAGE SPACE**

HALYARD* H100 STERILIZATION WRAP

Compact Packaging

Stock No.	Size†	Sheets /Bag	Bags /Case	Sheets /Case
37159	9in. x 9in. / 22cm x 22cm	100	10	1000
37046	12in. x 12in. / 30cm x 30cm	100	10	1000
37047	15in. x 15in. / 38cm x 38cm	100	10	1000
37049	18in. x 18in. / 45cm x 45cm	100	10	1000
37050	20in. x 20in. / 50cm x 50cm	100	10	1000
37051	24in. x 24in. / 60cm x 60cm	100	5	500

† Larger sizes available.

Medical practices: www.halyardhealth.com

Dental practices: www.halyarddental.com



KNOWLEDGE NETWORK* Accredited Education

Ongoing Customer Support

Expert Sales Force

Tools & Best Practices

Clinical Research

Commitment to Excellence

www.halyardhealth.com



*Registered Trademark or Trademark of Halyard Health, Inc. or its affiliates. © 2015 HYH. All rights reserved.

SS63 C15476