



ArcRoyal | Surgical Solutions

Contents	SURGICAL PRODUCTS	
	Standard Light Handle Cover	5
	Camera Light Handle Cover	5
	Skin Marker	6
	Esmark Bandage	6
	Tube Positioner	7
	SAFETY PRODUCTS	
	Standard Needle Counter	8
	Premium Needle Counter with Blade Remover	9
	Magnetic Instrument Pad	10
	Sponge Counter Bag	11
	ENDOSCOPY, LAPAROSCOPY & ELECTROSURGICAL PRODUCTS	
	Electrosurgical Tip Cleaner	12
	Antifog Solution Kit	12
	TEMPERATURE MANAGEMENT	
	Patient Warming Blankets	13
	Ordering Information	14
	Mission Vision Values	15

Surgical Products

ArcRoyal provides an essential range of single use Surgical Solution Products used for a variety of clinical needs including patient safety, sharps and other surgical essentials in the Operating Room.



Key



Latex free



Phthalates free



PVC free



Non-sterile



Contains PVC

Light Handle Covers

Disposable light handle covers facilitate aseptic adjustment of surgical lamps and lower the procedure cost compared to reprocessing of reusable handles. They allow for secure attachment of adaptor to minimise risk of contamination and remain securely in place for the entire duration of procedure.



Standard Light Handle Cover

Product Specifications:

- Ensures peace of mind for clinicians that the sterile field is protected
- Needs separate adaptor that fits specific operating room lights



Camera Light Handle Cover

Product Specifications:

- Instantly usable in the Operating Room on OR Light Cameras
- Available in square large or round small



Ordering Information

Light Handle Covers

Item No.	Description	Case Quantity
AR-LHCB01	Standard Light Handle Cover, 1 pack - Blue, Sterile	100
AR-LHCB02	Standard Light Handle Cover, 2 pack - Blue, Sterile	120
AR-LHCAM01S	Camera Light Handle Cover, Square, Large, 1 pack, Sterile	100
AR-LHCAM02R	Camera Light Handle Cover, Round, Small, 1 pack, Sterile	100

Skin Marker

The prep-resistant gentian violet ink marker allows the surgeon to create smooth, clear markings on patient's skin

Product Specifications:

- Includes flexible 15cm ruler
- Convenient packaging minimises storage requirements
- Non-toxic
- Non-tattooing



Ordering Information

Skin Markers

Item No.	Description	Pack Quantity	Case Quantity
AR-26001	Skin Marker, Standard Tip w/Ruler (blue barrel), Sterile	25 per box	150
AR-26002	Skin Marker, Fine Tip w/Ruler (white barrel), Sterile	25 per box	150

Esmark Bandage

Designed with excellent compression and elastic properties to facilitate a bloodless field.

Product Specifications:

- Advanced washing technology results in minimum powder



Ordering Information

Esmark Bandages

Item No.	Description	Case Quantity
AR-35309	Esmark Bandage 7.5cm x 275cm, Sterile	20
AR-35409	Esmark Bandage 10cm x 275cm, Sterile	20
AR-35609	Esmark Bandage 15cm x 275cm, Sterile	20
AR-35612	Esmark Bandage 15cm x 365cm, Sterile	20

Tube Positioner

Adhesive backing allows for secure placement and facilitates convenient management of tubes, lines and cords in the sterile field



Product Specifications:

- Easily detachable by gloved personnel
- With integrated radiopaque element to ensure x-ray detectability




Ordering Information

Tube Positioner

Item No.	Description	Case Quantity
AR-30802	Tube Positioner, Sterile	50

Needle Counters

Designed to provide a safe and easy way to store and count various types of sharps used during a procedure



Standard Needle Counter

Product Specifications:

- Surface magnet allows for safe storage of metallic components
- Foam facilitates secure parking of needles, covering the sharp tip
- Clearly visible markings on magnets and foam allow easy sharps placement and positioning
- Numerous ergonomic considerations such as simple latch opening mechanism and high contrast numbers to facilitate correct count of blades and needles
- Available with a transparent lid allowing clear visualisation of contents even when closed (AR-250103)



Ordering Information

Standard Needle Counter

Item No.	Description	Specifications	Case Quantity
AR-250101	Needle Counter Red, Double Magnet, Sterile	20 count, 40 capacity	100
AR-250103	Needle Counter Red, Single Magnet, Clear Top, Sterile	10 count, 20 capacity	100
AR-250104	Needle Counter Red, Double Magnet, Horizontal, Sterile	30 count, 60 capacity	100
AR-250401	Needle Counter Red, Foam Block/Magnet, Sterile	20 count, 40 capacity	100
AR-250802	Needle Counter Red, Foam strip, Sterile	15 count, 15 capacity	100
AR-250504	Needle Counter Red, Foam Adh., Sterile	40 count, 80 capacity	50



Premium Needle Counter - with Blade Remover

Product Specifications:

- Incorporates a blade removal feature for quick and easy removal of contaminated blades
- Overlapping box edges prevent sharps from escaping when closed
- Passive locking latch engages automatically and securely
- Three interlocking hinges prevents inadvertent opening




Ordering Information

Premium Needle Counter - with Blade Remover

Item No.	Description	Specifications	Case Quantity
AR-251104	Needle Counter Red, Double Magnet, Blade Remover, Sterile	30 count, 60 capacity	100
AR-251108	Needle Counter Red, Double Magnet, Blade Remover, Sterile	60 count, 120 capacity	108
AR-251401	Needle Counter Red, Foam Block/Magnet, Blade Remover, Sterile	20 count, 40 capacity	80
AR-251502	Needle Counter Red, Foam Block/Adhesive, Blade Remover, Sterile	20 count, 40 capacity	80
AR-251701	Needle Counter Red, Magnet Foam/Strips, Sterile	15 count, 30 capacity	100

Instrument Pads

Designed to provide a safe surface to pass and store instruments and to reduce the risk of accidental sharps injuries. The products' non-slip foam backing adheres securely to patient drape and encourages safe sharps' passing practices



Magnetic Instrument Pad

Product Specifications:

- Instruments cling to special magnetic surface



Ordering Information

Instrument Pad

Item No.	Description	Case Quantity
AR-25001	Instrument Pad, Magnetic, 25 x 40cm, Sterile	30

Sponge Counter Bag

Designed to provide a secure containment for used sponges to reduce risk by minimising exposure to blood-borne pathogens. Facilitates easy and continuous sponge counting.

Product Specifications:

- Pocket divider seal easily separates to allow for larger sponges to blood-borne pathogens
- Clear polyethylene pockets allow clear content visualisation, promoting accurate counts and also allows sponges to be viewed for visual estimate of patient blood loss
- Additional bags can be placed over filled bags as necessary
- Pockets accommodate five laparotomy sponges or ten pieces of 10x10cm gauze
- Additional bags can be placed over filled bags as necessary



Ordering Information

Sponge Counter Bag

Item No.	Description	Pack Quantity	Case Quantity
AR-25500 NS*	Sponge Counter Bag	50 per box	250

* Non-Sterile

Electrosurgical Tip Cleaner

The coated scratch pad eliminates disruptive scraping of electrodes and serves as a convenient debris removal system for use during electrosurgical procedures

Product Specifications:

- Extended adhesive peel-off backing enables easy removal by gloved staff
- Adhesive backing provides for placement on mayo or patient drape
- Facilitates quick and easy instrument cleaning



Ordering Information

Electrosurgical Tip Cleaner

Item No.	Description	Pack Quantity	Case Quantity
AR-32200	Cautery Tip Cleaner, 5cm x 5cm, Sterile	50 per box	300

Anti-fog Solution Kit

The unique anti-fog solution is designed to prevent fogging of the endoscope lens during procedures

Product Specifications:

- Bottle capacity: 6ml
- Adhesive foam application pad
- X-ray detectable



Ordering Information

Anti-fog Solution Kit

Item No.	Description	Case Quantity
AR-7304412	Anti-fog Solution Kit, Sterile	48

Patient Warming Blankets

The Cozee™ range of single-use patient warming blankets is a cost-effective choice to prevent hypothermia during operation. Cozee™ blankets are available in two material options, and are designed to ensure maximum patient safety and comfort as well as maximum convenience for clinical staff.



Product Specifications:

- Three Layer Warmth 20% better insulated than cotton blankets
- Shrink-wrapped: Needs less space than traditional folded blankets
- Cost-Effective: No electricity used, lower cost compared to active warming systems



Note: For reflective blanket - green side to be placed towards the body

Ordering Information

Patient Warming Blankets


Item No.	Description	Case Quantity
AR-4115175	Standard Patient Warming Blanket/ Large 130 x 190 cm	50
AR-480161	Reflective Warming Blanket/ Large 122 x 213 cm	25

Non-Sterile

Ordering Information

To order, please contact us directly by email, telephone or fax, with the Item No (code) and the case quantity required.

Item No.	Description	Quantity per Case	Page Number
AR-LHCB01	Standard Light Handle Cover, 1 pack - Blue, Sterile	100	5
AR-LHCB02	Standard Light Handle Cover, 2 pack - Blue, Sterile	120	5
AR-LHCAM01S	Camera Light Handle Cover, Square Large 1 pack, Sterile	100	5
ARLHCAM-02R	Camera Light Handle Cover, Round Small 1 pack, Sterile	100	5
AR-26001	Skin Marker, Standard Tip w/ruler, Sterile	150	6
AR-26002	Skin Marker, Fine Tip w/ruler, Sterile	150	6
AR-35309	Esmark Bandage 7.5cm x 275cm, Sterile	20	6
AR-35409	Esmark Bandage 10cm x 275cm, Sterile	20	6
AR-35609	Esmark Bandage 15cm x 275cm, Sterile	20	6
AR-35612	Esmark Bandage 15cm x 365cm, Sterile	20	6
AR-30802	Tube Positioner, Sterile	50	7
AR-250101	Standard Needle Counter, Red, 20 count, 40 capacity, Sterile	100	8
AR-250103	Standard Needle Counter, Red, 10 count, 20 capacity, Sterile	100	8
AR-250104	Standard Needle Counter, Red, 30 count, 60 capacity, Sterile	100	8
AR-250401	Standard Needle Counter, Red, 20 count, 40 capacity, Sterile	100	8
AR-250802	Standard Needle Counter, Red, 15 count, 15 capacity, Sterile	100	8
AR-250504	Standard Needle Counter, Red, 40 count, 80 capacity, Sterile	50	8
AR-251104	Premium Needle Counter with Blade Remover, Red, 30 count, 60 capacity, Sterile	100	9
AR-251108	Premium Needle Counter with Blade Remover, Red, 60 count, 120 capacity, Sterile	108	9
AR-251401	Premium Needle Counter with Blade Remover, Red, 20 count, 40 capacity, Sterile	80	9
AR-251502	Premium Needle Counter with Blade Remover, Red, 20 count, 40 capacity, Sterile	80	9
AR-251701	Premium Needle Counter with Blade Remover, Red, 15 count, 30 capacity, Sterile	100	9
AR-25001	Magnetic Instrument Pad, Sterile	30	10
AR-25500NS	Sponge Counter Bag	250	11
AR-32200	Electrosurgical Tip Cleaner, Sterile	300	12
AR-7304412	Anti-Fog Solution Kit, Sterile	48	12
AR-4115175	Standard Patient Warming Blanket	50	13
AR-480161	Reflective Warming Blanket	25	13



MISSION



We are a global healthcare services company delivering exceptional value to our customers.





VISION



Together, we make healing happen by connecting the world of medical products to the point of care.

VALUES



-  We Seek Simplicity.
-  Credibility is our Currency.
-  We are Connected by Purpose.
-  Outcomes are our Difference.



“ Together, we make healing happen by connecting the world of medical products to the point of care. ”

arcroyal.ie



ArcRoyal
Virginia Road, Kells
Co.Meath, Ireland

May, 2019 v.3.0.



Printed on FSC Paper