



Every day your team drives success

HALYARD* surgical gowns support your every need.

Our
Commitment
to you

Better
Infection
Protection



A Surgical gown for every job

In the OR (Operating Room) you are routinely at risk of exposure to blood and bodily fluids. Having the right level of protection for the task at hand can make a critical difference. Halyard offers a full range of surgical gowns, each offering various levels of protection and geared towards different types of surgical procedures. With Halyard, no matter what the surgical task, our gowns have you covered.

Maximum protection – as easy as 1-2-3

Knowing that different surgical tasks require different levels of protection, Halyard uses a colour coded system so you can easily and quickly identify the right surgical gown for your immediate need.

Our gowns offer three levels of protection, each identified by a coloured neckband:

Yellow

Standard performance

Green

Standard performance with fabric reinforcement

Red

High-performance film-reinforced gowns




Finding the right **surgical gown** as easy as 1-2-3

Level of Protection

- High Performance
- Standard Performance Fabric-Reinforced
- Standard Performance


Fabric features

-  **Breathability**
-  **Fluid protection**
-  **Comfort**



● ● ● "Standard and High Performance" as defined in the European Standard EN13795-1:2019

HALYARD* core competencies




SMS

Featuring two Spunbond outer layers of polypropylene for superior strength and durability and a Meltblown middle layer to protect against bacteria and liquids, our SMS technology means low lint levels and reduced risk of contamination.

Cool shield core

This proprietary copolymer core, featured in our AeroChrome* gowns, allows moisture vapour to escape for cool comfort while also preventing fluid penetration through the fabric – even under extreme high-pressure testing (per ASTM 1670 test method).



VALUESELECT*


Economical surgical gown

More value, less money.

Halyard's VALUESELECT* surgical gowns offer optimal comfort at an economical price. Available in both standard and high performance options, our EN 13795 compliant VALUESELECT* surgical gowns are the perfect complement to the HALYARD* range. Providing an economical, quality gown option, with VALUESELECT* you can get a significant reduction to your bottom line.



VALUESELECT* Standard Surgical Gowns



STERILE CODE	NON-STERILE CODE	SIZE	LENGTH*	ITEMS/CASE STERILE	ITEMS/CASE NON STERILE
46665	46669	Medium	108 cm	26	513
46666	46670	Large	117 cm	24	486
46667	46671	X-Large	126 cm	22	486
46668	46672	XX-Large	135 cm	22	486
47285	47281	L-X Long	150 cm	20	480
47286	47282	XL-XLong	150 cm	20	480
70433	70432	XXL-XLong	155 cm	20	420

+ length back (neck to bottom)



VALUESELECT* High-Performance PE-Reinforced Surgical Gowns



STERILE CODE	NON-STERILE CODE	SIZE	LENGTH*	ITEMS/CASE STERILE	ITEMS/CASE NON STERILE
46673	46677	Medium	108 cm	22	459
46674	46678	Large	117 cm	22	432
46675	46679	X-Large	126 cm	22	432
46676	46680	XX-Large	135 cm	22	478
47287	47283	L-XLong	150 cm	20	420
47288	47284	XL-XLong	150 cm	20	360
70335	70434	XXL-XLong	155 cm	18	360

+ length back (neck to bottom)

VALUESELECT* Standard Surgical Gowns:

Our lightweight standard performance surgical gowns put an emphasis on comfort.

- SMS
- Book-fold
- Set-in-sleeves
- Low-lint fabric
- Sterile gown packaged with two towels
- Comfortable
- Lightweight
- Soft cuffs
- Long Velcro adjustable collar for maximum flexibility
- Alcohol-repellent

VALUESELECT* High-Performance PE-Reinforced Surgical Gowns:

The front and arms of our high-performance surgical gowns come reinforced with polyethylene (PE), ensuring that the critical zones are liquid-resistant.


- SMS
- Book-fold
- Set-in-sleeves
- Low-lint fabric
- Sterile gown packaged with two towels
- Comfortable
- Lightweight
- Soft cuffs
- Long Velcro adjustable collar for maximum flexibility
- Alcohol-repellent

HALYARD BASICS*

Surgical gowns



HALYARD BASICS * Standard Surgical Gowns



STERILE CODE	NON - STERILE CODE	SIZE	LENGTH*	ITEMS/ CASE STERILE	ITEMS/ CASE NON STERILE
13337	-	X-Large, XX-Long	175 cm	40	-
11783	33916	Large - Towel in Overwrap (handi-bin)	125 cm	40	440
11784	33917	X-Large - Towel in Overwrap (handi-bin)	131 cm	40	400
99294	32109	Large, X-Long - Towel in Overwrap (handi-bin)	145 cm	40	460
99295	32113	X-Large, X-Long - Towel in Overwrap (handi-bin)	150 cm	40	380
-	79284	Large (handi-bin)	125 cm	-	640
-	79294	Large, X-Long (handi-bin)	145 cm	-	560
-	79285	X-Large (handi-bin)	130 cm	-	560
-	79295	X-Large, X-Long (handi-bin)	150 cm	-	560


+ length back (neck to bottom)

HALYARD BASICS* Surgical Gowns

- Colour Coded Neck Binding: White
- SMS Fabric
- Overlapping Back
- Raglan Sleeves
- Low Linting Fabric
- Sterile Gown packaged with 2 Towels
- Meets EN13795-1:2019 Standard Performance



HALYARD BASICS * High-Performance Fabric-Reinforced Surgical Gowns



STERILE CODE	NON - STERILE CODE	SIZE	LENGTH*	ITEMS/ CASE STERILE	ITEMS/ CASE NON STERILE
13815	-	X-Large, XX-Long	175 cm	32	-

+ length back (neck to bottom)

HALYARD BASICS*

Fabric-Reinforced Surgical Gowns:

- Colour Coded Neck Binding: Green
- SMS Fabric
- Fabric-Reinforced Chest and Sleeves
- Overlapping Back
- Raglan Sleeves
- Low Linting Fabric
- Sterile Gown packaged with 2 Towels

EVOLUTION* 4

Standard Surgical Gown

Raising the bar for reliability and quality


Our soft, single-use EVOLUTION* 4 surgical gowns are derived from Halyard's breathable SMS fabric. Providing the right balance of comfort and protection during procedures involving low to moderate fluid contact, these gowns meet standard performance per EN13795-1:2019 requirements.

EVOLUTION* 4 fabric-reinforced surgical gowns:

- Ideal for shorter and mid-size procedures
- EVOLUTION* 4 SMS fabric
- Low lint
- Flame resistant
- Adjustable neckline: hook and loop
- Book fold
- Sterile gown packaged with 2 towels

EVOLUTION* 4 standard surgical gowns:


- Ideal for shorter, minimally invasive procedures
- Available in both set-in-sleeve and raglan sleeve versions



EVOLUTION * 4
Standard Surgical Gowns
Set-in Sleeves

STERILE CODE	NON-STERILE CODE	SIZE	LENGTH*	ITEMS/ CASE STERILE	ITEMS/ CASE NON STERILE
90016	72978	Large	117 cm	36	640
90046	72988	X-Large	127 cm	32	640
90019	70391	Large, X-Long	132 cm	34	66
90049	70392	X-Large, X-Long	137 cm	34	60
-	79419	Large – Towel in Overwrap	111 cm	-	512
-	79434	X-Large – Towel in Overwrap	121 cm	-	448


+ length back (neck to bottom)



EVOLUTION * 4
Fabric-Reinforced Surgical Gowns
Raglan Sleeves

STERILE CODE	NON-STERILE CODE	SIZE	LENGTH*	ITEMS/ CASE STERILE	ITEMS/ CASE NON STERILE
90112	78911	Large	123 cm	32	36
90113	78913	Large, X-Long	143 cm	30	54
90142	78921	X-Large	127 cm	30	36
90143	78923	X-Large, X-Long	148 cm	32	48
-	44763	Large, X-Long	145 cm	-	432
-	44764	X-Large, X-Long	151 cm	-	384

+ length back (neck to bottom)



EVOLUTION * 4
Standard Surgical Gowns
Raglan Sleeves

STERILE CODE	NON-STERILE CODE	SIZE	LENGTH*	ITEMS/ CASE STERILE	ITEMS/ CASE NON STERILE
90002	-	Small	115 cm	36	-
90012	79970	Large	123 cm	36	40
90042	79980	X-Large	127 cm	34	36
90013	79973	Large, X-Long	143 cm	34	60
90043	79983	X-Large, X-Long	148 cm	34	48

+ length back (neck to bottom)

ULTRA

Enhanced Standard Performance Gown


More comfort.

Halyard's ULTRA surgical gowns are designed to provide an ideal balance of comfort and protection for mid-size procedures.

- Raglan sleeves
- Low lint fabric
- Flame resistant
- Adjustable neckline: hook and loop
- Book fold
- Sterile gown packaged with 1 towel
- Meet the standard performance per EN13795-1:2019 requirements



ULTRA Standard Surgical Gowns



STERILE CODE	NON-STERILE CODE	SIZE	LENGTH ⁺	ITEMS/CASE STERILE	ITEMS/CASE NON STERILE
-	47247	Small	115 cm	-	48
95113	-	Large	123 cm	32	40
95123	-	X-Large	128 cm	30	36
95131	71331	XX-Large	135 cm	28	36
95114	71626	Large, X-Long	138 cm	32	60
95124	71627	X-Large, X-Long	144 cm	30	48

+ length back (neck to bottom)

ULTRA


Fabric-reinforced surgical gowns with neck protection & thumb loops

Speciality Surgical Gowns

- Proprietary 4th generation SMS fabric
- Extra wide fabric reinforcements
- Raglan sleeves
- Low lint fabric
- Flame resistant
- Adjustable neckline: hook and loop
- Neck protection
- Book fold
- Sterile gown packaged with 1 towel



ULTRA Fabric-reinforced surgical gowns



STERILE CODE	NON-STERILE CODE	SIZE	LENGTH ⁺	ITEMS/CASE STERILE	ITEMS/CASE NON STERILE
90363	77337	Large, X-Long	143 cm	30	48
90373	77336	X-Large, X-Long	148 cm	28	36

+ length back (neck to bottom)

AERO CHROME*

High-Performance Surgical Gown

Say goodbye to hot & heavy Surgical Gowns!

- **Red collar:** Easily identify a high-performance gown.
- **Breathable back panel:** Maximizes airflow for cool comfort.
- **Soft, light, and breathable fabric:** Keeps staff protected yet cool and comfortable.
- **Low linting:** Helps reduce the risk of wound contamination and infection.
- **Meets the High Performance per EN13795-1:2019:** Ensures superior protection.
- **Meets AAMI Level 4 standard:** Provides the highest level of fluid and microbial protection in the critical zones.
- **Sterile gown packaged with 2 towels:** Added convenience for your team.
- **Raglan Sleeves**



AERO CHROME*

Surgical Gown



STERILE CODE	NON-STERILE CODE	SIZE	LENGTH*	ITEMS/ CASE STERILE	ITEMS/ CASE NON-STERILE
48395	48392NS	Large	125 cm	30	40
48397	48399NS	Large, X-Long	145 cm	30	36
48396	48393NS	X-Large	131 cm	28	32
48398	48400NS	X-Large, X-Long	151 cm	30	32
48402	48401NS	XX-Large, X-Long	137 cm	28	30
48403	48394NS	XX-Large	137 cm	28	32

+ length back (neck to bottom)



Our Commitment to you

Better Infection Protection

For more information, please contact your local O&M HALYARD Account Manager or Customer Service on:

UK & IE (toll-free)

- 📞 0800 917 65 85
- 📞 0800 169 02 35
- ✉️ customerservice.uk.ie@hyh.com

Other countries

- 📞 +32 2 700 68 51
- 📞 +32 2 711 26 91
- ✉️ customerservice.export@hyh.com