

PROMOTING STERILITY ENHANCING PROTECTION

High-Efficiency Materials for Sterilisation.



BUILDING TRUST IN STERILISATION

**When it comes to your sterilisation packaging,
preserving sterility means protecting your patient.**

At Halyard, we believe that there should be a sterilisation wrap available for every need and for every user.

Within the HALYARD* sterilisation wrap portfolio we offer:

- HALYARD* Sequential and Quick Check INTERLEAVED* Sterilisation Wrap: single sheets of POWERGUARD*-treated SMS for sequential wrapping
- HALYARD* Smart-Fold* Wrap: Designed with your toughest tasks in mind. Its impact-resistant fabric and triple-layer reinforcement zones are engineered to help protect your heaviest trays and loaner sets from tears and cuts in handling.
- HALYARD* QUICK CHECK* Sterilisation Wrap: a POWERGUARD* - treated ONE-STEP* product with a blue and a white layer allowing you to detect more easily if sterility has been breached.



H100



H200



H300



H400



H500

Next to the HALYARD* Sequential and Quick Check INTERLEAVED* and Smart-Fold* Wraps, Halyard is offering you a range of high quality accessories and infection prevention solutions that you need for enhanced efficiency and protection.



PURPLE NITRILE* XTRA and PURPLE NITRILE* MAX: These thick, extended-length gloves are specifically designed to provide CSSD technicians with robust protection during instrument decontamination.



FLUIDSHIELD* LEVEL 3 Procedure Masks: High quality masks that provide clear, reliable performance with no distractions, helping you perform at your best.

PACKAGING PRODUCTS



Heat Seal Pouches and Reels:

A comprehensive range of pouches and reels produced with a 68gsm high burst resistance medical grade paper, which is sealed to tear resistant film to provide excellent seal strength and good peel characteristics.

- Suitable for steam and ETO sterilisation
- Compliant with ISO 11607 & EN 868-5



HALYARD* HEAT-SEAL PAPER/FILM POUCHES

Code	Description	Case Count
48516	HALYARD* Heat-Seal Paper/Film Pouches 50 x 250 mm	2.000
48518	HALYARD* Heat-Seal Paper/Film Pouches 75 x 200 mm	2.000
48520	HALYARD* Heat-Seal Paper/Film Pouches 200 x 400 mm	2.000
48495	HALYARD* Heat-Seal Paper/Film Pouches 100 x 150 mm	2.000
48498	HALYARD* Heat-Seal Paper/Film Pouches 100 x 270 mm	2.000
48499	HALYARD* Heat-Seal Paper/Film Pouches 100 x 300 mm	1.000
48500	HALYARD* Heat-Seal Paper/Film Pouches 100 x 380 mm	2.000
48501	HALYARD* Heat-Seal Paper/Film Pouches 100 x 400 mm	2.000
48503	HALYARD* Heat-Seal Paper/Film Pouches 120 x 400 mm	1.000
48506	HALYARD* Heat-Seal Paper/Film Pouches 150 x 300 mm	1.000
48507	HALYARD* Heat-Seal Paper/Film Pouches 150 x 380 mm	2.000
48558	HALYARD* Heat-Seal Paper/Film Pouches 150 x 350 mm	1.000
48508	HALYARD* Heat-Seal Paper/Film Pouches 150 x 400 mm	2.000
48509	HALYARD* Heat-Seal Paper/Film Pouches 200 x 610 mm	1.000
48510	HALYARD* Heat-Seal Paper/Film Pouches 205 x 400 mm	1.000
48511	HALYARD* Heat-Seal Paper/Film Pouches 220 x 300 mm	1.000
48512	HALYARD* Heat-Seal Paper/Film Pouches 250 x 380 mm	1.000
48513	HALYARD* Heat-Seal Paper/Film Pouches 250 x 500 mm	500
48514	HALYARD* Heat-Seal Paper/Film Pouches 300 x 500 mm	500
48515	HALYARD* Heat-Seal Paper/Film Pouches 420 x 500 mm	500



HALYARD* HEAT-SEAL PAPER/FILM REELS

Code	Description	Case Count
48526	HALYARD* Heat-Seal Paper/Film Reels 75 mm x 200 m	2
48521	HALYARD* Heat-Seal Paper/Film Reels 100 mm x 200 m	1
48522	HALYARD* Heat-Seal Paper/Film Reels 150 mm x 200 m	1
48523	HALYARD* Heat-Seal Paper/Film Reels 200 mm x 200 m	1
48524	HALYARD* Heat-Seal Paper/Film Reels 250 mm x 200 m	1
48555	HALYARD* Heat-Seal Paper/Film Reels 300 mm x 200 m	1

HALYARD* HEAT-SEAL PAPER/FILM GUSSETED POUCHES/REELS

Code	Description	Case Count
48492	HALYARD* Heat-Seal Paper/Film Gusseted Pouches/Reels 150 x 50 x 400 mm	1.000
48494	HALYARD* Heat-Seal Paper/Film Gusseted Pouches/Reels 150 mm x 50 mm x 100 m	2

PROTECTING YOUR INSTRUMENT TRAY FROM THE INSIDE



ABSORBENT TRAY LINERS

Absorbent tray liners are designed to effectively absorb and distribute condensate during the sterilization process. They can also be used as general-purpose towels in healthcare environments.

Code	Description	Case Count
48139	HALYARD® Absorbent Trayliner 51 x 64 cm	400



MILLENNIUM TRAY LINERS (40GSM)

The Millenium blue liners are lightweight, user-friendly, and designed with a matte blue finish to reduce glare. Made from a 90% viscose and 10% polyester blend, these liners can absorb several times their own weight in water, locking it into their porous structure to prevent moisture from impacting the efficiency of the sterilization process. Typically placed inside instrument baskets, Millenium tray liners provide cushioning and absorption while also reducing glare from the instruments.

Code	Description	Case Count
E10011-26	Millenium Tray liner 25 x 25 cm	2.000
E10011-27	Millenium Tray liner 25 x 32 cm	2.000
E10011-31	Millenium Tray liner 28 x 50 cm	1.000



PROTECTIVE TRAY LINERS (60GSM)

Code	Description	Case Count
K100111	Protective trayliner 28 x 50 cm	500
K100112	Protective trayliner 25 x 33 cm	500



Tip from our clinical specialist: how do I choose my tray liner?

	Where	Indication
Millenium tray liners (40gsm)	Mainly used inside instrument basket	Absorption/distribution of wetness ++ Protection of instruments ++
HALYARD® absorbent tray liners (60gsm)	Can be used in or under instrument basket	Absorption/distribution of wetness +++ Reinforces bottom of packaging +
Protective tray liners (60gsm)	Can be used in or under instrument basket	Absorption/distribution of wetness + Reinforces bottom of packaging +++

PROTECTING YOUR INSTRUMENT TRAY FROM THE OUTSIDE



PURPLE TRANSPORT BAGS

Purple transport bags into which barrier-wrapped trays are inserted. The transport bag can be closed with the lip provided and with sterilizer tape.

The material is fully steam permeable.

Code	Description	Case Count
E10011-23	Transport Bag 40 x 70 cm	100
E10011-24	Transport Bag 50 x 70 cm	120
E10011-25	Transport Bag 60 x 70 cm	150



SILICONE TRANSPORT CARPET

HALYARD® Silicon Transport Carpet ensures shock absorption and reduces friction between instrument trays and transport baskets, keeping trays secure during transport. Its anti-slip surface prevents sliding or dragging and is designed for use in steam sterilizers, with an open structure that avoids wet packs.

Code	Description	Case Count
30010-00004	Silicone transport carpet 28 x 54 cm	10



TRANSPORTATION TRAY

HALYARD® Transport Trays are an aid to proper handling techniques of sterile packs. The Transport Trays are used to lift and support the sterile pack whilst it's in storage between sterilisation and use. HALYARD® Transport Trays are reusable after cleaning and decontamination. The XL transport tray is used to lift, and support wrapped robotic/heavy trays.

Code	Description	Case Count
KL-1500	Transportation Tray 29 x 55 cm	10
KL-1500XL	XL Transportation Tray 43 x 69.7 cm	5

INSPECTION



Tip from our clinical specialist: sometimes my instruments remain wet after automated cleaning.

You can use our low-linting single-use wipes: they can be used in the wash room to remove moisture from manually-cleaned instruments. Keep some in the assembly and packaging room so that you can remove the excess moisture droplets from instruments after automated cleaning and disinfection.

Product Codes	Description	Units/package
CHUX-WIPE01	Chux low linting absorbent towels, turquoise Folded dimensions: 155 x 170 mm Opened out: 300 x 337 mm	400
CHUX-VERA02A	Chux low-linting absorbent towels noncrepe, white Folded dimensions: 155 x 170 mm Opened out: 300 x 337 mm	300



MARKER PENS

Permanent marker pens, effective for writing on the surface of tray tags, pouches and reels as their ink is suitable for sterilizer tape, non-toxic and quick-drying.

Code	Description	Case Count
RT209P	Marker pen black	12



STERILISER TAPE

A complete range of unique, manufactured tapes which demonstrate excellent binding characteristics and comply with various types of sterilisation wraps.

Code	Description	Case Count
A3341CT0025B	Masking tape without indicator 25mm wide	36
A3352CT0019B	Latex free sterilizer tape with indicator 19 mm wide	48
A3352CT0025B	Latex free sterilizer tape with indicator 25 mm wide	36

References 1. Study Kurt Salmon Associates.



For more information,
please contact your O&M Halyard Customer Care
Team:

UK & IE (toll-free)

Tel: 0800 085 82 04 - Fax: 0800 066 91 06
customerservice.uk.ie@owens-minor.com

Other countries

Tel: +32 2 585 51 14 - Fax: +32 2 897 00 95
customerservice.export@owens-minor.com



www.halyardhealth.eu

*Registered Trademark or Trademark of O&M Halyard or its affiliates.
©2025. All rights reserved. Confidential and Proprietary to Owens &
Minor, Inc. and/or its affiliated companies. HC488-02-EN GB
The picture on the front cover is used with the permission of BELINTRA
www.belintra.com