

SIMPLE SOLUTIONS FOR HIGH-TECH ROBOTIC INSTRUMENTS



Once looked at as futuristic, today, robot-assisted surgery is now part of the surgical landscape. Packaging and storage of high-tech cameras represents a new challenge for the CSSD.

This issue can be solved combining the HALYARD* XL-transportation with SMART-FOLD* or QUICK CHECK* sterilisation wrap.

Due to its size (OD: 43 x 69.7cm; ID: 41.9 x 67.6cm), the XL transportation tray can accommodate camera trays used in robot-assisted surgery.

The XL polypropylene transport tray is used to lift, and support wrapped robotic trays. The transport trays are reusable after cleaning and decontamination.



The robotic tray is supported over its entire length, reducing the risk of damage to the sterile packaging.

Ordering Information

PRODUCT DESCRIPTION		PRODUCT CODE	
	XL transportation tray	5 trays/case	KL-1500XL
PARALLEL PACKAGING			
	137x182cm HALYARD* QUICK CHECK* sterilisation wrap	H400	34191
		H500	34195
ENVELOPE PACKAGING			
	137x137cm HALYARD* QUICK CHECK* sterilisation wrap	H400	34146
		H500	34143
SMART-FOLD*			
	HALYARD* SMART-FOLD* sterilisation wrap	H650	14309

For more information, please send an email to customerservice.uk.ie@hyh.com or visit

www.halyardhealth.eu



O&M Halyard, Inc., 9120 Lockwood Blvd., Mechanicsville, VA 23116, USA
 [EC] [RE] ArcRoyal Unlimited Company, Virginia Road, Kells, Co Meath, Ireland
 O&M Halyard Belgium BV, Berkenlaan 8B, 1831 Machelen (Brab.), Belgium
 UK RP • RLD Quality Limited Centenary House, Peninsula Park, Rydon Lane, Exeter, EX2 7XE, United Kingdom



*Registered Trademark or Trademark of O&M Halyard or its affiliates.
 ©2024 O&M Halyard. All rights reserved. HC 506-01-UK