

HALYARD* AND BELINTRA® SMART-FOLD* STERISYSTEM®
THIS CHANGES EVERYTHING.

A two-touch sterilization packaging, storage and transport system.



CATALOG OF PRODUCTS

CONTENTS

4 WRAP

SMART-FOLD* SELECT
SMART-FOLD* HIGH PERFORMANCE

7 STORAGE

STATIONARY

UBEFLEX®

UBEFLEX® ACCESSORIES

HIGH DENSITY MOBILE

UBETRACK™

UBETRACK™ ACCESSORIES

OPEN SHELVING

19 DRY BASE® INSTRUMENT TRAYS

21 TRANSPORTATION

OPEN CARTS

OPEN STERICART (CART WASHABLE)

MODULAR CART (NOT CART WASHABLE)

CLOSED CART

CLOSED STERICART (CART WASHABLE)

CART ACCESSORIES



PRESENTING THE SMART-FOLD* STERISYSTEM®

Change is good! That's why Halyard and Belintra have teamed up to bring you the SMART-FOLD* STERISYSTEM®, a unique system that combines HALYARD* SMART-FOLD* Sterilization Wrap with BELINTRAS storage and transport system. This two-touch solution not only ensures that sterility has been maintained from sterilization through to use in the OR, it confirms it! And THIS changes everything!

That's because the long-held belief that rigid containers are infallible just isn't true. They cause untold problems and headaches. Everything from failure to maintain sterility to inefficiencies in labor, design and space.

Created to solve these problems, the SMART-FOLD* STERISYSTEM® eliminates stacking, which reduces the chance for tears, cuts and holes. This enhances clinical confidence and reduces costs, time and labor associated with reprocessing.

It also increases efficiencies, is more ergonomic and features aseptic opening. And now with our new SMART-FOLD* Select Sterilization Wrap combined with BELINTRAS storage and transport system; we have an economical solution that is much preferred over problematic containers.

For the utmost in patient safety and the ultimate in long-term savings, this is the system for you. That's because it addresses four critical areas of focus:



2-Touch System
and
Aseptic Opening



Better
Throughput



Better
Ergonomics



Save on Space
and Storage

SMART-FOLD* STERILIZATION WRAP



SMART-FOLD* STERILIZATION WRAP

HALYARD* STERILIZATION WRAP HAS BEEN TESTED AND VALIDATED FOR USE WITH THE STERISYSTEM.®

SMART-FOLD* provides unique time saving and clinically beneficial features that make it the ideal wrap to be used with the efficient Sterisystem.® SMART-FOLD* allows Sterile Processing Professionals to process their heaviest trays with the utmost confidence that sterility will be maintained to help protect what matters most, the patient.

SMART-FOLD* HIGH PERFORMANCE

SMART-FOLD* HIGH PERFORMANCE Sterilization Wrap is Halyard's highest grade, and most durable wrap. We recommend this for your heaviest sets.

SMART-FOLD* SELECT

SMART-FOLD* SELECT Sterilization Wrap is expressly optimized at a lighter grade to work within the SMART-FOLD* STERISYSTEM® for trays weighing up to 25lbs.

In a recent user study, CS and OR clinicians found SMART-FOLD* HIGH PERFORMANCE Sterilization Wrap to be much better when compared to traditional wrap¹.

**ERGONOMIC WRAPPING DESIGN**

82% of the CS clinicians find SMART-FOLD* more ergonomic with less physical stress during wrapping.¹

**ASEPTIC UNWRAPPING**

88% of clinicians find the design of SMART-FOLD* better for aseptic opening.¹

**CONFIDENCE IN STERILITY**

78% of clinicians are more confident that sterility has been maintained with SMART-FOLD*.¹

SMART-FOLD* STERILIZATION WRAP

THE ASSURANCE YOU NEED.
THE EFFICIENCY YOU APPRECIATE.

Depend On SMART-FOLD* for Multifaceted Performance.

STRONGER

Reinforced zones protect tray edges and corners.

FASTER WRAPPING

Positioning Strips hold material for easier more ergonomic wrapping.

MORE ACCURATE

Reference Lines serve as guides for tray placement.

OPTIMIZED AIRFLOW

Patented Daisy Bond Pattern helps optimize available surface area for sterilant penetration and evacuation. Porosity allows for optimum airflow in and out of package, reducing time-consuming reprocessing and helping to keep your OR on schedule.








FASTER INSPECTION

Patented QUICK CHECK* Technology featuring dual contrasting layers bonded together make inspecting trays easier and faster than ever!

SMARTER

Pull Tabs improve aseptic presentation.

SMART-FOLD* Sterilization HIGH PERFORMANCE and SELECT Wrap					
HIGH PERFORMANCE Product Codes	Sheets per case	SELECT Product Codes	Sheets per case	Color Code by Size	Dimensions
14271	48	14268	48		40 inches wide x 47 inches long
14277	48	14275	48		40 inches wide x 55 inches long
14281	48	14279	72		28 inches wide x 46 inches long
14309	36	14313	48		48 inches wide x 61 inches long
14320	48	14324	72		22 inches wide x 45 inches long

STORAGE SOLUTIONS




STORAGE SOLUTIONS

TWO STORAGE SYSTEMS TO CHOOSE FROM:



UBeFLEX® Stationary Storage System



UBeTrack™ Mobile Storage System



UBEFLEX® STATIONARY STORAGE SYSTEM




Designed for quick assembly, functionality, durability and stability, the unique UBeFlex® modular stationary storage system is made from anodized aluminum and other high-quality, anti-static synthetic materials, which means it can be used in both sterile storage and non-sterile environments.

Its modular construction, in various column configurations, may be combined with open shelving so it can be utilized in multiple areas of your facility. And this ergonomic and hygienic storage solution can also be used as a "pass-through" rack.

For savings on space and storage and enhanced efficiency, this unit is ideal.



For assembly, the UBeFlex® modular stationary storage system includes just 2 posts and a spacer set. The spacer set features 4 connectors that connect the posts. To expand on this configuration, you simply need an extra post and spacer set. This flexibility allows it to be adapted to the size and shape of any room.

GUIDE FOR NUMBER OF SETS PER COLUMN, DETERMINED BY TRAY HEIGHT		
Up to 4.5" 10 Trays	Between 4.5" and 7" 8 Trays	Between 7" and 9.5" 6 Trays
		

COMPONENTS NEEDED TO BUILD UBEFLEX® COLUMN TYPES			
Component	DIN Deep Type Column	ISO Deep Type Column	ISO Wide Type Column
Posts	47360	47360	47489
Spacer Set	47365	47364	47490
Roof/Bottom Panel	47367	47366	47491
Runners	47369	47368	47495
Perfo-Safe® Shelves	47404	47405	47405

UBEFLEX® COLUMNS



POSTS

Code	Description	Case Quantity
47360	Deep Type (77"H x 25"D, 24" Depth)	2
47489	Wide Type (77"H x 18"D, 16" Depth)	2



SPACER SETS

Code	Description	Case Quantity
47365	DIN Deep Type (12" Wide)	1 Set
47364	ISO Deep Type (16" Wide)	1 Set
47490	ISO Wide Type (24" Wide)	1 Set




ROOF OR BOTTOM PANEL

Code	Description	Case Quantity
47367	DIN Deep Type (24"D x 12"W)	2
47366	ISO Deep Type (24"D x 16"W)	2
47491	ISO Wide Type (16"D x 24"W)	2



PERFO-SAFE® STAINLESS STEEL WIRE SHELVES

Code	Description	Case Quantity
47404	DIN Deep Type (24"D x 12"W)	10
47405	ISO Deep Type (24"D x 16"W) or ISO Wide Type (16"D x 24"W)	10



SILICONE MAT

Code	Description	Case Quantity
47408	Silicone Mat	10

UBEFLEX® COLUMNS (CONT'D)



HIGH PRESSURE LAMINATE SHELVES (HPL)

Code	Description	Case Quantity
48153	DIN Deep Type (24"D x 12"W)	5
48154	ISO Deep Type (24"D x 16"W) or ISO Wide Type (16"D x 24"W)	5



RUNNERS

Code	Description	Case Quantity
47369	DIN Deep Type	12 Pairs
47368	ISO Deep Type	12 Pairs
47495	ISO Wide Type	12 Pairs

47368, 47495



47369

UBEFLEX® ACCESSORIES



LABEL HOLDER AND INSERTS

HYH	Description	Case Quantity
47370	Label Holder	100
47372	Label Insert (Plastic White)	100



COLOR IDENTIFICATION STRIPS FOR POSTS

HYH	Description	Case Quantity
47492	Red	10
47493	Blue	10
47494	Green	10
47496	Orange	10



PERFO-SAFE® STAINLESS STEEL BASKETS


Code	Description	Case Quantity
47409	DIN Deep Low (24"D x 12"W x 4.5"H)	6
47410	DIN Deep High (24"D x 12"W x 6.5"H)	6
47411	ISO Deep Medium (24"D x 16"W x 4"H) or ISO Wide Medium (16"D x 24"W x 4"H)	6
47412	ISO Deep Medium High (24"D x 16"W x 5"H) or ISO Wide Medium High (16"D x 24"W x 5"H)	6



POLYCARBONATE BASKETS AND DIVIDERS

All polycarbonate baskets are ISO sized (24"D x 16"W) or (24"W x 16"D)

Code	Description	Case Quantity
47393	2"H Basket (Completely Closed)	10
47394	2"H Long Divider 24", Long	40
47395	2"H Short Divider 16", Short	40
47436	4"H Basket (Closed Base, Open Sides)	10
47437	4"H Long Divider	20
47438	4"H Short Divider	20
47439	8"H Basket (Closed Base, Open Sides)	6
47440	8"H Long Divider	10
47441	8"H Short Divider	10



UBETRACK™ MOBILE STORAGE SYSTEM

The dynamic UBeTrack™ is a high-density, space-saving mobile storage system with adjustable shelving.

Efficient and user-friendly this system is very lightweight, making the mobile racks quick and easy to move. Made from anodized aluminum and high-quality, antistatic synthetic materials, the unique UBeTRACK™ can be used in both sterile storage and non-sterile environments.

By placing the runner system at the top, cleanliness is maintained, and the unit is easy to pass with transport trolleys.

UBETRACK® RACKS



Mobile Post
pictured above

POSTS

Code	Description	Case Quantity
47383	Deep Type Stationary Post for Outer Rack (80"H x 25"D)	2
47388	Deep Type Mobile Post for Inner Rack (80"H x 27"D)	2

1 case = 2 per mobile unit



ASSEMBLY KITS

Code	Description	Case Quantity
47389	Kit (4"H x 5"D x 16"W)	1 Pair

1 case per stationary rack



TRACK KITS

Code	Description	Case Quantity
47391	2 Inner Rack Kit (3"H x 4.5"D x 59"W)	1
47392	3 Inner Rack Kit (3"H x 4.5"D x 73"W)	1

2 cases per UbeTrack™ unit (one for each side)



GUIDE ROLLER SETS

Code	Description	Case Quantity
47390	Roller set	1

1 case per mobile rack



SPACER SETS

Code	Description	Case Quantity
47365	DIN Deep Type (12" Wide)	1 Set
47364	ISO Deep Type (16" Wide)	1 Set

UBeTrack™ | Mobile Storage



ROOF OR BOTTOM PANEL

Code	Description	Case Quantity
47367	DIN Deep Type (24"D x 12"W)	2
47366	ISO Deep Type (24"D x 16"W)	2



47368

RUNNERS

Code	Description	Case Quantity
47369	DIN Deep Type	12 Pairs
47368	ISO Deep Type	12 Pairs




47369



PERFO-SAFE® STAINLESS STEEL WIRE SHELVES AND SILICONE MATS

Code	Description	Case Quantity
47404	DIN Deep Type (24"D x 12"W)	10
47405	ISO Deep Type (24"D x 16"W)	10



SILICONE MAT

Code	Description	Case Quantity
47408	Silicone Mat	10



HIGH PRESSURE LAMINATE SHELVES (HPL)

Code	Description	Case Quantity
48153	DIN Deep Type (24"D x 12"W)	5
48154	ISO Deep Type (24"D x 16"W)	5

UBETRACK™ ACCESSORIES



VERTICAL HANDLES FOR INNER RACKS

HYH	Description	Case Quantity
47505	Handles	1 Pair
48134	Handle Assembly Kit	1

1 Assembly Kit needed for each pair of handles.



LABEL HOLDER AND INSERTS

HYH	Description	Case Quantity
47370	Label Holder	100
47372	Label Insert (Plastic White)	100



COLOR IDENTIFICATION STRIPS FOR POSTS

HYH	Description	Case Quantity
47492	Red	10
47493	Blue	10
47494	Green	10
47496	Orange	10



PERFO-SAFE® STAINLESS STEEL BASKETS

Code	Description	Case Quantity
47409	DIN Deep Low (24"D x 12"W x 4.5"H)	6
47410	DIN Deep High (24"D x 12"W x 6.5"H)	6
47411	ISO Deep Medium (24"D x 16"W x 4"H)	6
47412	ISO Deep Medium High (24"D x 16"W x 5"H)	6



POLYCARBONATE BASKETS AND DIVIDERS

All polycarbonate baskets are ISO sized (24"D x 16"W)

Code	Description	Case Quantity
47393	2"H Basket (Completely Closed)	10
47394	2"H Long Divider 24", Long	40
47395	2"H Short Divider 16", Short	40
47436	4"H Basket (Closed Base, Open Sides)	10
47437	4"H Long Divider	20
47438	4"H Short Divider	20
47439	8"H Basket (Closed Base, Open Sides)	6
47440	8"H Long Divider	10
47441	8"H Short Divider	10

Open Shelving



POSTS

Code	Description	Case Quantity
47360	Deep Type (77"H x 25"D, 24" Depth)	2
47489	Wide Type (77"H x 18"D, 16" Depth)	2



SPACER SETS

Code	Description	Case Quantity
47508	36"	1 Set
47509	48"	1 Set
47510	60"	1 Set
47511	72"	1 Set



POLYMER SHELVES

Code	Description	Case Quantity
47506	Deep Type (24"D x 12"W)	12
47507	Wide Type (24"D x 16"W)	12



CHROME WIRE SHELF & DIVIDER

Code	Description	Case Quantity
47363	Deep Type (24"D x 12"W)	12
47382	Divider (24"D x 9"H)	12



LABEL HOLDER AND INSERT

Code	Description	Case Quantity
47371	Label Holder	50
47372	Label Insert (Plastic White)	100

INSTRUMENT TRAYS & LIDS



Dry-Base® Stainless Steel Instrument Trays and Lids

Size Guide



STERISYSTEM™ DRY-BASE® INSTRUMENT TRAYS

Our Sterisystem™ Dry-Base® instrument trays are manufactured from highly polished grade stainless steel.

The smooth, dry base completely traps any condensation moisture that may stick to the wrap. This ensures that condensation moisture can heat up or evaporate under the influence of radiation, convection or conduction.

Our instrument trays have a smooth finished top wire. This reduces the chance of perforating or tearing the wrap.



DIN

TRAY AND LID SOLD SEPARATELY (BULK)

Code	Dimensions	Case Quantity
47413	Low/Small Size Base (4 3/4 x 9 13/16 x 2 1/2)	12
47428	Small Size Lid (4 3/4 x 9 13/16 x 13/16)	6
47414	Low/Mid Size Base (9 7/16 x 9 13/16 x 2 1/2)	12
47417	High/Mid Size Base (9 7/16 x 9 13/16 x 3 15/16)	6
47442	Mid Size Lid (9 7/16 x 9 13/16 x 13/16)	6
47416	Low/Full Size Base (18 7/8 x 9 13/16 x 2 1/2)	6
47418	High/Full Size Base (18 7/8 x 9 13/16 x 3 15/16)	6
47429	Full Size Lid (18 7/8 x 9 13/16 x 13/16)	6

TRAY AND LID SOLD TOGETHER (EACH)

Code	Description	Case Quantity
47914	Low/Small Size Lid & Base (5 3/4 x 10 1/8 x 2 7/8)	1
47920	Low/Mid Size Lid & Base (9 49/64 x 10 1/8 x 2 7/8)	1
47921	High/Mid Size Lid & Base (9 49/64 x 10 1/8 x 4 3/8)	1
47922	Low/Large Size Lid & Base (14 31/64 x 10 1/8 x 2 7/8)	1
47915	Low/Full Size Lid & Base (19 11/64 x 10 1/8 x 2 7/8)	1
47916	High/Full Size Lid & Base (19 11/64 x 10 1/8 x 4 3/8)	1

ISO

TRAY AND LID SOLD SEPARATELY (BULK)

Code	Description	Case Quantity
47419	Low/Small Size Base (4 3/4 x 11 13/16 x 2 1/2)	12
47430	Small Size Lid (4 3/4 x 11 13/16 x 13/16)	6
47425	Low/Mid Size Base (9 7/16 x 11 13/16 x 2 1/2)	12
47431	Mid Size Lid (9 7/16 x 11 13/16 x 13/16)	6
47426	Low/Full Size Base (18 7/8 x 11 13/16 x 2 1/2)	6
47432	Full Size Lid (18 7/8 x 11 13/16 x 13/16)	6

TRAY AND LID SOLD TOGETHER (EACH)

Code	Description	Case Quantity
47917	Low/Small Size Lid & Base (5 3/64 x 12 3/12 x 2 7/8)	1
47918	Low/Mid Size Lid & Base (9 49/64 x 12 3/12 x 2 7/8)	1
47919	Low/Full Size Lid & Base (19 11/64 x 12 3/12 x 2 7/8)	1

Dry-Base® Instrument Trays and Lids are not currently available for sale through Halyard in the USA.

TRANSPORTATION



Sterisystem® Transport Carts

THREE CARTS TO CHOOSE FROM:



Open Stericart




Modular Cart



Closed Stericart

Sterisystem® Open Carts



OPEN STERICARTS (CART WASHABLE)

The open Stericart in stainless steel is designed for indoor transport of both instrument sets and consumables between the SPD and the OR.

The cart is available in both an ISO size (24 x 16 inches) and a double DIN size (24 x 24 inches). When using the double DIN version, the capacity per level is doubled if your facility uses DIN transport baskets or trays.

The yellow push handles are a unique indicator that the cart may safely go through the cart washer.



OPEN STERICARTS STAINLESS STEEL (CART WASHABLE)

Code	Description	Case Quantity
47397	ISO Deep Type (62"H x 28"D x 24"W) (24"D x 16"W)	1
47469	ISO Wide Type (62"H x 20"D x 31"W) (16"D x 24"W)	1



Code	Description	Case Quantity
47396	Double DIN (62"H x 28"D x 31"W) (24"D x 24"W)	1

Sterisystem® Open Carts



MODULAR CARTS (NOT CART WASHABLE)

Our ISO transport carts, based on ISO 16 x 24 inch, have been specifically designed to transport goods between various units.

They allow logistics staff, pharmacists and nursing staff to easily transfer trays from the storage system to the carts and vice versa.

The runners consist of U and L profiles. They are compatible with all the common ISO trays, baskets and shelves. Modular baskets or trays simply slide into the U profile; the L profile has been designed to accommodate baskets with a high edge.



Code	Description	Case Quantity
47810	ISO Deep Type (71"H x 28"D x 21"W)	1



Code	Description	Case Quantity
47812	ISO Wide Type (71"H x 20"D x 29"W)	1



Code	Description	Case Quantity
47811	Double ISO Deep (71"H x 28"D x 38"W)	1

Sterisystem® Closed Carts



CLOSED STERICARTS (CART WASHABLE)

The closed Stericart features washing-tunnel-resistant stainless steel, is manufactured based on the DIN standard dimensions of 12 x 24 inches and can be used with 2 DIN positions one after the other. The 2-inch extra depth at the back allows this transport cart to also be used for sterile transport of slightly longer sets.

The lower model (43 inch) closed Stericart can double as a sterile field in the operating room, while the "Extra-Wide" model can serve to transport packages up to 37 x 24 inches.



Fixed durable yellow plastic runners, bumpers, push handles and latches set this cart apart from the rest. The well thought out design is ergonomic and allows for enhanced space savings when transporting sterile instrument sets to the OR and contaminated sets back to the SPD.

STERICART STAINLESS STEEL CARTS (CART WASHABLE)



LOW CART

Code	Description	Case Quantity
47399	Low (44"H x 29"D x 29"W) (24" x 24")	1



HIGH CART

Code	Description	Case Quantity
47398	High (60"H x 29"D x 29"W) (24" x 24")	1



EXTRA-WIDE

Code	Description	Case Quantity
47403	Extra-Wide High (60"H x 29"D x 43"W) (37" x 24")	1


Stericart Accessories

PERFO-SAFE® STAINLESS STEEL WIRE SHELVES AND SILICONE MATS



PERFO-SAFE® STAINLESS STEEL WIRE SHELVES AND SILICONE MATS

Code	Description	Case Quantity
47404	DIN Deep Type (24"D x 12"W)	10
47406	Double DIN Type (24"D x 24"D)	10
47405	ISO Deep Type (24"D x 16"W) or ISO Wide Type (16"D x 24"W)	10
47487	Extra Wide Type (24D" x 37"W) (Upper Compartment of 47403)	6



SILICONE MAT

Code	Description	Case Quantity
47408	Silicone Mat	10





For more information or an evaluation of how your facility can benefit from this unique Sterilization Packaging and Storage and Transport solution, contact:

your Halyard Health representative, call **1-844-425-9273**
or visit **www.halyardhealth.com/find-a-sales-rep**

1 Data on file with Halyard. Based on a 2018 survey of 12 facilities with 61 randomly-selected clinicians representing the Central Sterile and Operating Room departments by Halyard sales personnel. Existing Maintenance of Package Integrity (MPI) reports on file with Halyard support the passing results of the sterility for HALYARD* SMART-FOLD* wrap.



www.halyardhealth.com
www.belintra.com/us

* Registered Trademark or Trademark of O&M Halyard or its affiliates. © 2019. All rights reserved.
UBeFlex® & Sterisystem® are Registered Trademarks or Patent Pending Trademarks of Belintra NV
COPY-03607 008095